

Improved Middle-Down Characterization of Antibodies Using Multiple Ion Activation Techniques and Ion-Ion Proton Transfer Reactions on a Modified Orbitrap Mass Spectrometer

John E.P. Syka, Christopher Mullen, Stephane Houel, Romain Huguet, Joshua A. Silveira, Helene L. Cardasis, and Vlad Zabrouskov; Thermo Fisher Scientific 355 River Oaks Parkway, San Jose, CA 95134

ABSTRACT

This is a preliminary investigation of utilizing ion-ion proton transfer reactions (IIPPT) subsequent to the different modes of fragmentation used in middle-down MS/MS analysis of monoclonal antibodies (mAbs). This study demonstrates that use of IIPPT reactions for *m/z* selected ranges of product ions from short electron transfer dissociation (ETD) reactions enables observation of long mAb subunit product ions.

INTRODUCTION

The sub-unit mass analysis of monoclonal antibodies is a common assay. However obtaining critical sequence information via "middle-down" MS/MS analyses wherein intact mAb sub-units Fc/2, LC and Fd are directly dissociated in an MS/MS experiment is still considered a challenge. Last year, we described LC-MS/MS approaches utilizing a broad precursor *m/z* range to *m/z* select several charge states of eluting mAb subunits for collective MS/MS analysis. The described assays involved a multiplicity of ion dissociation types: ETD and Ultra Violet Photo-dissociation (UVPD) and ETD reaction followed by collision cell type collisional dissociation (ETnCD). The best result was obtained from the aggregation of data from twelve LC MS/MS (15 minutes) analyses covering 6 different reaction/activation times for both ETD and UVPD. Sequence coverages ranging from 83% to 92% for the subunits of NIST mAb standards were obtained. However, study of the sequence maps indicated that the sequence coverage was obtained with few long sequence ions and not many complementary N- and C-terminal pairs of product ions were observed.

Previous reports indicated that IIPPT reactions subsequent to MS² ETD activation leads to increased sequence coverage.² At this conference we and one of our collaborators in this area have presentations demonstrating similar advantage in applying IIPPT reactions subsequent to UVPD^{3,4}.

As the molecular weight of apolypeptide precursor ions increases, so too does the number of possible product ions, and considerable spectral congestion occurs in the vicinity of the precursor ion *m/z* and intact charge reduced precursor ions (ETnCD products) preventing the observation of any discrete product ion isotopic clusters from large product ions species. For ETD and UVPD, if the reaction/activation times are extended, multiple generations of product ions are generated, leading to a population of relatively short and low charge sequence ions of that are more widely distributed in the *m/z* domain. IIPPT reactions differentially spread product ions across the *m/z* range without the necessity of over reacting/activating, thus conserving the large product fragments. Highly charged product ions undergo multiple sequential IIPPT reactions (ion-ion reaction rate constants vary as *z*²) and "diffuse" away in *m/z* space from the adjacent lower *m/z* product ions. Here we present results from a preliminary investigation demonstrating the potential for IIPPT of selected *m/z* windows of mAb subunit product ions to increase mAb sequence coverage of mAb subunits and allow observation of large sequence ion.

MATERIALS AND METHODS

Sample Preparation: The NIST mAb standard (RM 8671) was used for all analyses. It was first digested with IdeS protease and then reduced under denaturing conditions (Guanidine/DTT) to generate three ~25 kDa subunits (Fc/2, LC and Fd).

LC Separations: A Thermo Scientific™ UltiMate™ 3000 HPLC system was used for all separations (Solvent A: Water with 0.1% Formic Acid (v/v); Solvent B: Acetonitrile with 0.1% Formic Acid (v/v); flow rate: 1.5µl/min; Column 25 cm × 100 µm ID RP-4H; Sample load: 100 ng). The LC column was deliberately overloaded to increase analyte signal abundance and extend the time for acquisition of MS² and MS³ spectra.

MS Instrumentation and Methods: Targeted LC-MS experiments were performed using a modified Thermo Scientific™ Orbitrap Fusion Lumos™ Tribid™ MS with ETD capability. A second reagent inlet was added to the ETD source so that the IIPPT reagent, perfluoroperhydrophenanthrene (623 *m/z*, C₁₆F₂₄), as well as the standard ETD reagent, fluoranthene, (*m/z* 202, C₁₆H₁₂) could be introduced simultaneously. The reaction *q* (reagent) during the ion-ion reaction was 0.4 for all experiments. The high pressure cell of the dual cell linear ion trap (where the ion-ion reactions are performed) was an experimental device with a front section length extended to 35 mm from 12.5 mm. This gives the device a ~3 fold higher reagent ion capacity. The MS³ Orbitrap spectra were acquired in full profile mode (no thresholding data compression) with Precursor Ion Target: 1e6, Precursor Max Injection Time: 800 ms, ETD Reagent Target: 2E6, IIPPT Reagent Target: 2E6. See Figures 1–3 and Table 1 for scheduling and definitions of the MS³ scans and overall LC-MS³ experiments.

Data Analysis: All LC-MS³ spectra in each time window for the three mAb subunits were averaged in Thermo Scientific™ QualBrowser™ software and single (averaged) spectra were generated. These averaged LC-MS raw files were processed through Thermo Scientific™ BioPharma Finder™ 3.0 software and the MS³ spectra were automatically *m/z* to mass "deconvoluted" using the Xtract algorithm. The signal-to-noise ratio (SNR) threshold for fragment peak picking was set to 3. The resulting neutral monoisotopic mass peak lists MS³ were exported to Excel where the mass lists (5 for ETD and 2 for CID and HCD) for each subunit where summed, then fed to ProSight Lite software (the mass tolerance of fragment ions: 10 ppm) to produce sequence coverage maps. To account for the possibility that Xtract may error in monoisotopic mass assignments by exactly ±1 Da were, in certain instances, summed into the aggregate mass lists and re-searched.

Figure 1. Preparation of NIST mAb into LC, Fd and Fc/2 Subunits Suitable for "Middle Down" LC/MS³ and LC/MS² Analyses.




Figure 2. LC/MS Selected Ion Chromatograms of Charge States of the NIST mAb GOF Isoform LC, Fd and Fc Subunits Selected for MS² and MS³ Analyses. The Liquid Chromatography Column was Deliberately Overloaded to Maximize Precursor Ion Abundances and to Extend the Time for Acquisition of MS² and MS³ Spectra.




Figure 3. Averaged MS Spectra Showing the Precursor Ion Charge States and *m/z* Selection Window Width of the NIST mAb GOF Isoform LC, Fd and Fc Subunits Selected for MS² and MS³ Analyses.




Table 1. Description Analyzer MS Selection and Dissociation/Reaction Settings for the "Middle-Down" LC-MS³ Experiments Targeting the Fc/2 (23+), LC (23+) and Fd (24+) Sub Units of the NIST mAb Standard (5 ea. MS/ETD/MS/IIPPT/MS, 2 ea. MS/CID/MS/IIPPT/MS and 2 ea. MS/HCD/MS/IIPPT/MS).

LC-MS ³ Experiment No.	NIST mAb Sub Unit Targeted	MS ³ Precursor <i>m/z</i> (Th)	MS ³ Precursor <i>m/z</i> Isolation Width (Th)	MS ³ Precursor Activation or Reaction Type	MS ³ Precursor Normalized Energy Setting (Reaction Time (ms))	Targeted MS ³ Precursor (MS ³ Product) <i>m/z</i> Selection Window	MS ³ Precursor Window (Th)	MS ³ Precursor Isolation Width (Th)	Targeted MS ³ <i>m/z</i> (Th)	MS ³ Precursor Activation or Reaction Type	MS ³ Precursor Reaction Time (ms)
1	Fc/2	1098.1	2	ETD	3	Between Low <i>m/z</i> and Intact MS ³ Precursor <i>m/z</i>	780-1080	300	930	IIPPT	6
	LC	1006.5	2	ETD	3		780-1080	300	850	IIPPT	6
	Fd	1071	2	ETD	3		780-1080	300	910	IIPPT	6
	Fc/2	1098.1	2	ETD	3	Between Precursor <i>m/z</i> and 1st Charge Reduced Intact Precursor <i>m/z</i>	1100.5-1140.6	40	1120.6	IIPPT	6
	LC	1006.5	2	ETD	3		1000-1060	40	1020	IIPPT	6
2	Fc/2	1098.1	2	ETD	3	Between 1st Charge Reduced Intact Precursor <i>m/z</i> and 2nd Charge Reduced Intact Precursor <i>m/z</i>	1073.5-1113.5	40	1093.5	IIPPT	6
	LC	1006.5	2	ETD	3		1100-1110	45	1120.5	IIPPT	6
	Fd	1071	2	ETD	3		1050-1100	45	1072.5	IIPPT	6
	Fc/2	1098.1	2	ETD	3	Between 2nd Charge Reduced Intact Precursor <i>m/z</i> and 3rd Charge Reduced Intact Precursor <i>m/z</i>	1200-1220	50	1200	IIPPT	6
	LC	1006.5	2	ETD	3		1100-1115	40	1130	IIPPT	6
3	Fc/2	1098.1	2	ETD	3	Between Low <i>m/z</i> and MS ³ Precursor <i>m/z</i>	1085-1088	300	938	IIPPT	6
	LC	1006.5	2	ETD	3		696.5-896.5	300	846.5	IIPPT	6
	Fd	1071	2	ETD	3		761-1061	300	911	IIPPT	6
	Fc/2	1098.1	2	ETD	3	Between MS ³ Precursor <i>m/z</i> and High <i>m/z</i>	1108-1408	300	1208	IIPPT	6
	LC	1006.5	2	ETD	3		1016.5-1316.5	300	1169.5	IIPPT	6
4	Fc/2	1098.1	2	ETD	3	Between Low <i>m/z</i> and MS ³ Precursor <i>m/z</i>	1085-1088	300	938	IIPPT	6
	LC	1006.5	2	ETD	3		696.5-896.5	300	846.5	IIPPT	6
	Fd	1071	2	ETD	3		761-1061	300	911	IIPPT	6
	Fc/2	1098.1	2	ETD	3	Between MS ³ Precursor <i>m/z</i> and High <i>m/z</i>	1108-1408	300	1208	IIPPT	6
	LC	1006.5	2	ETD	3		1016.5-1316.5	300	1169.5	IIPPT	6
5	Fc/2	1098.1	2	ETD	3	Between Low <i>m/z</i> and MS ³ Precursor <i>m/z</i>	1085-1088	300	938	IIPPT	6
	LC	1006.5	2	ETD	3		696.5-896.5	300	846.5	IIPPT	6
	Fd	1071	2	ETD	3		761-1061	300	911	IIPPT	6
	Fc/2	1098.1	2	ETD	3	Between MS ³ Precursor <i>m/z</i> and High <i>m/z</i>	1108-1408	300	1208	IIPPT	6
	LC	1006.5	2	ETD	3		1016.5-1316.5	300	1169.5	IIPPT	6
6	Fc/2	1098.1	2	HCD (Coll. Cell)	10	Between Low <i>m/z</i> and MS ³ Precursor <i>m/z</i>	1085-1088	300	938	IIPPT	6
	LC	1006.5	2	HCD (Coll. Cell)	10		696.5-896.5	300	846.5	IIPPT	6
	Fd	1071	2	HCD (Coll. Cell)	10		761-1061	300	911	IIPPT	6
	Fc/2	1098.1	2	HCD (Coll. Cell)	10	Between MS ³ Precursor <i>m/z</i> and High <i>m/z</i>	1108-1408	300	1208	IIPPT	6
	LC	1006.5	2	HCD (Coll. Cell)	10		1016.5-1316.5	300	1169.5	IIPPT	6
7	Fc/2	1098.1	2	HCD (Coll. Cell)	10	Between Low <i>m/z</i> and MS ³ Precursor <i>m/z</i>	1085-1088	300	938	IIPPT	6
	LC	1006.5	2	HCD (Coll. Cell)	10		696.5-896.5	300	846.5	IIPPT	6
	Fd	1071	2	HCD (Coll. Cell)	10		761-1061	300	911	IIPPT	6
	Fc/2	1098.1	2	HCD (Coll. Cell)	10	Between MS ³ Precursor <i>m/z</i> and High <i>m/z</i>	1108-1408	300	1208	IIPPT	6
	LC	1006.5	2	HCD (Coll. Cell)	10		1016.5-1316.5	300	1169.5	IIPPT	6
8	Fc/2	1098.1	2	HCD (Coll. Cell)	10	Between Low <i>m/z</i> and MS ³ Precursor <i>m/z</i>	1085-1088	300	938	IIPPT	6
	LC	1006.5	2	HCD (Coll. Cell)	10		696.5-896.5	300	846.5	IIPPT	6
	Fd	1071	2	HCD (Coll. Cell)	10		761-1061	300	911	IIPPT	6
	Fc/2	1098.1	2	HCD (Coll. Cell)	10	Between MS ³ Precursor <i>m/z</i> and High <i>m/z</i>	1108-1408	300	1208	IIPPT	6
	LC	1006.5	2	HCD (Coll. Cell)	10		1016.5-1316.5	300	1169.5	IIPPT	6
9	Fc/2	1098.1	2	HCD (Coll. Cell)	10	Between Low <i>m/z</i> and MS ³ Precursor <i>m/z</i>	1085-1088	300	938	IIPPT	6
	LC	1006.5	2	HCD (Coll. Cell)	10		696.5-896.5	300	846.5	IIPPT	6
	Fd	1071	2	HCD (Coll. Cell)	10		761-1061	300	911	IIPPT	6
	Fc/2	1098.1	2	HCD (Coll. Cell)	10	Between MS ³ Precursor <i>m/z</i> and High <i>m/z</i>	1108-1408	300	1208	IIPPT	6
	LC	1006.5	2	HCD (Coll. Cell)	10		1016.5-1316.5	300	1169.5	IIPPT	6

RESULTS

The aggregated results from experiments 1-5 (MS/ETD/MS/IPT) were by far the most promising with many long sequence ions observed (see Figures 5-8) for all the NIST mAb subunits. The results from experiments 6-7 (MS/CID/MS/IPT/MS) and Experiments 8-9 (MS/HCD/MS/IPT/MS) were less promising. The pair of CID experiments provided the following aggregate sequence coverages: 19% (Fc/2), 31% (LC) and 31% (Fd/2). The pair of HCD experiments provided the following aggregate sequence coverages: 13% (Fc/2), 17% (LC) and 17% (Fd/2). Very few sequence ions having lengths more than 50% of the subunits were observed.

Figure 5. Sequence Coverage Map NIST mAb Fc/2 Subunit from the Product Ion Mass List From Experiments 1-5 (LC-MS/ETD/MS/IPT/MS Analyses).




Figure 6. Sequence Coverage Map NIST mAb Fd Subunit from the Collective Product Ion Mass List From Experiments 1-5 (LC-MS/ETD/MS/IPT/MS Analyses).



Figure 4. Example MS/ETD/MS Spectrum of the Targeted Fd Precursor Ion. A) Full Vertical Scale: Illustrating the 5 ea. MS³ IIPPT Precursor Ion *m/z* Selection Windows for the 5 Targeted MS/ETD/MS/IIPPT/MS Analyses. B) 20×Expanded Vertical Scale: Illustrating the *m/z* Distribution of ETD and ETnCD Product Ions.





Figure 7. Sequence Coverage Map NIST mAb LC Subunit from the Collective Product Ion Mass List From Experiments 1-5 (LC-MS/ETD/MS/IPT/MS Analyses).



Figure 8. Sequence Coverage Map NIST mAb LC Subunit from the Collective Product Ion Mass List From Experiments 1-5 (LC-MS/ETD/MS/IPT/MS Analyses) Where Additional Masses Exactly +1 Da and -1 Da from the Masses Obtained by *m/z* to Mass "Deconvolution" of the Averaged Spectra Were Included in the Search.



Figure 9. Averaged MS/ETD/MS/IIPPT Spectrum of Targeted Fd Precursor Ion Obtained in LC-MS³ Experiment #4 A) With *m/z* Range Chosen to Entire Range of *m/z* of Product Ions. B) 50 Th Window (Highlighted in Red) With Mostly Lower Charge State (5+ to 7+) Product Ions. C) 50 Th Window (Highlighted in Blue) With Mostly Higher Charge State (8+ to 10+) Product Ions.



CONCLUSIONS

This preliminary study demonstrates that IIPPT of selected *m/z* ranges of ETD product ions greatly enhances the observation of large sequence ions for all of the NIST mAb subunits. A similar increase of large sequence ions for the IIPPT subsequent to CID and HCD was not demonstrated. We believe this was due to the relatively large (300 Th) *m/z* windows of CID and HCD product ions selected for IIPPT. We anticipate that observation of large UVPD product ions of mAb sub units will also be enhanced by a subsequent IIPPT of selected *m/z* windows as was done in the ETD experiments.

REFERENCES

- Houel, S.; Fornelli, L.; Srzentic, K.; et al. Middle-down analyses of monoclonal antibody on Orbitrap Fusion Lumos Tribid. *Proceedings of the 16th Annual World Congress of the Human Proteome Organization (HUPO) in Dublin, Ireland, 2017*.
- Anderson, L. C.; English, M.A.; Wang, W. H.; et al. Protein derivatization and sequential ion/ion reactions to enhance sequence coverage produced by electron transfer dissociation mass spectrometry. *JMS* 2015, 377, 717–724.
- Mullen, C.; Schwartz, J. C.; Syka, J. E. P.; et al. *Improved top-down sequence coverage on a Tribid mass spectrometer by ion-ion proton transfer (IIPPT) reactions subsequent to ETD and UVPD*. Proc. 66th ASMS Conf. Mass Spectrometry and Allied Topics, San Diego, VA [TP741].
- Sanders J, Holden D, et al. Reducing Congestion in Ultraviolet Photodissociation spectra of large protein ions using Proton Transfer Reactions. Proc. 66th ASMS Conf. Mass Spectrometry and Allied Topics, San Diego, VA [TOH pm 02:30].

TRADEMARKS/LICENSING

© 2018 Thermo Fisher Scientific Inc. All rights reserved. ProSight is a trademark of Proteinaceous, Inc. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.