



ThermoFisher
S C I E N T I F I C

Introducing the LC-MS/MS Portfolio – TSQ Quantis and Altis

Osama Abu-Nimreh
CMD Sales Support Specialist
MECEC , Dubai

The world leader in serving science

Introduction to TSQ Altis and TSQ Quantis

Performance: Sensitivity, Selectivity (H-SRM)



	Thermo Scientific™ TSQ Altis™ <i>High-end</i>	Thermo Scientific™ TSQ Quantis™ <i>Mid-tier</i>
Mass Range	5-2000	5-3000
SRM/sec	600	600
Selectivity (H-SRM)	0.2 Da FWHM	0.4 Da FWHM
Sensitivity (HESI Reserpine 1 pg)	500,000:1	150,000:1
Targeted Market	Omics, Research, Pharma/Biopharma, Clinical Research and Forensic Toxicology	Environmental and Food Safety, Clinical Research, and Forensic Toxicology

Robustness, Reproducibility, Speed, Ease-of-Use, Flexibility

TSQ Altis: Sensitivity with Robustness, No Compromises

AIM+
TECHNOLOGY

Active Ion Management Plus (AIM+) - The next step in precision design delivers the ultimate in ion management, inception to detection, from the OptaMax™ ion source housing to the enhanced electron multiplier. Incorporates segmented quadrupoles with hyperbolic surface and enhanced RF Electronics to further optimize ion management precision, reliability, speed, and reproducibility.

Ion beam guide with neutral blocker
Reduces chemical background

High capacity ion transfer tube (HCTT)
Increases ion flux

Electrodynamic ion funnel (EDIF)
Increases ion flux

OptaMax™ NG
APCI ready

Segmented Quadrupoles
with hyperbolic surface for enhanced performance with both SRM and H-SRM (0.2 FWHM)

Active collision cell with axial DC field
facilitates more SRMs/sec

Enhanced dual-mode electron multiplier detector
Ensures excellent linearity and dynamic range

NEW!

NEW!

NEW!

TSQ Quantis: Unprecedented Robustness, Day After Day



Active Ion Management Plus (AIM+) - The next step in precision design delivers the ultimate in ion management, inception to detection, from the OptaMax™ ion source housing to the enhanced electron multiplier. Incorporates segmented quadrupoles with hyperbolic surfaces and enhanced RF Electronics to further optimize ion management precision, reliability, speed, and reproducibility.

Enhanced dual-mode electron multiplier detector
ensures excellent linearity and dynamic range



Stacked ring ion guide (SRIG)
Increases ion flux

Segmented Quadrupoles
with hyperbolic surfaces for enhanced performance with both SRM and H-SRM (0.4 FWHM)



OptaMax™ NG
APCI ready



Ion beam guide with neutral blocker
Reduces chemical background

Active collision cell with axial DC field
facilitates more SRMs/sec

OptaMax NG Source Housing

Benefits: Reliable and consistent performance with improved usability!

• Re-designed APCI discharge assembly

- Built-in to every source (separate APCI sprayer required for APCI mode)
- Re-designed on/off switch (to improve usability)

• Re-designed HESI Sprayer

- Needle adjustment is no longer possible during acquisition (locked position)
- Tool available to help the user to correctly set needle protrusion

• Usability and Consistency

- Vertical adjustment moved to the side for easier access
- New drain insert with improved latching and locating pin to prevent rotation
- Improved sprayer alignment and stability
- New finer threads on HESI and APCI sprayers to make installation easier



Segmented Quadrupoles


Benefits: Increased Sensitivity (more significant at higher mass range)
Flat tuning for consistent and robust performance

- The use of RF only pre-filters (segments) between the entrance lens and the quadrupole minimizes the effects of fringe fields, leading to improved transmission (and therefore sensitivity) at unit and higher resolution.
- With the RF only pre-filter, the tuning of several lenses is flat across mass range allowing the voltage to be set and not tuned. This helps reducing the complexity of the tune and making the systems more consistent.




Benefits: Increased electron multiplier lifetime. Increased Uptime!

- Increased number of dynodes (21) for extended lifetime.
- Improved electron multiplier calibration routine.
- Excellent linearity and dynamic range across the mass range.
- Reduced number of service visits leading to more uptime.



Benefits: More compounds in the same run or longer dwells on existing method

- New main RF/DC electronics
- Analyze more compounds in the same time window or better Quantitation results with better ion statistics (more scans across your chromatographic peak)
- Up to **600 SRM/sec**



What Makes the new Triple Quads Robust?

Robustness

The ability to perform at the expected level (LOD, LOQ, MRL) under adverse conditions (complex samples, limited sample preparation) for the desired period of time without maintenance.

OptaMax NG



Re-designed source housing

NEW!

Reliable and consistent performance with improved usability!

ITT/Sweep Cone



Ion transfer tube and sweep cone

No need to break vacuum for basic maintenance!

Ion Beam Guide





Ion beam guide with neutral blocker

Blocks neutral keeping the ion path clean while reducing noise!

What Makes the new Triple Quads Consistent?

Consistency


The ability to report equivalent results every time independently of the system or the user (day-to-day and system-to-system)

OptaMax NG	Optics	Segmented Quads
		
Re-designed Source Housing	Improved Optics Alignment	Segmented Quads with hyperbolic surfaces
		
<i>Reliable and consistent performance with improved usability!</i>	<i>More consistent performance from system-to-system!</i>	<i>Flat tuning across the mass range!</i>

What Makes the new Triple Quads Reliable?


Reliability

Increased time between service calls (mean time between failure- MTBF)
Increased uptime







Demonstration of Robustness – Food Safety

Atrazine QC monitored in leek for more than 400 injections with 4.5% RSD . Red lines represent $\pm 20\%$ response at $10 \mu\text{g}/\text{Kg}$. Yellow lines show the time the system was placed in standby mode for 12h to demonstrate consistent performance after standby period





Demonstration of Robustness – Clinical Research

Immunosuppressants (Cyclosporin A, Everolimu, Sirolimus and Tacrolimus) monitored in crashed whole blood for more than 1500 injections. Red lines represent $\pm 20\%$ of calculated amounts (ng/mL)



Sensitivity with Reproducibility



Active Ion Management Plus (AIM+)


The next step in precision design delivers the ultimate in ion management, inception to detection, from the OptaMax NG source housing to the enhanced electron multiplier. Incorporates segmented quadrupoles and enhanced RF Electronics to further optimize ion management precision, reliability, speed, and reproducibility.

TSQ Altis: Collaboration with Beatson Institute for Cancer Research


LPA 17:0 @ 0.01 ng/mL

TSQ Quantiva

TSQ Altis




2x Greater!



RSD 3.4% (n=3)

RSD 0.9% (n=3)

Glutamic Acid




2.5 x Greater!

LPA at 0.01 ng/mL with RSD below 1%! 2 x greater response on TSQ Altis! Confirming Ion detected on TSQ Altis! A 2.5 fold response increase was observed when transferring the metabolomics assay from TSQ Quantiva to TSQ Altis.

TSQ Altis: Quantitation of Therapeutic Drugs in Plasma

Compound Name	TSQ Quantiva LOQ (pg/mL)	TSQ Altis LOQ (pg/mL)
Desomorphine	5	5
Desmethyldoxepin	10	2.5
Flecainide	2.5	1
Midazolam	5	2.5
Imipramine	10	2.5
Amitriptyline	10	2.5
Fluoxetine	5	5
Diazepam	5	2.5

Quantitation of Desmethyldoxepin in plasma
2.5 pg/mL



~ 2.5 X average sensitivity improvement over TSQ Quantiva

TSQ Altis: Confident Quantitation of Challenging Analytes in Environmental Matrices


IC-MS/MS solutions for environmental analysis Quantitation of **Glyphosate** at 5 ng/L, 3% RSD

Experimental Details

ICS 5000+
Flow rate: 0.3 mL/min
Eluent Source: Eluent Generator
Mobile Phase: KOH


MS: TSQ Altis

Software: TraceFinder Software 4.1





TSQ Quantis: Robust Quantitation of ETG and ETS

Excellent Precision for the most challenging quantitative assays




50 pg on column for Ethyl Sulfate in Urine



Experimental Details

UHPLC : Vanquish Flex Binary
Flow rate: 0.4mL/min
Solvent A: 0.1%FA in H₂O
Solvent B: 0.1%FA in MeOH


TSQ Quantis: Sensitivity Improvement for Pesticide Residue Analysis



3 x average sensitivity improvement across mass range!

TSQ Quantis: Sensitivity in Regulatory Methodology (EPA 544)

TSQ Endura





In 2014, half-a-million people were without drinking water in the state of Ohio due to high toxin level!



Algae blooms in Western Lake Erie

EPA 544 released in February 2015

TSQ Quantis



~ 4 x average sensitivity improvement over TSQ Endura


Speed

- **Active Collision Cell with axial DC field**

- 90° cell design for noise reduction


- **New main RF/DC electronics**

- Analyze more compounds in the same time window or better Quantitation results with better ion statistics (more scans across your chromatographic peak)
- Up to **600 SRM/sec**




Robust, Reliable, Fast Quantitation Workflows

Excellent Quantitative Performance at Lower Dwell Times!




SRM Visualization from Instrument Control Software 3.0 displaying the number of transitions per unit time



~ 160 Transitions Monitored Simultaneously with
Polarity Switching. Excellent Reproducibility (% RSD 2.3) below the MRL

Performance at Extreme Speeds – 500 SRMs/sec vs 600 SRMs/sec

Equivalent Quantitative performance between 500 and 600 SRMs/sec!
Good Performance at extremely low dwell times!



SRMs/Sec	Total Number of Transitions	Dwell Time (mSec)
500	1075	0.769
600	1291	0.437

Atrazine Concentration (ppb)	500 SRM/Second		600 SRM/Second	
	Average Area	%CV	Average Area	%CV
1	21682	9	18090	9
10	475465	4	369612	5
100	4326117	1	4296555	1

Selectivity with High Resolution SRM (H-SRM)



*New segmented quadrupoles
with hyperbolic surfaces*




TSQ Altis – 0.2 Da FWHM

TSQ Quantis – 0.4 Da FWHM

*Reduced noise when
analyzing complex matrices –
better S/N – lower LOD/LOQ*

Superior Sensitivity with H-SRM (0.2 Da FWHM) – GPSVFPLAPSSK

GPSVFPLAPSSK - Peptide from monoclonal antibody hinge region



Experimental Details

LC : Ultimate NCS-3500RS
 Flow rate: 25 μ L/min
 Solvent A: 2% ACN in H₂O w/0.1% FA
 Solvent B: 5% H₂O in ACN w/0.1% FA

25 ng GPSVFPLAPSSK and IS

- **Market Specific Software:**
Thermo Scientific™ TraceFinder™ software

- SRM databases for Environmental, Food Safety, Clinical Research and Forensic Toxicology applications: easy instrument method setup
- Compatible with GC and LC high resolution systems. One software package for different platforms
- Customized reports
- Flexibility on Data Review – compound or sample centric modes for easy data review. Customized flagging system for easy review of important criteria.



TraceFinder

- **Instrument Control Software: Tune 3.0**

- **User**

- Automated Compound Optimization
(via infusion or Autosampler)
- Dynamic Retention Time (dRT)
- Method Conversion Tool (Agilent/Sciex)
- Improved Auto tune and Calibration Routines
(Faster with improved reporting)
- Favorites – Easy to setup source settings for Calibration



- **Service**

- New/Smarter Diagnostics – better serviceability of the system – improved reporting for service.

Basic Maintenance



5 min operation!
No need to break vacuum!

OptaMax NG




**Re-designed
source housing**

NEW!

*New tool available to ensure
consistent positioning of
HESI needle*

ITT/Sweep Cone



**Ion transfer tube
and sweep cone**

*No need to break vaccum for
basic maintenance
5 minute Operation!*

Flexibility

- **LC Options**

- **High Sample Throughput**

- Multi channel – Transcend and Prelude
- CTC PAL 3 – reduced time between injections, barcode reading

- **High Performance LC**

- UHPLC – Vanquish Horizon and Vanquish Flex

- **MS Ion Sources**

- H-ESI
- APCI
- OptaMax – APCI ready (only need to change the sprayer – 2 min operation)
- APPI
- Easy spray source and flex for nanoLC applications

- **Software Options**

- TraceFinder Software


- Method development to automated report generation – offering comprehensive solution for every target application

- Thermo Scientific™ Chromeleon™ Software

- Data acquisition and processing

- AB Sciex Skyline Software

- Compatibility with TraceFinder



Critical Resources Worth Looking At

thermoscientific

APPLICATION BRIEF 6495

Quantification of cyanotoxins in drinking water according to EPA 544 guidelines

Authors Goal

thermoscientific

APPLICATION BRIEF 6497

Robustness, reproducibility, reliability with best-in-class sensitivity; Increased confidence in targeted quantitation of pesticides in food matrices

Pesticide Explorer Collection – Standard Quantitation

Authors Goal
Katerina Bousova,* Michal Godula,*

thermoscientific

APPLICATION BRIEF 6492*

Reduced injection volume applied to the quantification of cylindrospermopsin and anatoxin-a in drinking water according to EPA Method 545

Authors Goal
Al Higgins,* Andy Eaton,*
Nelson Wjanatha,* Claudia Martinez*
*Eurofins Eaton Analytical.

thermoscientific

TECHNICAL NOTE 64969

Testing robustness: Immunosuppress blood with a TSQ Quantis MS for cli

Authors Goal
Kristine Van Natta, Nelson Wijesinghe,
Claudia Martinez
Thermo Fisher Scientific,
San Jose, CA

Introduction



Better Data
Reliable, Reproducible, Robust
Increased Productivity



thermoscientific

TECHNICAL NOTE 64970

Analysis of ethyl glucuronide and ethyl sulfate in urine with a TSQ Quantis MS

thermoscientific

Pesticides in Food: Targeted Quantitation Solution Guide
Your lab, your challenges—our solutions to help you address them

thermoscientific

Biopharma: Targeted Quantitation Solution Guide
Your lab, your challenges—our solutions to help you address them



thermoscientific

Confident Quantitation
Robust and reliable performance day after day

Thermo Scientific TSQ Quantis triple-stage quadrupole mass spectrometer

ThermoFisher SCIENTIFIC

ThermoFisher SCIENTIFIC Confident Quantitation

Confident Quantitation
Any compound, any matrix, any user.

Learn how to make a difference in your lab today

Confident Quantitation
Every molecule, every matrix, every user.

