

Trifluoroacetic Acid

Product Specification

Trifluoroacetic acid (TFA) is a very versatile reagent. Some of the more popular uses for TFA are as a silyl catalyst when derivatizing carbohydrates, as a reagent for purifying large peptides, as an ion pair reagent, and as an LC mobile phase additive.

Features/Benefits

Silyl Catalyst Addition of a small amount of acidic catalyst usually increases the rate or degree of silylation. In acid catalysis, protonation of the silyl donor weakens the Si-X bond (X is the leaving group). TFA derivatives are stable and volatile. Use of TFA in combination with HMDS avoids the formation of ammonium chloride.

Ion Pair Reagent Use of TFA as an ion pair reagent in reversed phase ion-pair chromatography on a nonionic column can maximize polarity differences between proteins or peptides in a sample mixture and improve chromatographic separation. (To avoid protein denaturation, use 0.1% or less TFA.)

Mobile Phase Additive In reversed phase HPLC, inclusion of TFA in the aqueous mobile phase provides an ion with detectable properties and with affinity for the solid phase. Analyte injection into the system gives rise to equilibrium disturbances and influences the distribution of the detectable ionic component.

Typical Procedures

These procedures are intended to be guidelines and may be adapted as necessary to meet the needs of a specific application. Always take proper safety precautions when using a silylating reagent – consult MSDS for specific handling information. TFA is extremely sensitive to moisture and should be handled under dry conditions.

Carbohydrates in Syrups

1. Place 60-70mg 80% soluble syrup in a vial and dissolve in 1 mL pyridine.
2. Add 900µL HMDS, then 100µL TFA. Shake 30 seconds, then let stand 15 minutes, with occasional shaking. Inject aliquot onto GC column for analysis.

The reagent and the reactivity of the compounds to be derivatized play equally important parts. Derivatization times will vary widely, depending upon the compound being analyzed. To determine when derivatization is complete, analyze aliquots of the sample at selected time intervals until no further increase in product peak(s) is observed.


Aflatoxins

(AOAC Method 990.33, Official Methods of Analysis, 16th Ed. 1995.)

1. Place extracted sample containing aflatoxins in screw cap vial. Evaporate to dryness using *clean* nitrogen.
2. Add 200µL hexane to re-dissolve aflatoxins.
3. Add 50µL TFA, cap, and vortex 30 seconds. Let stand 5 minutes.
4. Add 2.0mL deionized water/acetonitrile (9:1). Vortex 30 seconds, then allow layers to separate.
5. Remove aqueous (bottom) layer containing aflatoxins. Filter through 0.45µm syringe-tip filter and inject aliquot onto LC column.

Properties

Structure:



CAS Number:

76-05-1

Molecular Formula:

CF₃COOH

Formula Weight: 114.02

bp: 72.4°C

d: 1.480

Appearance:

clear, colorless liquid
moisture sensitive

796-0342

Toxicity – Hazards – Storage – Stability

TFA is a moisture-sensitive liquid. It may irritate eyes, skin, and/or the respiratory system. Store in a brown bottle or amber ampul at room temperature, in a dry, well ventilated area away from ignition sources. Use only in a well ventilated area and keep away from ignition sources.

Properly stored, this reagent is stable indefinitely. Recommended storage conditions for the unopened product are stated on the label. Moisture will decompose both the reagent and its derivatives. To exclude moisture, Supelco packages this product under nitrogen. If you store an opened container or transfer the contents to another container for later reuse, add desiccant. **Before reuse, validate that your storage conditions adequately protected the reagent.**

Ordering Information

Description	Cat. No.
TFA	
10 x 1mL	33077
25mL	33075
100mL	33076

Micoreaction Vessels with Hole Caps and Septa

1mL, pk. of 12	33293
3mL, pk. of 12	33297
5mL, pk. of 12	33299

Books

<i>Handbook of Derivatives for Chromatography</i> K. Blau and J. Halket	Z246220
<i>Handbook of Analytical Derivatization Reactions</i> D.R. Knapp	23561

Additional Reading

- K. Blau and J. Halket *Handbook of Derivatives for Chromatography* (2nd ed.) John Wiley & Sons, New York, 1993.
- D.R. Knapp *Handbook of Analytical Derivatization Reactions* John Wiley & Sons, New York, 1979.
- M. Morvai, P. I. Monlnar *Simultaneous Gas Chromatographic Quantitation of Sugars and Acids in Citrus Fruits, Pears, Bananas, Grapes, Apples and Tomatoes* *Chromatographia*, 34 (9-10): 502-504 (1992).
- G.W. Chapman, R.J. Horvat *Determination of Non-volatile (Organic) Acids and Sugars from Fruits and Sweet Potato Extracts by Capillary GLC and GLC-MS* *J. Agric. Food Chem.*, 37 (4): 947-950 (1989).
- P. Englmaier *High Resolution GLC of Carbohydrates as Their Dithioacetal-Trimethylsilylates and Trifluoroacetates* *J. High Res. Chromatogr.*, 13 (2): 121-125 (1990).

Contact our Technical Service Department (phone 800-359-3041 or 814-359-3041, FAX 800-359-3044 or 814-359-5468) for expert answers to your questions.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA · Sigma-Aldrich de Argentina, S.A. · Buenos Aires 1119 **AUSTRALIA** · Sigma-Aldrich Pty. Ltd. · Castle Hill NSW 2154 **AUSTRIA** · Sigma-Aldrich Handels GmbH · A-1110 Wien
BELGIUM · Sigma-Aldrich N.V./S.A. · B-2880 Bornem **BRAZIL** · Sigma-Aldrich Quimica Brasil Ltda. · 01239-010 São Paulo, SP **CANADA** · Sigma-Aldrich Canada, Ltd. · 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC · Sigma-Aldrich s.r.o. · 186 00 Praha 8 **DENMARK** · Sigma-Aldrich Denmark A/S · DK-2665 Vallensbaek Strand **FINLAND** · Sigma-Aldrich Finland/YA-Kemia Oy · FIN-00700 Helsinki
FRANCE · Sigma-Aldrich Chimie · 38297 Saint-Quentin-Fallavier Cedex **GERMANY** · Sigma-Aldrich Chemie GmbH · D-82041 Deisenhofen **GREECE** · Sigma-Aldrich (o.m.) Ltd. · Ilioupoli 16346, Athens
HUNGARY · Sigma-Aldrich Kft. · H-1067 Budapest **INDIA** · Sigma-Aldrich Co. · Bangalore 560 048 **IRELAND** · Sigma-Aldrich Ireland Ltd. · Dublin 24 **ISRAEL** · Sigma Israel Chemicals Ltd. · Rehovot 76100
ITALY · Sigma-Aldrich s.r.l. · 20151 Milano **JAPAN** · Sigma-Aldrich Japan K.K. · Chuo-ku, Tokyo 103 **KOREA** · Sigma-Aldrich Korea · Seoul **MALAYSIA** · Sigma-Aldrich (M) Sdn. Bhd. · 58200 Kuala Lumpur
MEXICO · Sigma-Aldrich Quimica S.A. de C.V. · 50200 Toluca **NETHERLANDS** · Sigma-Aldrich Chemie BV · 3330 AA Zwijndrecht **NORWAY** · Sigma-Aldrich Norway · Torshov · N-0401 Oslo
POLAND · Sigma-Aldrich Sp. z o.o. · 61-663 Poznań **PORTUGAL** · Sigma-Aldrich Quimica, S.A. · Sintra 2710 **RUSSIA** · Sigma-Aldrich Russia · Moscow 103062 **SINGAPORE** · Sigma-Aldrich Pte. Ltd.
SOUTH AFRICA · Sigma-Aldrich (pty) Ltd. · Jet Park 1459 **SPAIN** · Sigma-Aldrich Quimica, S.A. · 28100 Alcobendas, Madrid **SWEDEN** · Sigma-Aldrich Sweden AB · 135 70 Stockholm
SWITZERLAND · Supelco · CH-9471 Buchs **UNITED KINGDOM** · Sigma-Aldrich Company Ltd. · Poole, Dorset BH12 4QH
UNITED STATES · Supelco · Supelco Park · Bellefonte, PA 16823-0048 · Phone 800-247-6628 or 814-359-3441 · Fax 800-447-3044 or 814-359-3044 · email: supelco@sial.com

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

H

BFD