



Tips, Tricks, and Troubleshooting for Successful Headspace Analysis

Mark Sinnott
GC Application Scientist
September 17, 2024



What Is Headspace?



$$K = \frac{C_S}{C_G}$$

$$\beta = \frac{V_G}{V_S}$$

$$C_S = \frac{W_S}{V_S}$$

$$C_G = \frac{W_G}{V_G}$$

$$W_S + W_G = W_0$$

$$C_0 \times V_S = (C_G \times V_G) + (C_S \times V_S) = (C_G \times V_G) + (K \times C_G \times V_S) = C_G \times [K \times V_S + V_G]$$

$$C_0 = C_G \left[\frac{K \times V_S}{V_S} + \frac{V_G}{V_S} \right] = C_G (K + \beta)$$

$$C_G = \frac{C_0}{K + \beta}$$

$$A \propto C_G = \frac{C_0}{K + \beta}$$

HS Advanced Operation user guide:

<https://www.agilent.com/cs/library/usermanuals/public/G4556-90016.pdf>

Why Headspace?

- Offers clean injections into GC systems
 - Less GC maintenance – only the volatile vapors are injected into the system
 - Non-volatiles/semi-volatiles/acids etc. stay in the liquid phase of the HS vial
- Less sample preparation
- Ideal for analysis of volatile analytes in dirty matrices that can't be directly injected into the GC.

*Not suitable for some applications



Some Math to Make it Fun

$$CoVs = CgVg + CsVs$$

Partition coefficient: $K = \frac{Cs}{Cg}$

Phase ratio: $\beta = \frac{Vg}{Vs}$

$$Cg = \frac{Co}{(\beta + K)}$$



What Should We Focus On?

Partition coefficient: $K = \frac{C_s}{C_g}$

The smaller the “ K ”, the greater the concentration of the analyte in the gas phase.

Like dissolves like. The greater the solubility or affinity that an analyte has for the matrix, the larger the K .

What drives K ?

Temperature:

- Higher temperatures drive K down = better sensitivity

Reduce Solubility = Reduce “ K ” = Improved response

- Add salt
- Add another solvent to the matrix

What Parameters Drive Success?

- Incubation temperature
 - Typically, 10-20 °C below the solvent BP (need to avoid boiling the solvent!)
- Incubation time
 - Determined experimentally
 - Need to ensure equilibrium is reached
 - Once equilibrium is reached, no need to incubate longer
- Shaking
 - Shortens incubation time
- Efficient transfer of the sample from the vial to the column
- Use of salts



Things to Consider

- You will need to have at least 5 mL of headspace in the vial.
 - 15 ml maximum for 20 mL vials / 10 mL maximum for 10 mL vials
- Keep the incubation temperature 10 to 20 °C below the BP of the solvent/matrix.
- Long incubation times 'generally' only delay the first sample.
 - The HS can over-lap samples
- Higher split ratios help get the sample onto the column more efficiently; this results in sharper peaks.
 - Lower splits are 'OK' with larger id columns. Higher volumetric flow transfers sample faster.
- Shake, but try to keep the sample from splashing the vial septum.
 - Sample can get into the sample probe and contaminate the loop
- Think about the temperature limitation of vial septa
 - Be considerate of sample/analyte degradation

Headspace Parameters

Temperatures

- Oven
- Sample loop
- Transfer line
- Transfer line interface

Times

- Vial equilibration
- Injection duration
- GC cycle time

Vial and loop

- Vial size
- Shake vials while in oven
- Vial fill mode
- Loop fill mode

Starting Parameters

Temperatures

- Oven **20 °C below the BP of the matrix**
- Sample loop **Same temperature as the oven or higher (+10 c)**
- Transfer line **Hot enough not to have anything condense (+ 10 c)**
- Transfer line interface **Same as the inlet**

Times


- Vial equilibration **10 minutes, but use method development**
- Injection duration **0.5 minutes**
- GC cycle time **Run time + cool down to ready**

Vial and Loop


- Vial size **20 mL**
- Shake vials while in oven **3 (low)**
- Vial fill mode **Default 15 psi**
- Loop fill mode **Default**

Incubation Temperature Increase

20 minutes
K decreases with T
Not equal for all analytes



Incubation Time




Change in Vial Size


4 mL sample, changing β

10 mL vial $\beta = 1.5$

20 mL vial $\beta = 4$



Change in Sample Volume in a 20 mL Vial




What Else Can Affect Signal?

- Loop size
- Split ratio
- Liner type




Change in Loop Size


40:1 split (64 mL/min)



Change in Split Ratio



Liner Size and Type



Use of Salts

Decreases the solubility of polar analytes in aqueous samples

Decreases K favoring the gas (headspace) phase

Potassium carbonate (K_2CO_3)

Ammonium chloride (NH_4Cl)

Ammonium sulfate ($(NH_4)_2SO_4$)

Sodium chloride ($NaCl$)

Sodium citrate ($Na_3C_6H_5O_7$)


Sodium sulfate (Na_2SO_4)

Use high quality, low impurity salts

Be sure to check blanks with salt for contaminants


Effect of salt

20 mL vial
80 °C oven temperature
20 minute incubation





Can I Inject Multiple Times?

Only for developing a chromatographic method




Method Development Tools


 Edit Method Parameters




Temperatures




Times




Vial and Loop




Carrier



Advanced Functions



Sequence Actions



Method Development


Method Development

[Manual](#)


Would you like to increment a method setting over subsequent runs?

None


[Assisted](#)



Create method based on a specific application




Convert an existing valve and loop Headspace method





Convert an existing pressure transfer Headspace method


Method Development Tool


Edit Method Parameters



Temperatures



Times


Vial and Loop


Carrier


Advanced Functions


Sequence Actions


Method Development

Method Development

Manual

Would you like to increment a method setting over subsequent runs?

Temperature

None

Temperature

Vial Equilibration

Shaking


Temperature increment:

10 °C


Maximum oven temperature:

80 °C

Assisted



Create method based on a specific application



Convert an existing valve and loop Headspace method

Choose what you want to increment

Choose your setpoints

Types of Vials



5182-0837




5188-2753



5067-0226



5182-0838



Consumables



Good for SPME



Safety cap
Tears at 45 psi



Max temp 125 °C
Butyl/PTFE



Max temp 180 °C
Silicone/PTFE

High Performance Septa

Maximum temperature is 300 °C

Reduce siloxane interferences at high temperatures



High-power crimpers are required for steel caps

Publication number: 5990-9385EN

High-Power Crimper



Crimper with 20 mm jaw set
5190-4067
Required when using steel camps

Standard Crimpers



5040-4669



5190-3189

How Tight is Right?



Good crimp



Too tight




Too loose



Troubleshooting

- No peaks/low peak response
- Retention times not repeatable
- Peak areas not repeatable/poor RSD
- Contamination or carryover




Troubleshooting – General

Always consider the entire system:

Follow the flow path...

- Carrier gas →
- Sample/sample preparation →
- Sample introduction/inlet →
- Detection →
- Data handling and...
- All points in between



Troubleshooting – General


Always check the logs:

- Does the problem correlate with any recent changes or events?
 - Maintenance
 - Power Failure
- Check HS sequence log
- Check maintenance logs
- Half-split the system
- Check the GC
 - Check error messages/logs
 - If possible, mount an ALS to see if the issue persists
- Run something simpler:
 - System blanks/non-injections
 - Solvent blanks
 - Standards



<https://www.agilent.com/cs/library/usermanuals/public/G4556-90018.pdf>

The Logbook Is in the Instrument Control Screen





Troubleshooting – Missing Peaks/Low Peak Response


- Check for proper vial crimp
- Check delivery gas pressures/tanks
 - Delivery pressure should be 20 PSI higher than the highest pressurization setting
- Enable dynamic leak checking
- Check for leaks in HS
 - Transfer line leaks/broken?
 - Probe, valve lines, and fittings





Vial Leaks



Temperatures



Times


Vial and Loop


Carrier







Advanced Functions


Sequence Actions


Method Development

Sequence Actions

What should the sequence do if it encounters the following:

Vial Missing	Wrong Vial Size	Leak Detected	System Not Ready	Wrong Cooling Plate Temp
				
Skip	Continue	Continue	Abort	Continue

The system always logs detected issues and the action taken.

Troubleshooting – Retention Time Instability

- Check the GC
 - Inlet leaks-septum, liner
 - Gas supply pressure
 - Is the GC given enough time to stabilize between runs?
- Confirm that the column installed matches the one in the method
- Check for leaks in HS
 - Transfer line leaks/broken?
 - Probe, valve lines, and fittings



Troubleshooting – Peak Area Instability

- Check for inconsistencies in sample preparation
- Check for carry-over
- Check for proper vial crimp
- Run standard replicates and calculate RSD
- Enable dynamic leak checking
- Sample not at equilibrium
- Check for leaks in HS
 - Transfer line leaks/broken?
 - Probe, valve lines, and fittings



Troubleshooting – Contamination or Carryover

- Check gas supply/traps
- Check the GC
 - Inlet consumables: septum, liner gold seal, vent trap
 - Check error messages/logs
- Run several instrument blanks
- Run several solvent blanks
 - If carryover does not decay, it probably condensed in the flow path (loop, probe, transfer line)
- Enable dynamic leak checking
- Check HS parameters
 - Loop purge flow/time
 - Enable standby flow
 - Check sample path temps proper: loop/valve, transfer line



Other Common Issues

Carryover/contamination

- Too much sample in the vial
- Shaking is set too high
- The sample is condensing in the loop

Contaminates the probe and loop

Septum or caps blowing off

- The oven temperature is too high (too much pressure in the vial)

High %RSD

- Vial leaks – check vial crimping; check sequence actions and logbook
- Condensation in the flow path
- Check temperatures
- The vial equilibration time is too short

Can run a leak check

Sequence makes it through first sample only

- GC cycle time is too short; check sequence actions and logbook

Summary

- Stay 10 to 20 °C below the boiling point of the solvent/matrix
- Keep a minimum of 5 mL of headspace in the vial
- Use the Method Development tools
 - Don't forget to turn off the functions when done!
- Optimize parameters based on compounds with the highest K
 - Not every compound responds/reacts in the same way
- To improve sensitivity:
 - Use 10 ml vials
 - Use salts
 - More sample
 - 2 mm liner
- Be consistent with crimping vials; set the crimper properly so that every user is successful
- When troubleshooting:
 - Think logically about what can/cannot cause problem
 - Half-split the system so you know where to focus your energy
- Contact technical support

Additional Resources

[7697A Headspace Sampler Troubleshooting \(PDF\)](#) G4556-90018

[7697A Headspace Sampler Advanced Operation \(PDF\)](#) G4556-90016

Search for 7697A Headspace Sampler on [Agilent.com](https://www.agilent.com)

Contact Agilent Chemistries and Supplies Technical Support



1-800-227-9770 option 3, option 3:

Option 1 for GC and GC/MS columns and supplies

Option 2 for LC and LC/MS columns and supplies

Option 3 for sample preparation products, filtration, and QuEChERS

Option 4 for spectroscopy supplies

Available in the U.S., 8-5 all time zones



gc-column-support@agilent.com




GC columns and supplies

lc-column-support@agilent.com


spp-support@agilent.com

spectro-supplies-support@agilent.com


Agilent 7697A/8697 Loop System



Vial pressurization



Loop fill




Injection


Change in Loop Pressure

First two eluting peaks


Vial fill pressure: 40 psi
Loop fill rate: 30 psi/min
Inlet pressure: 28.3 psi




Is That a Good Way to Increase Signal?



Change in Split Ratio





The Effect of Vial Pressure, Loop Pressure, and Fill Rate





Change the Loop Purge Time and Flow


Carryover issues



Temperatures



Times


Vial and Loop


Carrier



Advanced Functions



Sequence Actions



Method Development

Advanced Functions

Extraction Mode

☒ Single extraction


☐ Multiple extractions


☐ Concentrated extractions


Venting and Purging


☒ Vent vial pressure after extraction

Post-injection purge:

Default

Purge flow: 100 mL/min

Purge time: 1 min



Dynamic Leak Checking

Acceptable leak rate:

Default

Leak flow: 0.2 mL/min

Barcoding of Vials

Barcode symbology:

Enable All

☐ Vial barcodes include checksum