

Preclinical evaluation of a new stereoelectroencephalography - guided radiofrequency ablation system with real time temperature monitoring

Robert E. Gross², Gerald A. Grant⁹, Guy M. McKhann II⁵, Daniel E. Couture⁷, Joshua P. Aronson⁶, David P. Darrow⁸, Jamie J. Van Gompel⁴, Vishwadeep Ahluwalia³, Maria Porto Cruz¹, Maria Vomero¹, Hijaz Haris¹, Camilo Diaz-Botia¹, Mary McNeil¹, Aura Kullmann¹

¹NeuroOne Medical Technologies Corporation, USA; ²Emory University, USA; ³Georgia Institute of Technology, USA; ⁴Mayo Clinic Rochester, USA; ⁵NY-Presbyterian Hospital Columbia University Medical Center, ⁷Wake Forest School of Medicine, ⁸University of Minnesota, ⁹Duke University, ⁶Beth Israel Deaconess Medical Center

INTRODUCTION

sEEG-guided RF ablation is a procedure that uses sEEG electrodes implanted for seizure evaluation to create tissue lesions. The procedure has been performed primarily in Europe for >20 years and has shown to result in seizure free or reduction of seizure frequency with minimal risks¹.

OBJECTIVES

Previous studies have been performed without monitoring the temperature at which ablations were performed. Here we present preclinical in vivo data for the characterization of a new system capable of monitoring and controlling the ablation temperature during the procedure.

METHODS

Six neurosurgeons implanted 18 sEEG electrodes into the brain of three pigs (6/animal). Monopolar and bipolar RF ablations (n=34) were performed using minimum (50°C, 5s) default (80°C, 120s) and maximum (85-90°C, 600s) temperature and time settings (Table 1). MRI and histology were used for lesion characterization. Impedance measurements were performed pre- and post-ablation. User feedback was obtained using a questionnaire.

Fig 1. Swine sEEG implantation





Table 1. Ablation parameters

TC=temperature control; m=manual;

T=temperature; t=time; M=monopolar; B= bipolar


Settings: T (°C), t (s)	Config	N
85°C, 600s (Max)	M, TC	4
85°C, 600s (Max)	B, TC	3
~88-90°C, 600s (Max)	M, m	4
~88-90°C, 600s (Max)	B, m	5
80°C, 120s (Default)	M, TC	6
80°C, 120s (Default)	B, TC	4
50°C, 5s (Minimum)	M, TC	5
50°C, 5s (Minimum)	B, TC	3

sEEG guided RF ablation system (FDA cleared K231675)




RESULTS

Histology: lesions are well delineated



- Monopolar – spherical; Bipolar – ellipsoid
- Necrotic core: coagulative necrosis
- Transition zone: a band of variably vacuolated and degenerated but intact tissue
- Viable tissue: intact neuroparenchyma, no evidence of cellular degeneration

- Minimum settings: no visible lesion
- Default and maximum settings: lesion size correlates with ablation parameters
- Maximal lesion size:
 - width up to 10 mm;
 - height up to 12 mm (bipolar)
 - volume up to ~400 mm³




SUMMARY AND CONCLUSIONS

- The sEEG guided new ablation system is able to create lesions in the swine brain in vivo
- Lesion sizes are proportional to the ablation parameters: temperature and time
- Lesion characterization in MRI and histology:
 - Similar appearance and shape (spherical and ellipsoid)
 - Well demarcated lesions for ablations at default and worst-case scenario parameters
 - No detectable lesion for ablations performed at minimum parameters
 - Histologically distinct zones: necrotic core, transition zone, surrounded by intact parenchyma
- Intracranial EEG recordings are feasible pre- and post-ablation
- The system is intuitive and easy to use
- The only FDA cleared sEEG guided RF ablation system able to monitor and regulate temperature during ablation
- Potential to reduce the number of surgical procedures by providing treatment using the same sEEG electrodes for ablation and diagnostic purposes

RESULTS


In vivo: lesion size is proportional to ablation parameters




- Monopolar – spherical; Bipolar – ellipsoid
- Inner core and outer ring
- Minimum settings: no visible lesion
- Default and maximum settings: lesion size correlates with ablation parameters
- Maximal lesion size:
 - width up to 9 mm
 - height up to 10 mm (bipolar)
 - volume up to ~350 mm³

Monopolar Ablation Parameters	Average Width (mm) ± SD	Average Volume (mm ³) ± SD	n
80°C, 120 s, M-TC	6.04 ± 0.70	126.14 ± 32.49	3
85°C, 600 s M-TC	6.92 ± 1.74	213.04 ± 122.85	4
88°C-90°C, 600 s, M-m	7.38 ± 2.17	232.56 ± 122.18	4

Bipolar Ablation Parameters	Average Width (mm) ± SD	Average Height (mm) ± SD	Average Volume (mm ³) ± SD	n
80°C, 120 s, B-TC	6.27 ± 0.72	7.89 ± 0.71	174.27 ± 46.76	3
85°C, 600 s B-TC	7.32 ± 0.81	8.53 ± 1.78	251.19 ± 110.51	3
88°C-90°C, 600 s, B-m	6.09 ± 1.03	9.17 ± 0.85	202.71 ± 81.38	4



Overall customer feedback is very strong



Recordings pre and post-ablation are feasible

Monopolar Ablation Parameters	Impedance pre- and post-ablation		N*
	Begin Average ± SD (Ohm)	End Average ± SD (Ohm)	
50°C, 5 s, TC	577.04 ± 116.67	515.06 ± 86.25	5
80°C, 120 s, TC	573.82 ± 116.44	509.25 ± 167.68	4
85°C, 600 s TC	552.02 ± 90.88	501.05 ± 83.09	4
88°C-90°C, 600 s, m	421.15 ± 138.40	334.70 ± 113.03	4

Bipolar Ablation Parameters	Impedance pre- and post-ablation		N*
	Begin Average ± SD (Ohm)	End Average ± SD (Ohm)	
50°C, 5 s, TC	768.20 ± 102.73	703.23 ± 97.19	3
80°C, 120 s, TC	505.27 ± 126.62	483.42 ± 186.12	4
85°C, 600 s TC	669.50 ± 209.17	586.30 ± 142.18	3
88°C-90°C, 600 s, m	891.16 ± 126.06	711.94 ± 117.54	5

Funding: NeuroOne Medical Technology Corporation
 Acknowledgements: NeuroOne Team
 Contact: rg1246@rwjms.rutgers.edu; aurak@nmtc1.com