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Use of the birth pool “Made in water”

Content of the suitcase and/or the bags we give you:

- 1 inflatable purple pool
- 1 electric air pump to inflate the pool
- 1 plug adaptor
- 1 single-use liner to place in the pool
- Different rings to attach hose to the tap
- 1 pump to pump water out of the pool
- 1 single-use hose to fill up and empty the pool after use

It is important that you only inflate the pool once, ie. when labour has started. Every time the pool is inflated, counts as a one time use as the material wears out a tiny bit.

Set up:

- Place the plastic sheet on the floor. This is not only to protect your floor, but also to protect the pool from any sharp bits that might damage it.
- Place the pool on the sheet and inflate using the electric air pump.
- Reposition the pool on the sheet if necessary.
- **!!! Place the liner in the pool !!!**
- Fill the pool with water at a temperature that the woman in labour finds comfortable, but not over 38 degrees. Let her check the temperature from time to time.

Break down:

- Place the water pump in the water, attach the hose.
- Place the end of the hose in your bath or shower and turn on the pump.
- The pump is not capable of pumping out the last centimeters of water. You will have to take that out manually.
- Remove the liner and throw away.
- Clean the outside of the pool with a humid cloth and leave to dry.
- Deflate the pool.
- Fold the pool.
- Clean the plastic sheet with a humid cloth and leave to dry. Fold the sheet.

Content of the suitcase and/or bags when you return it:

- 1 plastic sheet to place under the pool
- 1 inflatable purple pool
- 1 electric air pump to inflate the pool
- 1 plug adaptor
- Different rings to attach hose to the tap
- 1 pump to pump water out of the pool