

FFP2 mask  
Ear straps



# FFP2

## Filtration face mask

> 98% filtration of particulate aerosols



### Technical information

**Product name:** Filtration face mask

**Classification:** FFP2 - NR (non-reusable)

**Model:** MSH

**BFE:** > 98 %

**PFE:** > 94 % (>98% NaCl aerosol sprays)

**Size:** 157 x 107 cm

**Standard:** EN 149:2001 + A1:2009

**Composition:** 66.7% nonwoven fabric, 33.3% melt-blown non-woven fabric

**Allergy warning:** Non-woven textile

**Expiration time:** Two years if stored under the specified conditions

### Specifications

Designed for **face perfect fit**

**High breathability**

**Doble welding** that guarantees the resistance of the cords up to 25 newtons

With cord stoppers for **size regulation**

**Ear saver** included

**Individual packaging including:**

- Zip closure to store the mask while not in use
- Hours of use indicator
- Owner identifier
- Antibacterial ventilation system

### Fitting instructions

1. Hold the mask with one hand, place it on your face, covering nose and chin. Pull the top strap over your head. and pass the lower strap behind your head and place it on the nape of the neck.
2. Pull the top strap over your head.
3. Pass the lower strap behind your head and place it on the nape of the neck. Adjust the cord stopper for perfect fit.
4. Using both hands, adjust the nose clip to the shape of the nose and ensure that the edges of the mask are in direct contact with the face.
5. Fit testing for leaks. Allows you to check that the mask is fitted correctly. To do this, cover the mask with both hands and exhale with force. If any leaks are detected at the edges, you will need to adjust the position of the mask against the face, readjust the nose clip and/or the tensile position of the elastic straps. Repeat this process until no leaks are detected.



CE 2163

EN 149:2001 + A1:2009

FFP2 mask  
Ear straps



# FFP2

## Filtration face mask

> 98% filtration of particulate aerosols



### Important information

- NR equipment must not be used for more than one shift.
- Use in accordance with the manufacturer's recommendations, following the instructions on the label.
- These masks do not provide oxygen.
- Not to be used in environments with low oxygen concentration or poor ventilation.
- The presence of special facial features, abundant facial hair or glasses could hinder correct fitting of the mask.
- Keep away from heat sources.
- Indoor storage with less than 80% humidity, free from corrosive gases and well ventilated.

### Individual bag specifications

Quantity: 1 unit/bag  
Dimensions: 130 x 205 mm  
Weight: 12 g

SKU: MASK-ST-00029.0001



SKU: MASK-ST-00013.0001



### 10 pcs box specifications

Quantity: 10 units (individual bags)/box  
Dimensions: 132 x 210 x 45 mm  
Weight: 130 g

SKU: MASK-ST-00029.0010



SKU: MASK-ST-00009.0010



### Master carton specifications

Quantity: 60 box/master carton  
Total quantity: 600 units  
Dimensions: 485 x 415 x 444 mm  
Weight: 9.5 kg

CE 2163

EN 149:2001 + A1:2009