

FFP3

Filtration face mask

> 99% filtration of particulate aerosols



Technical information

Product name: Filtration face mask

Classification: FFP3 - NR (non-reusable)

Model: MSH

Size: 156 x 118 mm

Standard: EN 149:2001 + A1:2009

Composition: 66.7% nonwoven fabric, 33.3% melt-blown non-woven fabric

Allergy warning: Non-woven textile

Expiration time: Two years if stored under the specified conditions

Specifications

Designed for **face perfect fit**

High breathability

Doble welding that guarantees the resistance of the cords up to 25 newtons

With cord stoppers for **size regulation**

Ear saver included

Individual packaging including:

- Zip closure to store the mask while not in use
- Hours of use indicator
- Owner identifier
- Antibacterial ventilation system

Fitting instructions

1. Hold the mask in one hand with the nose clip towards your fingers and the straps hanging free.
2. Place the mask over the face, fully covering the nose and mouth. The nose clip must be over the nose.
3. Using one hand, hold the mask against the chin. Pull the strap and place it behind your ears.
4. Using both hands, adjust the nose clip to the shape of the nose and ensure that the edges of the mask are in direct contact with the face.
5. Fit testing for leaks. Allows you to check that the mask is fitted correctly. To do this, cover the mask with both hands and exhale with force. If any leaks are detected at the edges, you will need to adjust the position of the mask against the face, readjust the nose clip and/or the tensile position of the elastic straps. Repeat this process until no leaks are detected.



CE 0370

EN 149:2001 + A1:2009

FFP3

Filtration face mask

> 99% filtration of particulate aerosols



Important information

- NR equipment must not be used for more than one shift.
- Use in accordance with the manufacturer's recommendations, following the instructions on the label.
- These masks do not provide oxygen.
- Not to be used in environments with low oxygen concentration or poor ventilation.
- The presence of special facial features, abundant facial hair or glasses could hinder correct fitting of the mask.
- Keep away from heat sources.
- Indoor storage with less than 80% humidity, free from corrosive gases and well ventilated.

Individual bag specifications

Quantity: 1 unit/bag
 Dimensions: 205 x 130 mm
 weight: 12 g

10 pcs bag specifications

Quantity: 10 units (individual bag)/bag
 Dimensions: 210 x 131 x 60 mm
 Weight: 135 g

SKU: MASK-ST-00024.0001



SKU: MASK-ST-00024.0010



Master carton specifications

Quantity: 60 box/master carton
 Total quantity: 600 units
 Dimensión: 635 x 415 x 444 mm
 Weight: 11.6 Kg

CE 0370

EN 149:2001 + A1:2009