Filtering half mask FFP2 NR Ear straps

> 98% filtration of particulate aerosols

Technical information

Product name: Mascarilla de filtración **Classification:** FFP2 - NR (No reutilizable)

Model: MSH

Size: 157 mm x 107 mm

Standard: EN 149:2001 + A1:2009

Composition: 66.7% Tela no tejida, 33.3% tela fundida

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Designed for face perfect fit

High breathability

Doble welding that guarantees the resistance of the

cords up to 25 newtons

With cord stoppers for size regulation

Ear saver included

Packed in individual bag

Fitting instructions

- 1. Hold the mask in one hand with the nose clip towards your fingers and the elastic straps hanging free.
- 2. Place the mask over the face, fully covering the nose and mouth. The nose clip must be over the nose.
- 3. Using one hand, hold the mask against the chin. Pull the strap and place it behind your ears or attach it to the ear saver for added comfort.
- 4. Using both hands, adjust the nose clip to the shape of the nose and ensure that the edges of the mask are in direct contact with the face.
- 5. Fit testing for leaks. Allows you to check that the mask is fitted correctly. To do this, cover the mask with both hands and exhale with force. If any leaks are detected at the edges, you will need to adjust the position of the mask against the face, readjust the nose clip and/or the tensile position of the elastic straps. Repeat this process until no leaks are detected.

Filtering half mask FFP2 NR Ear straps

> 98% filtration of particulate aerosols

Important information

- NR equipment must not be used for more than one shift.
- Use in accordance with the manufacturer's recommendations, following the instructions on the label.
- These masks do not provide oxygen.
- Not to be used in environments with low oxygen concentration or poor ventilation.
- The presence of special facial features, abundant facial hair or glasses could hinder correct fitting of the mask.
- Keep away from heat sources.
- Indoor storage with less than 80% humidity, free from corrosive gases and well ventilated.

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EU TYPE EXAMINATION CERTIFICATE

Certificate No: 2163-PPE-707

Respiratory protective devices, filtering half masks to protect against particles manufactured by

Huizhou Hengda Innovation Communication Equipment Co., Ltd.

Building A, Wanli Industrial Co., Ltd., Dalongkeng, Ganpo, Zhenlong Town, Huiyang District, Huizhou City, Guangdong Province, China

are tested and evaluated according to

EN 149:2001 + A1:2009 Respiratory Protective Devices -Filtering Half Masks to Protect Against Particles -Requirements, Testing, Marking

Based on the type examination conducted with the evaluation of test reports, technical file according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 5, it is approved that the product meets the requirements of the regulation.

Product Definition

Brand Name: KSA Model: MSH Filtering half mask Classification: FFP2 NR

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as shown below, on the Category III product models given above, with;

- Issuing an appropriate EU Declaration of Conformity according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 9.
- Ongoing successful performance in fulfilment of the requirements set out in Personal Protective Equipment Regulation (EU) 2016/425 and harmonised standards, ensured by assessments based on Annex 7 (Module C2) or Annex 8 (Module D) of the regulation no later than 1 year from the beginning of serial production

This certificate is initially issued on 04/06/2020 and will be valid for 5 years, if there is no change in the relevant harmonised standard affecting the essential health and safety requirements.

Suat KAÇMAZ
UNIVERSAL CERTIFICATION
Director

Mascarilla de filtración FFP2 - NR

Certificado CE

CERTIFICATE OF CONFORMANCE

Certificate No: 2163-PPE-707/01

Respiratory protective devices, filtering half masks to protect against particles manufactured by

Huizhou Hengda Innovation Communication Equipment Co., Ltd.

Building A, Wanli Industrial Co., Ltd., Dalongkeng, Ganpo, Zhenlong Town, Huiyang District, Huizhou City, Guangdong Province, China

Continues to fulfil the requirements of

EN 149:2001 + A1:2009 Respiratory Protective Devices -Filtering Half Masks to Protect Against Particles -Requirements, Testing, Marking

Based on the evaluation of test reports and internal quality control audit reports according to EN 149+A1:2009 and Personal Protective Equipment Regulation (EU) 2016/425 Annex VII (Module C2). This certificate implies that the manufactured products show below are in conformance with the approved EU Type Examination model and meets the requirements of the regulation.

Product Definition

Model	Class	EU Type Examination Certificate		
		Serial No	Date	Issuing NB No
KSA/MSH	FFP2 NR	2163-PPE-707	04.06.2020	2163

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as shown below, on the Category III product models given above, with;

- Issuing an appropriate EU Declaration of Conformity according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 9.
- Taking all measures necessary so that the manufacturing process and its monitoring ensure the homogeneity of production and conformity of the manufactured PPE with the type described in the EU type examination certificate.

This certificate is issued on 21/06/2020 and will be valid for one year, until 20/06/2021 if the manufacturer makes no major change in the product designs and manufacturing processes affecting the product performance on the essential health and safety requirements.