



UNSURPASSED PRODUCTIVITY

Volatile Metabolites as Early Diagnostics for Blood-Borne Pathogens



Dr. Yatin J. Mange
Syft Technologies Ltd
September 2, 2016
yatin.mange@syft.com



www.syft.com



www.syft.com

Summary

Introduction to SIFT-MS

SIFT-MS for the detection of blood-borne pathogens

Conclusion

Selected Ion Flow Tube Mass Spectrometry (SIFT-MS) – Key Attributes



Analytical technique for measuring Volatile Organic Compounds (VOCs)

Rapid and real-time measurement


Selectivity, sensitivity and speed for a diverse range of chemical species

Soft ionization resulting in less fragmentation

Detection limit up to pptv

Easy to operate

SIFT-MS – how this soft chemical ionization technique works



Multiple reagent ions

Pure reagent ion delivery

Ultra-soft sample ionization

Mass spectrometer


SIFT-MS can measure most VOCs and inorganic compounds

hydrocarbons	alkanes, alkenes, aromatics, monoterpenes
oxygenates	alcohols, aldehydes, ketones , esters, ethers, carboxylic acids, formaldehyde
nitrogen compounds	amines, amides, nitriles , nitrated organics
sulfur compounds	mercaptans, thioethers, carbonyl sulfide
halogenated compounds	aliphatic and aromatic fluorides, chlorides, bromides and iodides
inorganics	ammonia, hydrogen cyanide, hydrogen sulfide , hydrogen chloride, hydrogen fluoride, carbon dioxide , sulfur dioxide,


Who is using Syft Instruments



SIFT-MS is extremely sensitive and has a dynamic range of 6 orders of magnitude



Concentrations can be monitored in real time



SIFT-MS for the detection of blood-borne pathogens

SIFT-MS - Biomedical applications



Breath diagnosis – non-invasive, point-of-care diagnostics and monitoring:

- Diagnosis of Chronic Obstructive Pulmonary Disease (COPD) and asthma.
- Determination of dialysis efficacy
- Diagnosis of acute renal failure and transplant rejection

Laboratory – high throughput, low operator skill:

- Detection and identification of bacteremia in blood culture
- Diagnosis of urinary tract infections
- Detection of solvents and toxic industrial chemicals in urine

SIFT-MS application – VOCs from Blood-Borne Pathogens



VOCs present as a result of metabolic or disease processes

Present in breath, or headspace over body fluids or cultures


Conventional automated systems detect CO₂ liberated by growing microbes.
(Eg: BacT/ALERT techniques)

Conventional techniques – Time consuming, trained technician

Timely detection enables more accurate treatment

In the current study, 5 bacterial species detected from 10 CFUs

Identification of bacterial pathogens from blood cultures



Bacterial VOCs at 6 hours in BacT/ALERT medium

Bacterial Species

- Pseudomonas aeruginosa*
- Streptococcus pneumoniae*
- Escherichia coli*
- Staphylococcus aureus*
- Neisseria meningitidis*

SIFT-MS decreased time-to-positive result by as much as one half


Enhanced sensitivity enables detection of bacterial growth in blood culture in 6 hrs

VOCs from each bacterial culture headspace completed in 30 s without prior sample preparations

Profiles of MVOCs allows differentiation between different species of bacteria

Potential for integration into automated systems

Automated SIFT-MS analysis – VOCs from headspace



Summary

SIFT-MS technology is adaptable, testing a range of body samples (including breath, blood and urine)

Volatile compounds produced by microbes are measured quickly and accurately

Unique set of volatile metabolites indicates particular species of bacteria

Shorter incubation time and rapid analysis – Powerful early screening tool



UNSURPASSED PRODUCTIVITY



Questions?

Email: yatin.mange@syft.com



Find out more at
Exhibition Hall 1
(1F)
Poster No: PA-94

