boste
purea

SPECTRA **LED Face Mask**

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English

Product Specifications

Product name: Spectra LED Face Mask
with Remote Control Black
Model: SM2309
Dimensions: 324x214x4.3mm
Net weigh: 170g
Included in the set:

- LED face mask
- Remote control
- TYPE-C charging cable
- Silicone safety goggles
- Adjustable nylon strap

Charging voltage: 5V/1A
Operating current: ≤ 3.5A
Rated power: ≤ 15W
Battery capacity: 2600 mAh (lithium)
Light colors: 7 colors
Water resistance rating: IPX7
Country of origin: Chiny

Product Description

The LED photorejuvenation mask is a modern, silicone-based device designed for home use. It is lightweight, convenient, and easy to operate. Thanks to its flexible material, it comfortably adapts to the shape of the face, and the adjustable strap ensures a perfect fit for any user.

The mask is equipped with 472 LED beads, divided into 118 sets. Each LED emits light in one of four specific wavelengths:

- Red (630 nm) – improves skin elasticity
- Blue (460 nm) – helps reduce blemishes and inflammation
- Yellow (590 nm) – evens out skin tone and brightens
- Infrared (850 nm) – penetrates deeper skin layers

The set includes a remote controller, which allows you to select 3 energy levels and customize the lighting time. The LED light supports natural regenerative processes in the skin, stimulates collagen and elastin production, and with regular use, helps visibly improve skin tone, texture, and radiance.

Photon Rejuvenation instrument

Controller

Type-C Cable

Silicone safety goggles

Adjustable nylon buckle

Light Modes

-  Red Light (630 nm): Level 1, 2, 3
-  Blue Light (460 nm): Level 1, 2, 3
-  Yellow Light (590 nm): Level 1, 2, 3
-  Mixed Light (Combination): Level 1, 2, 3

Infrared light can be controlled independently or used in combination with other light modes.

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Red Light – 630 nm
Anti-aging effect,
supports collagen production,
firms the skin, and improves elasticity.

Blue Light – 460 nm
Sterilizes the skin, eliminates
acne-causing bacteria, reduces inflammation
and scars, and increases skin oxygenation.

Yellow Light – 590 nm
Evens skin tone,
boosts blood circulation, reduces redness,
puffiness, and fine lines.

Infrared – 850 nm
Regenerates the skin, i
mproves blood circulation, enhances hydration,
and supports deep skin repair.

Important Information

Before first use, make sure to fully charge the mask's battery. During charging, the LED indicator on the controller will flash. Once fully charged, the display will remain lit. Estimated charging time: approximately 4 hours.

The mask is equipped with an adjustable hook-and-loop strap, allowing comfortable and stable fit on the head.

How to put on the mask:

Thread the strap through the openings on opposite sides of the mask.
Wrap the strap around your head.
Press and fasten the strap securely.

Device Controller

Device Operation

- 1. Power button:** Press and hold for 2.5 seconds to turn the device on or off. When turned on, the controller indicator light will turn red.
- 2. "E" Button (Energy Level):** Short press to change the energy level (up to 3 levels of output).
- 3. Timer setting:**
Press the T+ button to increase the session time.
Press the T- button to reduce the time in 5-minute intervals.
Available durations: 10 / 15 / 20 / 25 / 30 minutes
It is recommended to use the device for 15–25 minutes per session.
- 4. "M" Button (Light Mode Selection):**
Short press to cycle through the available light modes:
blue light / yellow light / combination / red light / light cycle

The display will show the corresponding symbols:

Red light (R)
Blue light (B)
Yellow light (Y)
Combination (G)
- 5. "NIR" Button (Infrared Light):** Short press to turn infrared light on or off.
- 6. Charging Status:** The battery indicator flashes while charging. When fully charged, the indicator turns off.

Precautions

- 1. Use the device strictly according to the manual.**
Improper use may lead to equipment damage or unexpected reactions.
- 2. Do not use the device for purposes other than those intended.**
This product is designed exclusively for cosmetic skincare treatments.
- 3. Not suitable for pregnant or breastfeeding women, or children.**
- 4. Avoid using if you have skin conditions or other serious health issues.**
If you suffer from autoimmune diseases, light sensitivity (photosensitivity), lupus, eczema, albinism, or neurological/mental disorders — the use of this device may be contraindicated. Always consult a doctor (e.g. dermatologist) before use.
- 5. Do not use if you have light sensitivity (photosensitivity).**
It may cause severe skin reactions in individuals with this condition.
- 6. Always protect your eyes during use.**
Use protective goggles and keep your eyes closed. If you have a genetic eye disease — do not use LED masks.
- 7. If you feel any discomfort — stop using the device immediately.**
- 8. The device is not waterproof.**
Do not immerse in water or rinse. If water contact is suspected — do not use.
- 9. Keep out of reach of children and infants.**
Store safely.
- 10. Do not use if the product is damaged.**
Damaged components can cause injury or malfunction.
- 11. For home use only.**
Do not disassemble, modify, or repair the device yourself — it may lead to malfunction or damage.
- 12. Recommended usage frequency:**
3–4 times per week, for 15–25 minutes per session.

Cleaning and Maintenance

- After each use, wipe the surface of the mask with a dry, soft cloth (e.g. microfiber) to remove sweat, dust, and residue.
- Do not use water, detergents, alcohol, or chemicals — they may damage the silicone and LED surface.
- Do not immerse the controller or mask in water. The product is not waterproof.
- Store the product in a cool, dry place, away from heat sources, sunlight, or moisture.
- Do not bend the mask sharply — it may damage internal LED circuits.
- Do not pull on cables or press hard on buttons or display components.
- If the device is not in use for a longer time — recharge the battery at least once every 1–2 months to avoid deep discharge.

Frequently Asked Questions

Can I use the mask every day?

– It is recommended to use it 3–4 times a week, for 15–25 minutes per session. Using it too frequently may irritate the skin.

Can I use the mask after applying cosmetics?

– Do not use the mask directly after applying products containing acids, retinol, or other photosensitizing ingredients. Wait a few hours or perform the session on clean skin.

Can I use the mask around the eyes?

– It is recommended to close your eyes during the session and use the included protective goggles. Do not shine the light directly into your eyes.

Is the device waterproof?

– No. The mask and controller are not water-resistant. Do not immerse in water or rinse under running water.

How long does charging take?

– A full charge takes approximately 4 hours. The LED indicator will stop blinking when the battery is fully charged.

Can pregnant or breastfeeding women use the mask?

– Use is not recommended for pregnant or breastfeeding women without prior consultation with a doctor.

Can children use the mask?

– No. The device is intended for adult use only.

What do the light colors mean?

– Each color has a different effect: Red – anti-aging, Blue – anti-acne, Yellow – brightening, Infrared – regenerating

This product complies with the requirements of the following EU directives:

Directive 2014/30/EU – Electromagnetic compatibility

RoHS Directive 2011/65/EU – Restriction of hazardous substances

The crossed-out wheeled bin symbol means

that within the European Union, used electrical and electronic equipment should not be disposed of with household waste.

It should be taken to appropriate collection points for recycling.

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