ideal

IZ9102

INTELLIGENTNA WAGA DLA NIEMOWLĄT
SMART BABY SCALE
INTELLIGENTE BABYWAAGE
ИНТЕЛЛЕКТУАЛЬНЫЕ ДЕТСКИЕ ВЕСЫ
BALANZA INTELLIGENTE PARA BEBÉS
BALANCE INTELLIGENTE POUR BÉBÉ
BILANCIA INTELLIGENTE PER BAMBINI

Eldom Sp. z o.o. • ul. Pawła Chromika 5a • 40-238 Katowice, POLAND
tel: +48 32 2553340 • fax: +48 32 2530412 • www.eldom.eu

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/EN/ The disposal of used electrical and electronic equipment (applies to European Union countries and other European countries with separate waste-collection systems).

This symbol on the product or its packaging indicates that it should not be classified as household waste. It should be handed over to an appropriate company dealing with the collection and recycling of electrical and electronic equipment. The correct disposal of the product will prevent potential negative consequences for the environment and human health resulting from hazardous substances present in the product. Electrical devices must be handed over to restrict their re-use and further treatment. If the device contains batteries, remove them, and hand them over to a storage point separately. DO NOT THROW EQUIPMENT INTO THE MUNICIPAL WASTE BIN. Material recycling helps to preserve natural resources. For detailed information on how to recycle this product, please contact your local authority, the recycling company, or the shop where you bought it.

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SAFETY INSTRUCTIONS

Carefully read the entire contents of this manual before first use.

Keep all warnings and safety instructions for future reference.

1. The equipment is intended for domestic use only.
2. This equipment may be used by children of at least 8 years of age and by persons with reduced physical and mental capabilities and persons with inexperience and lack of familiarity with the equipment, if supervision or instruction is provided on how to use the equipment in a safe manner so that the associated risks are understood. Children should not play with the equipment. Unsupervised children should not carry out cleaning or maintenance of the equipment.
3. Do not exceed the maximum allowable weight (25kg) - this may damage the scale or cause it to calibrate.
4. Handle the scale with care. Do not drop the scale.
5. Use a damp, soft cloth without chemicals to clean, then wipe dry.
6. Do not immerse the equipment in water.
7. Repairs to the equipment may only be carried out by an authorised service centre. Any upgrades or the use of other other than original spare parts or components of the device is prohibited and jeopardises the safety of use.
8. Store the scale in a shaded and dry place.
9. Do not press on the scale with other objects during storage.
10. The accuracy of the scale indication depends on the way it is used.
11. Eldom sp. z o.o. shall not be held liable for any damage resulting from improper use of the device.

WARNING: Improper use of the appliance may cause injury.

CAUTION: The plastic bag can be a hazard - to avoid suffocation with the bag, keep it away from babies and small children.

Hereby, Eldom Sp. z o.o., 40-238 Katowice, ul. Pawła Chromika 5a declares that the radio equipment type **Ideal IZ9102** is in compliance with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at the following internet address: **www.eldom.eu/wsparcie/deklaracje**

OPERATING INSTRUCTION

SMART BABY SCALE IZ9102

EN

• GENERAL DESCRIPTION - Fig. 1

1. Button  ZERO
2. Display 
3. Button 
4. Weighing platform
5. Measuring scale
6. Measuring scale lock
7. Battery cover
8. Feet

• Display (2) – Fig. 2

- A. Weight measurement
- B. Comparison marks
- C. Measurement units

• TECHNICAL DATA

- maximum weight: 25kg
- accuracy of measurement: 5g
- power supply: 4xAAA battery (excluded)

• INTENDED USE

The **Ideal IZ9102** scale allows for measuring the body weight of infants and small children (up to 25 kg), saving the results on the device, and transferring them to an app. In the app, you can also record the child's height and head circumference. Measurements can be taken using the built-in measuring tape. The device is safe for children and easy to use.

NOTE.

The images, colors, accessories, user interface, and other illustrations in the manual are for reference only. Due to product improvements, the actual device may differ slightly from the illustrations. Please refer to the actual product.

• TIPS

- It is recommended that you take measurements at the same time daily (e.g. between 7:00-9:00 p.m.).
- Place the scales on a firm and flat surface.

• USE

- Before use, check that the unit is complete and not damaged.
- Remove all transport protection.
- Insert the batteries.

• BATTERIES

The scale is powered by 4 AAA batteries.

Battery replacement - Fig. 3

- Open the battery cover located on the bottom of the unit.
- Gently lift up and remove the batteries.
- Insert new batteries. **Pay attention to the polarity.**
- Close the cover.

Used batteries should be disposed of in special containers.

• OPERATION

Taking measurements without using the app – standard weighing mode (Fig. 4)

- Place the scale on a stable, dry surface.
- Turn on the scale using the button (1).
- Place the child on the platform (4).

WARNING: During weighing, always ensure the child does not slip off the platform (4).

- Wait a few seconds, then press the button (3). The weight reading will stabilize, and the measurement will be saved in the device's memory.

WARNING: Measurement results are stored in the device until the batteries are replaced.

- Read the display (2).

- On the display (2), in addition to the measurement result, one of the indicators "**B**" will appear. The comparison marks "**B**" show whether the child's body weight has increased or decreased compared to the previous measurement.

Built-in measuring scale (5)

- The measuring scale can be used to measure the child's height and head circumference.
- To take a measurement, open the silicone cover and pull out the measuring scale (5).
- The measurement result can be saved in the dedicated app.
- To retract the measuring scale (5), press and hold the lock button (6). The measuring scale (5) will automatically return to its place. Then close the silicone cover.

Taking measurements without using the app

NOTICE. Due to the continuous development of our products, the app and its appearance may change over time.

The scale connects with devices that support wireless communication with at least Bluetooth 4.0 and operating systems iOS 12.0/Android 6.0 or newer.

1. Download the app using the code located on the cover of the manual or search for the keyword "iMoni Health" in the "App Store" or "Google Play".
2. For the first use, create a new user account.
3. Ensure that Bluetooth is enabled on your smartphone.
4. Launch the app.
5. Measurements will appear on the scale's display and in the app.
6. All measurements will be automatically saved in the app's data history.
7. In the app, you can manually enter the child's head circumference and height data.

The app allows you to monitor your child's development by tracking changes in body weight, height, and head circumference. Detailed explanations of the indicators are available in the app.

Automatic switch-off

The scale will automatically turn off after approximately 2 minutes of inactivity.

Changing the weighing units

You can change the weighing unit (kg/lb) by pressing and holding the button.

Low battery indicator – signals the need to replace the battery.

Weight overload message - do not exceed the maximum allowable weight (25 kg). This may damage the device.

• CLEANING AND MAINTENANCE

- Do not immerse in water.
- Use a damp, soft cloth without chemicals to clean, then wipe dry.
- Wipe the housing with a moistened cloth.

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