molti

ZM001

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EN | USER MANUAL

1. Unpacking the set:

Open the package, ensuring that all components of the set (roller and balls) are included.

Check for any visible damage.

2. Preparing for massage:

Choose a suitable surface for massage: hard but comfortable to ensure effectiveness and comfort.

Decide which body parts you want to massage (e.g., back, thighs, calves, feet).

3. Using the roller:

Position: Sit or lie on the floor.

Massage: Place the roller under the muscles you want to massage. Slowly move your body back and forth to roll out the muscles.

Pressure: Adjust the pressure on the roller according to your comfort level. You can increase the pressure for a deeper massage.

4. Using the balls:

Position: Sit or lie on the floor.

Massage: Place the ball under the muscles you want to massage (e.g., feet, glutes, back). Gently move around to apply pressure on tight areas.

Adjustment: You can use different balls in the set to tailor the massage to specific body parts.

5. After the massage:

Make sure to rest and drink water to hydrate your body after the massage.

EN | WARNING CARD

- **Do not use in case of injury:** Do not use the roller or balls if you have open wounds, bruises, fractures, or other injuries.
- **Consult a doctor:** If you have health issues, chronic muscle pain, or other complaints, consult a doctor before starting the massage.
- **Avoid excessive pressure:** Do not apply too much pressure on body areas to avoid bruising and discomfort.
- **Do not use on the spine:** Do not apply the roller or balls directly on the spine to avoid injury.
- **Storage:** Store the set in a dry place, away from heat sources, to maintain its quality.

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