ZTB6010W

zelmer

ZTB6010W

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WE WOULD LIKE TO THANKS YOU FOR CHOOSING ZELMER, WE WISH THE PRODUCT PERFORMS TO YOUR SATISFACTION AND PLEASURE

WARNING

PLEASE READ THE INSTRUCTIONS FOR USE CAREFULLY PRIOR TO USING THE PRODUCT. STORE THESE IN A SAFE PLACE FOR FUTURE REFERENCE.

DESCRIPTION

1. Protective cap
2. Brush head
3. Main body
4. ON/OFF button
5. Clean mode indicator (CLEAN)
6. Massage mode indicator (REFRESH)
7. Whitening mode indicator (WHITE)
8. Removal of plaque and tartar mode indicator (GUM CARE)
9. Sensitive mode indicator (SENSITIVE)
10. Charging indicator
11. Charging base
12. USB power wire

SAFETY INSTRUCTIONS

The sonic toothbrush is intended to be used solely by adults for oral care. It is intended for private indoor use only. It must not be used for professional purposes. It is not intended for industrial or commercial use. It must not be used at healthcare or veterinary facilities. Following the instructions in this manual will ensure safe installation and use of the product.

The device should not be used or maintained by children younger than 8 or by individuals with limited mental, intellectual, or physical abilities, unless such individuals are supervised or have

been instructed how to use the device safely, in compliance with the instruction manual, and understand the hazards involved.

The device is not a toy. Children must not play with the device. Children must not clean or perform maintenance of the device unless they are supervised by an adult.

The device can only be powered with low voltage electric current matching the markings on the device.

Unplug the power supply plug prior to assembly, cleaning, or maintenance of the device components.

The charging base power cord is not replaceable. If the power cord becomes damaged, do not attempt to repair it, but replace the charging base with a new one instead. Remember to replace the charging base with an original one to avoid hazard.

Warning! Do not use any other charging base to charge the device. Moreover, do not use the charging base included for charging other products, otherwise burns or a fire due to a short circuit may occur.

Warning! Clean the power supply plug

regularly to avoid dust accumulation. If this happens, moisture

can damage the insulation, which can result in a fire. Disconnect the power supply, and then wipe the plug with clean and dry cloth.

Do not use the device for any purpose other than its intended use.

The device is intended solely for private indoor use. It is not suitable for commercial, industrial, medical, or veterinary applications.

Do not use the device in locations where the use of electrical devices is forbidden, and where it could result in harmful interference or other hazards.

Do not turn on or use the device in the vicinity of flammable materials.

Do not use the device for cleaning other surfaces than teeth, gums, and tongue.

Warning! Keep the device out of the reach of children and infants. Do not allow children and infants to use the product on their own. Otherwise, it may result in an accident or injury caused by accidental swallowing or aspiration of moving and

non-fixed parts.

Danger! Regularly check the toothbrush head for cracks. A broken toothbrush head may pose a risk of choking.

Warning! Discontinue using the product if any of its parts are damaged (the head, the device body, the charging base).

Warning! If a failure occurs during operation, turn off and unplug the device. Otherwise, a fire or electric shock may occur.

Warning! If one of the following faults or failures occurs:

- the device or the charging base are deformed or excessively hot
- the device or the charging base give off a burning odor
- the device or the charging base emit abnormal sounds during operation or charging, discontinue using and contact the service facility.

Prior to using the product, consult a dentist if any of the following circumstances applies to you:

- an oral surgery performed within the last 2 months,
- an implanted pacemaker or another

device in the body,

- gum bleeding when using the product,
- any other medical conditions.

Caution! Individuals with severe periodontal disease, history of tooth disorders, or suspected oral cavity disease should consult a doctor prior to using the device, as the toothbrush may cause discomfort to the teeth and gums.

Caution! Individuals with a sensitive oral cavity should not use the product.

Caution! Do not press the toothbrush head against the teeth or gums too strongly, do not clean for a prolonged period of time.

Caution! Do not share the toothbrush with your family members or other persons, as it may result in injuries to the teeth and gums.

When using the device, make sure that next to the brush location there are no items which may get damaged, broken, or overturned.

The charging base power cord should be routed so that it does not touch very hot objects. The cord should not be exposed to crushing, pulling, or other damages.

The charging base power cord should not hinder the movement of people and pets, as one might trip over it, fall down, and suffer bodily injuries.

Danger! Do not clean the charging base when it is connected to the power source.

Do not expose the charging base and its power cord to water and other liquids, as it may negatively affect the operation and life of the device components, result in a short circuit, or even electric shock.

Danger! Do not immerse the device or the charging base in water.

Caution! When the charging base remains unused, disconnect it from the power source. Otherwise, a leaking battery can hinder the insulation properties and result in electric shock or fire.

Danger! Prior to plugging in the charging base, make sure it is completely dry.

Warning! Prior to plugging in, make sure the charging base operating voltage matches its rated voltage. Otherwise, fire or electric shock may occur.

Danger! The product includes a built-in battery. Do not damage or destroy

the battery, do not put it into fire. Keep the battery away from sources of high temperature. When exposed to high temperature, the product may overheat, catch fire, or explode.

Avoid exposing the battery to very low or very high temperatures (below 0°C/32°F or above 40°C/104°F). Extreme temperatures can negatively affect the battery capacity and lifespan.

Avoid exposing the battery to contact with liquids and metal items, as it may partially or completely damage the battery.

Warning! Do not remove or open the battery, otherwise fire, electric shock, or injuries may occur.

Warning! The battery may only be removed prior to the final disposal of the product. Put the worn or damaged battery in a special container.

Warning! Do not disassemble the device except for final disposal, as it may result in fire, electric shock, or injuries.

Protect the device against vibration and shock during transport.

Keep the user's manual and the packaging,

if possible.

If the device is passed on to another user, it must be accompanied by this manual.

IMPORTANT WARNINGS

Any incorrect use or improper handling of the product shall render the warranty null and void. Prior to plugging in the product, check that your mains voltage is the same as the one indicated on the product label. The mains connection cable must not be tangled or wrapped around the product during use. Do not use the device, or connect and disconnect to the supply mains with the hands and / or feet wet.

Do not pull on the connection cord in order to unplug it or use it as a handle.

In order to prevent any risk of danger, do not open the device. Only qualified technical personnel from the brand's Official technical support service may carry out repairs or procedures on the device.

B&B TRENDS SL disclaims all liability for damages that may occur to people, animals or objects, for the non-observance of these warnings.

BEFORE THE FIRST USE (Important)

1. Remove the product from packaging and check all the accessories according to the list given in the "DESCRIPTION" section. Caution! For the sake of children's safety, do not leave any packaging materials accessible and unattended (plastic bags, cardboard boxes, Styrofoam, etc.). Risk of suffocation!

2. Clean all the detachable items.

3. Charge the device.

CHARGING THE DEVICE

Before using it for the first time, please turn on the device, if you feel that the battery is low or if the device was not used for more than 6 months, please charge it for more than 6 hours. Then another charging (after about 45 days) continue for 4 - 6 hours.

CAUTION! Prior to handling the charging base and charging the battery in the device, make sure your hands, the device, and the charging base are completely dry.

Gently unwind the charging base power cord and put the base on a stable, level surface, so that the cord can reach the charging socket. Connect the charging base USB cord plug to a charging socket, e.g., in a computer, or to a mains adapter, such as one typically included with a phone / smartphone / tablet and equipped with a USB socket providing the recommended charging parameters: 5V  1A.

Put the device on the charging base. It will take about 6 hours to fully charge the device. The optimal ambient conditions for charging are 0-40°C. When the ambient humidity or temperature is too low or too high, the battery may not be charged correctly. The device heats up during charging and operation, it is normal. The LED charging indicator flashes slowly and rhythmically during charging. Once the battery is completely charged, the LED indicator will illuminate with steady light. When the battery is about to be discharged, the

LED indicator will flash rapidly to remind about the need to charge. To keep the battery sufficiently charged, you may put the toothbrush on the charging base whenever the brush is not being used. This will not result in overcharging.

OPERATION MANUAL

1. Attaching or Replacing the Brush Head

Attaching the Brush Head:

Align the brush head (2) so that the bristles face the front of the ON/OFF switch (4).

Firmly insert the brush head vertically onto the metal connecting shaft until it is securely in place.

Replacing the Brush Head:

To remove the brush head (2), pull it straight out vertically. Do not twist the brush head when removing it, as twisting may damage key parts of the toothbrush.

When replacing the brush head, ensure you purchase one from this product series. Brush heads from different electric toothbrush models are not interchangeable, so please pay attention when purchasing.

2. After inserting the toothbrush head (2) into the toothbrush handle (3), there will be a gap of 1–2 mm. This gap is designed to allow the toothbrush head to move freely. Please avoid applying excessive force on the toothbrush head, as this may damage the toothbrush or its components.

Note: the brush bristle must be at the same face with ON/OFF switch (4). CAUTION! It is only allowed to assemble and disassemble the set components when the device is turned off.

3. Put the small amount of the tooth paste onto the brush head (2).

4. Soak it first and put it into the mouth, switch the device on by pressing the ON/OFF button (4). Choose the appropriate mode according to your needs:

(5) Clean mode indicator (CLEAN)

Applicable for: Most people

For the problem: Designed to address various oral health issues by cleaning teeth, enhancing their self-protection, reducing gum bleeding and odor, and effectively removing plaque.

(6) Massage mode indicator (REFRESH)

Applicable for: People with bleeding gums.

For the problem: Provides a professional and safe massage experience for the gums, promoting microcirculation, preventing gum atrophy, and reducing gum bleeding and mouth ulcers. It ensures that teeth are not sensitive to hot and cold temperatures.

(7) White mode indicator (WHITE)

Applicable for: People with teeth discoloration caused by pigments such as smoke and tea stains.

For the problem: Targets tooth discoloration caused by pigments such as smoke, tea, and coffee stains.

(8) Removal of plaque and tartar mode indicator (GUM CARE)

Applicable for: People with tartar.

For the Problem: Utilizes a targeted vibration mode to prevent plaque calcification, effec-

tively remove stubborn tartar, polish teeth, and reduce calculus.

(9). Sensitive mode indicator (SENSITIVE)

Applicable People: People with sensitive gums and dentin.

For the problem: Effectively cleans dental plaque and tartar while caring for and massaging periodontal tissue. Strengthens gums, making it especially suitable for users with sensitive teeth.

5. Brushing the teeth

- Point the brush bristles at an angle of 45 degrees, so that the bristles can clean the side surfaces of the teeth and gums.
- Run the brush tip slowly, pressing it very gently against the surface of the teeth and the line where the teeth connect with the gums.
- The standard brushing cycle lasts 4 x 30 s. It is recommended to divide the surfaces of the teeth to brush each of them for 30 s – the top (front and back) and the bottom (front and back). Then clean the contact surfaces of the upper and lower jaw teeth.
- The brush will remind you to change the surface every 30 seconds by pausing for a moment.
- For the first few times you use the brush, the gums may bleed slightly (this happens rarely), or you may feel discomfort (this is common in individuals who have not used mechanical toothbrushes before). It usually happens when the gums are not used to the new brushing technique and the high vibration speed. If this occurs, it is recommended to use the “Sensitive” program. If after 2 weeks of using the brush as recommended the gums keep bleeding, discontinue using the device and consult a dentist.
- Switch OFF the device by pressing and holding the ON/OFF button (4) for 2 seconds.
- Rinse your mouth with water after brushing.

CLEANING AND MAINTENANCE

Proper and regular cleaning will ensure safe operation and extend the lifespan of the device. Caution! Do not clean the device with any chemicals, alkali, abrasive, or disinfecting agents, as they may damage the surface of the device.

Do not immerse in water! Warning!

Cleaning of the device

1. The cleaning should not be made by children.
2. Before cleaning and maintenance, make sure the device is switched off.
3. Disconnect the brush head (2) from the main body (3).
4. Rinse the brush head (2) in running lukewarm water.
5. Do not wash it in hot water, the temperature of which is above 80°C!
6. Shake to remove excess water.
7. Do not pull by the brush bristles.
8. It is recommended to change the brush head (2) every two months.
9. Rinse the main body (3) with running lukewarm water.
10. Do not scratch or scrub the main body (3) to avoid damage to the rubber gasket, which provides water tightness.

11. Never immerse the whole device in water.
12. Wipe dry and put back on the charging base or in a safe place.
13. Prior to connecting for charging, make sure that the device is completely dry.

PRODUCT DISPOSAL

This product complies with European Directive 2012/19/EU on electrical and electronic devices, known as WEEE (Waste Electrical and Electronic Equipment), provides the legal framework applicable in the European Union for the disposal and reuse of waste electronic and electrical devices. Do not dispose of this product in the bin, instead going to the electrical and electronic waste collection centre closest to your home.

The product may contain batteries. Remove them before disposing of the product and dispose in special containers approved for this purpose.

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Telefon: 0769 224 338 / 0264 434 242

EN / WARRANTY REPORT

B&B TRENDS, S.L. guarantees compliance of this product for the use for which it is intended for a period of two years. In the case of breakdown during the term of this warranty, users are entitled to repair or else the replacement of the product at no charge if the former is unavailable, unless one of these options proves impossible to fulfill or is disproportionate. In this case, you can then opt for a reduction in price or cancellation of the sale, which must be dealt with directly with the sales vendor. This also covers replacement of spare parts provided that the product has been used according to the recommendations specified in this manual for both cases, and has not been tampered with by any third party that is not authorised by B&B TRENDS, S.L. The warranty will not cover any parts subject to wear and tear. This warranty does not affect your rights as a consumer in accordance with the provisions in Directive 1999/44/EC for member states of the European Union.

USE OF WARRANTY

Customers must contact a B&B TRENDS, S.L. authorised Technical Service for repair of the product. Since any tampering of the same by anyone not authorised by B&B TRENDS, S.L., or the careless or improper use of the same shall render this warranty null and void. The warranty must be fully completed and delivered along with the receipt or delivery docket for the effective exercise of rights under this warranty.

This warranty should be retained by the user as well as the invoice, receipt or the delivery docket to facilitate the exercise of these rights. For technical service and after-sales care outside the Polish territory, please submit your query to the point of sale where you purchased the item or check for further aftersales information at www.zelmer.com

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