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EN - Massage Pillow

User Manual:

1. The package includes a power adapter and an additional adapter for a car cigarette lighter. Choose the adapter suitable for your usage area and ensure it is firmly plugged into the massager's power socket.
2. To turn on the massager, use the switch. To set the direction of the massage heads' rotation, choose the appropriate position:
 - I – outward direction
 - II – inward direction.
3. To turn off the massager, set the switch to the 0 position.
4. Always unplug the adapter when the product is not in use.

Warnings:

Not suitable for use by children or elderly individuals. Do not use during pregnancy, immediately after childbirth, after surgical procedures, or in cases of spinal injury, heart disease, rheumatoid arthritis, circulation problems, nerve damage, sensory disorders, or any skin conditions. It is recommended to consult a doctor regarding the safe use of this product. Discontinue use immediately if discomfort occurs during operation. Use the massager in 30-minute cycles, with at least a one-hour break between uses. Do not use in bed or during sleep. Do not use in combination with heating devices or warming agents, including warming creams. Do not use on damaged or irritated skin. Do not cover the massager during use. Do not use the massager near flammable materials. Do not open the massager while it is in operation – moving parts may cause injury. Do not use while driving. Keep away from water or moisture during storage and use. Use only with the provided adapters. Improper use of this product may result in electric shock.

Maintenance:

Before cleaning, disconnect the device from the power supply. Wipe with a soft cloth and a small amount of neutral cleaning agent. Do not use corrosive cleaning agents. If the product or power adapter is damaged, stop use immediately.

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