BLOW

Elektrosymulator mięśni brzucha i pleców EMS - Abdominal and back muscles massager

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DISCLAIMER

PLEASE DO NOT USE THE EMS OR ANY VIBRATION DEVICE WITHOUT FIRST OBTAINING APPROVAL FROM YOUR DOCTOR IF ANY OF THE FOLLOWING APPLY:

Pregnancy, diabetes with complications such as neuropathy or retinal damage, wear of pace-makers, recent surgery, epilepsy or migraines, herniated disks, spondylolisthesis, spondylolysis, or spondylosis, recent joint replacements or IUD's, metal pins or plates or any concerns about your physical health. Frail individuals and children should be accompanied by an adult when using any vibration device. These contra indications do not mean that you are not able to use a vibration or exercise device but we advise you to consult a doctor first.

Ongoing research is being done to the effects of vibration exercise for specific medical disorders. This will likely result in the shortening of the list of contra indications as shown above. Practical experience has shown that integration of vibration exercise into a treatment plan is advisory in a number of cases. This must be done on the advice of and in accompaniment with a doctor, specialist or physiotherapist.

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USER MANUAL

Attention: Please set a lower intensity when using it for the first time and adjust it after adapting to the low intensity.

1. Please fully charge the massager with the Micro USB cable inside the package. While Charging, the red LED will stay on, when fully charged, red LED will turn off.
2. Press button ③ to turn on the massager, then choose the mode by button ⑤, then you can increase intensity by pressing button ③.
3. Press button ④ to decrease intensity, when reach to 0 intensity, press button ④ again to turn off the massager.
4. There are total 8 different massage modes, by pressing button ⑤ to change mode.
5. Also, you can make all above operation with the remote controller.

Specification

Voltage 3.7V
 Working Mode 8 Modes
 Auto Shut-Down Time 15mins
 Intensity Level 19 Levels
 Output Current 35mA Maximum
 Used for Body Part Neck, Waist, Leg, Arm,
 Power Normal Mobile Chargers
 Battery 200mAh
 Working Time 60-80Mins
 Charging Time 2H
 Accessory: Micro USB Cable*1,

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The product meets the requirements of the directives:

2014/30/EU, 2014/35/EU, 2011/65/EU, (EU)2015/863, (EU)2017/2102

The declaration of conformity is available for download from the website www.prolech.com.pl

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