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Electric Toothbrush MagniSonic 506

zenwire

Thank you for purchasing the Zenwire product. We wish you a pleasant experience!**Product**

Description:

The electric toothbrush is a high-frequency vibration toothbrush. High-frequency vibration effectively removes bacteria and impurities around teeth, and its cleaning range covers various areas around the teeth.

The sound wave energy generated by the motor is transmitted to the tooth surface and gum line through the bristles. This helps remove plaque, loosen tartar, and eliminate stains while improving adhesion to teeth. It also helps destroy bacteria hidden in gum pockets and various types of parasitic bacteria on the tooth surface. Additionally, it can reduce pulp inflammation, prevent gum bleeding, and help prevent gum recession.

Function Description:

1. Pressing the switch within five seconds after turning on the device will switch to the cleaning program change mode. If the button is held for more than 5 seconds, the machine will turn off.
2. After turning on the toothbrush, the function of reminding you to change the cleaning area every 30 seconds and automatic shut-off after 2 minutes of use will be activated.
 - During charging, the battery indicator flashes, and when the battery is fully charged, the light stops flashing.
 - It is important to monitor the battery level and recharge the toothbrush after approximately 10 cycles of use.
 - While charging, ensure the charging interface of the electric toothbrush remains clean.
3. The toothbrush is charged via a USB-C port. To charge the device, remove the protective cover located at the bottom of the toothbrush and then connect the USB-C cable.

Important Notes:

1. Before using the product for the first time, carefully read the user manual and keep it for future reference.
2. If you have undergone oral surgery or dental treatment in the last few months, consult a doctor before use.
3. If gum bleeding occurs after using the product and persists for more than a week, consult a dentist.
4. The product is intended solely for cleaning teeth, gums, and the tongue – do not use it for other purposes. If you experience discomfort or pain during use, stop immediately and contact a doctor.
5. The brush head should be replaced every 3 months or sooner if it becomes deformed or the bristles are worn out.
6. Do not use tips other than those recommended by the manufacturer.
7. Do not allow the battery to fully discharge. Do not store the device with a discharged battery, as this may shorten its lifespan. Charge the battery regularly to ensure optimal performance.
8. The product is not a toy – keep it away from animals and small children without adult supervision.
9. The toothbrush is intended for users aged 5 and older.
10. Always switch off the appliance before cleaning and maintenance, and if it is not working properly and after use.
11. Do not use with damaged tips, as this may lead to injury.
12. Do not leave the appliance switched on unattended.
13. Do not expose the device to extreme temperatures (above 40°C or below 0°C).
14. Store in a dry place, away from direct sunlight.
15. This unit has rechargeable batteries that cannot be replaced.

Important Safety Information:

Before using the product, carefully read this user manual and keep it for future reference.

1. After each use, rinse the toothbrush head and clean it thoroughly.
2. During charging, keep the charging port dry and free from water and dust.
3. Do not clean the charging station and cable with cleaning agents or immerse them in water. Do not charge the toothbrush in excessively high (above 40°C) or low temperatures (below 0°C).
4. If the toothpaste used contains peroxide, sodium hydride, or other strong whitening agents, clean the toothbrush head thoroughly with soapy water after each use, as toothpaste residue may cause plastic to crack.

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⚠ Please pay attention to the following precautions during use:
Do not use chargers larger than 2A.

RoHS