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SPECIFICATION

- Material: high-quality elastic material, covered with cooling gel, ensuring a comfortable fit to the shape of the face and head
- Dimensions: 22x17x0.5 cm (universal size, suitable for most head sizes)
- Color: black
- Cold therapy (for migraines, headaches, stress)
- Heat therapy (for muscle tension, relaxation)
- Benefits: light blocking, soundproofing, reduction of headache, migraine, tension, stress, and swelling
- Safety: safe for use during pregnancy, does not require medication
- Easy to clean: simply wipe with a damp cloth
- Brand: MEDISAN

SAFETY INSTRUCTIONS

1. The product is not a toy – The cooling/heating mask is intended for adult use only and is not a toy. It should be kept out of reach of children to prevent any accidents or improper use.
2. Safe use – The mask can be used for both cold and heat therapy, but before use, ensure that the temperature is appropriate (do not overheat the mask in the microwave).
3. Avoid contact with sharp objects – Although the mask is made of durable materials, avoid contact with sharp edges to prevent surface damage.
4. Do not use on open wounds or damaged skin – The product is designed for use on healthy skin.
5. Use during pregnancy – The mask is safe for pregnant women and provides a natural method of pain relief without the need for medication.
6. Regular inspection – Before each use, check the mask for any damage that may affect its effectiveness.

ASSEMBLY INSTRUCTIONS

1. Preparing the mask: For cold therapy, place the mask in the freezer for at least 1 hour to reach the optimal cooling temperature. For heat therapy, place the mask in the microwave for 10-15 seconds, depending on the microwave's power. Ensure the mask is not too hot before applying it to your face.
2. Wearing the mask: Place the mask on your face, making sure it fits well against the skin and blocks out light. This can help enhance comfort and therapy effectiveness.
3. Using the mask:
4. Cold therapy: Keep the mask on your face for 15-20 minutes to relieve migraines, headaches, or stress.
5. Heat therapy: Apply the mask to your face for 15-20 minutes to relax tense muscles and improve circulation.
6. Cleaning the mask: To maintain hygiene, wipe the mask with a damp cloth. Avoid using strong detergents that could damage the material.
7. Storage: After use, store the mask in a dry and cool place. Avoid contact with sharp objects that could damage the material.

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