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EN: USER'S MANUAL S150

Product features: Electrostimulator of the buttock muscles. Dimensions: 20 x 26 cm. Packaged in a box. The set includes a control panel and a wristband. Panel powered by 2x AAA batteries (not included). High stimulation power. Quick and visible results. High durability.

Service:

Before use, insert the battery: open the battery compartment and insert 2 AAA batteries (not included), remembering to observe the correct polarity, then close the battery compartment. Put the devices on your buttocks and turn them on using the control panel. The device can be used while exercising or walking, but also while cleaning, watching movies, reading books, etc.

Warnings:

Before using the product, read the instructions in their entirety and keep them for the life of the product. The owner is responsible for ensuring that all users of this product are properly informed about how to use it safely. Use the product as intended. This product is intended for home use only. Protect against moisture and dirt. Do not use or store the product in humid conditions. Protect the product from direct sunlight. Avoid contact with water and other liquids. Do not allow the product to heat up or be near sources of fire. Do not disassemble or modify the product. Any repairs may only be performed by an authorized person. The product may only be maintained by an authorized person. Never use a damaged product. If the product does not function as intended, turn it off immediately. If you are not sure whether stimulation will harm your health, consult your doctor before using it. Do not stimulate places with cancerous changes or their vicinity. Do not use the device in the bathtub or shower. Do not use the device while sleeping. Stimulation should not be used while driving, operating machinery or performing other activities that pose a risk of injury. The device should only be attached to normal, undamaged and clean healthy skin. Stimulation should be stopped or the intensity level reduced when you start to feel discomfort. If this feeling does not go away, contact your doctor. Do not use the stimulator near shortwave or microwave therapeutic devices as this may reduce its output power. The stimulator should not be used together with electronic health monitoring devices (e.g. pulse oximeter, EGG alarm unit), because in such a case their readings may be incorrect (disturbed). Skin irritation or burns may occur at the point of contact with the electrodes. This device must not be used if we have an implanted pacemaker, defibrillator or other metal or electric equipment. In this case, there is a risk of electric shock, burns, electrical noise or death. Do not mix old and new batteries. Batteries that are not used for a long time should be removed. Batteries must be inserted with correct polarity. The item is not a toy. Keep away from children. Information for users on disposal of electrical and electronic devices (applies to households) The symbol on the products or the documentation attached to them informs that faulty electrical or electronic devices cannot be disposed of with household waste. The correct course of action when it is necessary to dispose of electrical or electronic devices, disposal of them, reuse them or recover components is to take the device to a specialized collection point, where it will be accepted free of charge. In some countries, you can return the product to your local distributor when purchasing another device. By disposing of your device correctly, you can save valuable resources and avoid negative impacts on health and the environment that may be at risk due to inappropriate waste handling. Details of your nearest collection point can be obtained from your local authority. Improper waste disposal may result in penalties provided for in relevant local regulations. Users in European Union countries: If you need to dispose of electrical or electronic equipment, please contact your nearest point of sale or supplier who will provide further information. Waste disposal in countries outside the European Union: This symbol is only valid in the European Union. If you need to dispose of this product, please contact your local authorities or dealer for instructions on how to properly dispose of it.

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