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EN - Wireless Body Massager

10 intensity levels
12 massage modes

A multifunctional body massager helps reduce stress and muscle tension in the shoulders, arms, back, abdomen, thighs, and feet. The set includes 6 massage heads. You can choose from several intensity levels to adjust the massage strength to your needs. Made from high-quality materials to ensure durability, effectiveness, and maximum results. Designed for easy massage of the back and legs. This device allows you to perform an effective massage in the comfort of your home.

Maintenance

Before cleaning, disconnect the device from the power source. Wipe with a soft cloth and a small amount of neutral detergent. Do not use abrasive cleaners. If the adapter or product is damaged, stop using it immediately.

User Manual:

1. The set includes a power adapter.
2. Charge the product before first use.
3. To start the massager, press the power button.
4. ↑ - Increase speed.
5. ↓ - Decrease speed.
6. Auto - Change program.

The massager automatically turns off after 15 minutes. It is recommended not to exceed 15 minutes of use.

Charging:

1. Connect the charger to the power outlet.
2. The indicator will turn red during charging.
3. Once fully charged, the indicator will turn green, and the device will be ready for use.
4. Disconnect the massager from the power source before use; the device will not operate while charging.

Cleaning:

1. Before cleaning, disconnect the device from the power source and ensure it is turned off.
2. Wipe with a soft cloth and a small amount of neutral detergent.
3. Allow to dry completely before reuse.

Warnings:

Not suitable for use by children or the elderly. Do not use during pregnancy or immediately after childbirth, after surgeries, spinal injuries, heart diseases, rheumatoid arthritis, circulation problems, nerve damage, sensory disorders, or any skin diseases. It is recommended to consult a doctor regarding the safety of using this product. If you experience discomfort while using the product, stop using it immediately. Do not use in combination with other heat-generating devices or warming agents, including warming creams. Do not apply to damaged or irritated skin. Do not use near flammable materials. Do not open the massager while in use - moving parts may cause injury. Do not use while driving. Keep away from water and damp places during storage and use. Do not submerge the device in water. Use only with the adapter provided in the package. Improper use of the product may cause electric shock. Ensure hair is tied back during use; moving parts may cause entanglement. Do not pull the power cord. Do not use if the product shows signs of damage, as this may result in electric shock. If the power cord is damaged, it must be replaced with a new one. If the product emits unusual sounds, smells, or becomes excessively hot, stop using it immediately. When not in use, store the product in a place inaccessible to children. Not suitable for use with pacemakers.

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