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Cechy produktu:

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How to use:

Charge the battery before use. Before first and repeated use, the product should be thoroughly rinsed under running water. Select the appropriate nozzle and attach it by inserting it into the top hole shown in illustration 1. Fill the tank with water and close the cover. The default setting of the new razor is the first gear. Press the "on/off" button to turn on or off. Use the button shown in illustration 4 to select the appropriate operating mode. The first operating mode is suitable for gentle cleaning of the oral cavity. The second mode has a gentle water pressure, suitable for beginners and sensitive teeth. The third gear is a combination of the first and second, it is used for deep cleaning of the oral cavity, not suitable for people with sensitive gums. It is possible that after the first use, the gums will bleed because they are not adapted to this type of cleaning. When the indicator starts flashing, it means that the product is discharging, immediately connect the device to the port and leave it to fully charge.

Warnings:

Use the product only when fully charged, it is forbidden to use the product connected to the power source. Do not change, store and clean the oral cavity at high temperature, to avoid accidents due to scalding. Please keep the machine out of the reach of children. The product should be cleaned with clean water, do not use corrosive agents or detergents. Protect the product from mechanical damage. The product is not a toy. Not suitable for children. Do not expose the product to direct sunlight or harmful conditions. Store the product in a safe place. Information for users on the disposal of electrical and electronic devices (applies to households) The symbol shown on products or in the documentation accompanying them informs that faulty electrical or electronic devices must not be disposed of with household waste. The correct procedure for the disposal of electrical or electronic devices, recycling, reuse or recovery of components is to hand the device over to a specialised collection point, where it will be accepted free of charge. In some countries, the product can be returned to the local distributor when purchasing another device. Correct disposal of the device allows you to save valuable resources and avoid negative impacts on health and the environment, which may be endangered by inappropriate waste handling. Detailed information about the nearest collection point can be obtained from your local authorities. Incorrect waste disposal is subject to penalties provided for in the relevant local regulations. Users in the European Union: If you need to dispose of electrical or electronic equipment, please contact your nearest point of sale or supplier for further information. Disposal in countries outside the European Union: This symbol is only valid in the European Union. If you want to dispose of this product, please contact your local authorities or the retailer for information on the correct procedure

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