medivon.

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PL

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zmiana trybu

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Safety instructions

- Remove the adapter from the socket when not in use.
- Do not place or use the massage in the bathroom or wet / damp places.

- For safety, we equip the massager with protection against overheating. The heat generated in engines is subjected to heavy loads and excessive pressure. Therefore, use the device within the recommended time of 15 minutes. After this period, the massage device will turn off automatically. Resumption of work will take place after a few minutes break.
- We recommend that you do not exceed a continuous massage of 15 minutes. Prolonged use of the device may lead to muscle damage and fatigue.
- The massager meets all technical and safety rules.
- Do not get wet, never remove the cover.
- This item is NOT a TOY. Close adult supervision is imperative when used by children
- Never leave it unattended while connected to the mains.
- Any possible repairs may only be carried out by authorized specialist personnel.
- Improper use and unauthorized repairs are not permitted for safety reasons and may void the warranty.
- Never touch the plug with wet hands.
- Avoid contact of the device with water, high temperatures and direct sunlight.
- Do not use damaged cables, plugs or loose sockets.
- Do not use the device under a cover where the ventilation opening may be blocked.
- If the plugs or cords are damaged, they must be replaced by the manufacturer, service agent or qualified personnel.
- In case of failure, disconnect from the mains immediately.
- Do not use if you have skin problems, open wounds, or swelling
- Improper use or incorrect use excludes all liability for damage.
- Do not use this device while driving.
- Do not use while sleeping.
- Massage should not be performed during pregnancy, in places of body contusions, skin damage, all kinds of inflammation, swelling.
- Do not use if you have a pacemaker or are undergoing treatment for ailments and discomfort.

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The crossed-out wheeled bin symbol indicates that in the European Union after the end of the product's useful life it must be disposed of in a separate, dedicated collection point. Do not dispose of these products together with unsorted general household waste

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English

Hereby, Medivon Sp. z o. o. declares that the device equipment Medivon MD-7280 is in compliance with Directive 2014/35/UE. The full text of the EU declaration of conformity is available at the following internet address:

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TERMS OF WARRANTY

1. The product warranty is valid for a period of 24 months.
2. The warranty covers the period from the date of receipt of the device from Medivon Sp. z o. o. and covers defects resulting from causes inherent in the received device.
3. The device will be repaired as soon as possible, not exceeding 21 working days from the date of acceptance of the device for repair by an Authorized Service Center.
4. The customer has the right to apply for replacement of equipment free from defects, if the Authorized Service finds that the removal of the defect is impossible or requires excessive costs.
5. If only part of the product is defective and can be detached from the part of the product working properly, the Customer's rights resulting from this warranty limit only to the defective part of the product.

In order to benefit from the warranty conditions, please send a complaint notification together with the warranty card to the e-mail address serwis@rgmedia.pl. The exact warranty conditions are under the address:

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I represent that I have read and accepted the terms and conditions specified in the Warranty Certificate.

.....