

Dr.Finger VOLUMEBOOSTER

TDDS-Electroporation Scalp Massager



ROOTONIX



1.1 Dr. Finger Volume Booster Gold Package introduction

Dr. Finger Volume Booster Gold Edition

“Scalp-Hair loss care solution Proposal”



1.2 Dr. Finger Volume Booster Gold Package introduction – Before & After

4 DAYS

Experimental Data

"I felt a change in just 4 days."



1 month

Experimental Data

"It's amazing how much can change in just one month."



Result value used in

Result value used in

1.3 Dr. Finger Volume Booster Gold Package introduction


“Experience everything from scalp diagnosis to management at once!”

STEP 1



FHD high-resolution scalp scanner
Easily measure at home

STEP 2



Get diagnosed with A.I
expert data

STEP 3



Volume booster and ampoule
Start of scalp care solution

2.1 Dr.Finger Volume Booster Point

PATENT of Dr.Finger Volume Booster



ELECTROPORATION



MICRO VIBRATION



LED THERAPY



Voice Guide



Device Manufacturing

Health Care Product Manufacturing – Scalp Care Device



VOLUME BOOSTER

With ROOTONIX's own technology, anyone can easily manage it at home. Data based on real users, including UX/UI for user convenience based design is leading the home care market. This is the strong potential of ROOTONIX

R-TDDS (Rootonix Transfer Drugs Delivery System)
Electroporation

ROOTONIX's 'Electroporation' technology
The world's first patent for drug absorption using micro-current
for scalp care
It is a technology that delivers drugs to the dermis and hair follicles
by the power of electric current rather than needles.
Compared to the injection method or MTS method, it is painless
and is a very effective technique for drug absorption.


Methods such as MTS or injection are no longer needed.



No wounds



No wounds




High Absorbency




2.2 Dr.Finger Volume Booster function

Volume Booster strengthens the hair root volume and make the hair more plentiful.
Use this device comfortably for 10 minutes a day and then you will feel a miracle


Technique 01



Technique 02





Technique 03



R-TDDS(ROOTONIX-Transfer Drug Delivery System)

Through specific conditions and specific microcurrent values of ROOTONIX Technology that helps to reach the dermis layer by shaking the skin lipid film




ROOTONIX Co., Ltd. X Seoul National University/Technical University Berlin Experimental Results (2020~2021)



서울대학교 화학과
교수 이연





베를린 공대 전자공학 교수
Prof. Dr. Manfred Oppel



2.3 Dr.Finger Volume Booster absorbency test

“ Hair loss ampoule absorption rate 300%,
Amazing effects of Dr.Finger Volume Booster ”



<Report : Seoul National University Electroporation Penetration Experiment, 2020>

It will make good ingredients penetrated into follicles **3 times** intensively than applied by hand, so it can help you quickly solve your hair loss or scalp problems.

3.1 Scalp Scanner

“With a FHD camera,
Scan and Diagnose Problems”

Anywhere, Anytime Scan and
diagnose scalp troubles with FHD
precision scanner.

Easily check and care for your
never-before-seen scalp at home



Photo



Video



LED Light




HD



DIRECT

Wi Fi

4.1 monadoo A.I APP SERVICE



“Scalp care, hair loss care
Now with Monadoo.”



Download on the
App Store

GET IT ON
Google Play

4.2 monadoo A.I APP introduction



Scalp type analysis

ROOTONIX's AI analyzes your scalp based on the diagnosis of real experts.



Evaluation by item

Check the score for each item for scalp condition. You can easily find out what you need to focus on.




Expert Suggestions.

Get suggestions from experts on scalp care methods and effective ingredients based on comprehensive diagnosis.



Comprehensive scalp evaluation

We will inform you in detail about the current scalp condition and precautions.



Package introduction

Gold package box



Manual 3ea

(App guide, volume booster manual, scanner manual)



Dr.Finger
Volume Booster

SCALP SCANNER



Product



Dr.Finger Volume Booster

Name	scalp care device
Model Name	RT-01
Registration number	R-R-rio-RT-01
Rated voltage	3.7 V, 600mAh
Made by	ROOTONIX Co., Ltd. Korea
Battery model name	PN702535
Battery authentication number	XU101151-17001A-9
Return and Exchange	the place of purchase
Customer Service	070-7796-0101

CE / IC / PSE / FCC Certification No.
TMC190610102 / UNIA21090701ER-81
KR21-YEF0095 / DREKFCC2005-0132




SCALP SCANNER



Network standard	IEE 802.11 b/g/n
Image Sensor	High performance sensor
Sensor size	1/4
magnification	50-1000times
Lens structure	2G+IR
Camera light source	8 adjustable LED light
Operating temperature	32 to 113 F(0-45 °C)
Battery	800mAh lithium battery
Working time	1H
Power input	DC 5V 500 mAh/1A
Focal length	0~infinity
Camera resolution	1920*1080P
Picture format	JPEG
Video format	MP4



MONADOO



APP Name	MONADOO
Developer	ROOTONIX Co., Ltd.
Volume	59.5MB
language	KOR/ENG



Award & Certification

2023 CES
Innovation Award



Certification
completed



Korea
(KC)



USA
(FCC)



Japan
(PSE)



Australia
(RCM)



Canada
(IC)



Europe
(CE)



Scalp center in my hand

ROOTONIX

ROOTONIX's Volume Booster will be the most important partner for your steady scalp care.

Contact us

ROOTONIX Co., Ltd / CEO Daekwon Jeong (ceo@rootonix.com)

website : www.rootonix.kr Tera Tower2 B-1311 201, Songpa-daero, Songpa-gu, Seoul

2023 Copyright © ROOTONIX Co.,Ltd All right reserved