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CE 0044
Microlife Corporation Taipei - TW

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**EN : PACIFIER THERMOMETER**

In order to use the product properly, please read these notes carefully and keep them for reference if required.

- Always keep an eye on your child when using the lollipop thermometer. Do not leave your child unsupervised when using it.
- This product is not a toy and should not be used as a pacifier.
- The pacifier thermometer is suitable exclusively for oral measurement of the body temperature of children up to three years of age!
- The minimum measurement time until the signalling tone [beep] must be maintained without exception!
- The readings displayed by the pacifier thermometer may be too low under some conditions (e.g. if the child is crying or persistently breathing through his mouth). Therefore, in case of doubt, or if there are signs of fever (reddening of skin, sweating, agitation, unwillingness to drink), please always perform a follow-up rectal measurement if the result is more than 38,5 °C.
- Protect the thermometer from hitting and striking!
- Avoid ambient temperatures above 60 °C. NEVER boil the thermometer!
- The thermometer can/must not be opened! The batteries cannot be changed for safety reasons.
- Use only the commercial disinfectants listed in the section «cleaning and disinfecting» to clean the thermometer. The thermometer must be intact when immersed in liquid disinfectant.

• Warning: The thermometer contains a long-life silver oxide battery (not replaceable).

For disposal, please comply with local regulations for special waste.

- It is not necessary to calibrate the measurement accuracy within the limited 2-years battery life time.

• Ensure that children do not use the instrument unsupervised; some parts are small enough to be swallowed.

• The function of this device may be compromised when used close to strong electromagnetic fields such as mobile phones or radio installations and we recommend a distance of at least 1 m. In cases where you suspect this to be unavoidable, please verify if the device is working properly before use.

TURNING ON THE THERMOMETER

To turn on the thermometer, push the button over the display window: a short signalling tone [beep] signals «thermometer ON». A display test is performed. All segments should be displayed. Then at an ambient temperature of less than 32 °C, an «L» and a flashing «C» appear at the upper right of the display. The thermometer is now ready for use.

USING THE THERMOMETER

The pacifier thermometer is suitable only for oral measurement of the body temperature of children, and it must not be used as a pacifier. Place the thermometer in the child's mouth like a normal baby pacifier. Be sure the child keeps the thermometer in his mouth during the measurement, and keep the child's mouth closed as much as possible. The measurement time signalled is normally about 5 minutes. During a measurement, the current temperature is displayed continuously. The «C» symbol flashes during this period of time.

When the «C» symbol stops flashing, a signalling tone rings 10 times indicating temperature taking is finished. Wait for the alarm to finish before removing the thermometer.

Warning! Certain external factors can lead to a certain measurement inaccuracy in any measurement of body temperature, especially on small children. After a measurement, the maximum temperature measured is shown on the display until the thermometer shuts off. To prolong the life of the battery, turn it off after use by pressing the operating button briefly. Otherwise the thermometer will automatically turn off after about 10 minutes.

STORAGE OF MEASURED VALUES

If you hold the button down for more than 3 seconds when turning the thermometer on, the maximum temperature stored automatically at the last measurement will be displayed. At the same time, a small «M» for memory will appear at the right of the display.

CLEANING AND DISINFECTION

Name: Isopropyl alcohol 70%; immerse: max. 24 hours.

BATTERY EXHAUSTION

When the «▼» symbol (upside-down triangle) appears in the lower right of the display, the battery is exhausted and thermometer needs to be replaced. For battery disposal please refer to «safety instructions».

TECHNICAL SPECIFICATIONS

Type: maximum thermometer

Pacifier: silicone rubber

Measurement range: 32,0 °C to 43,9 °C

Temp. < 32,0 °C: display «L» for low (too low)

Temp. > 43,9 °C: display «H» for high (too high)

Measurement accuracy: ± 0,1 °C between 34 °C and 42 °C at an ambient temperature of 18°C to 28°C

Self-test: automatic internal check at a test value of 37 °C; if there is a deviation of > 0,1 °C, «ERR» (error) is displayed.

Display: liquid crystal display (LCD) with 3 digits

Smallest unit of display: 0,1 °C

Signalling tone: for signalling that the thermometer is ready to use and that the temperature increase is less than 0,1 °C / 32 seconds

Memory: for storing the last measured value

Storage conditions: 25°C +60°C; 15 - 95 % relative maximum humidity

Operating conditions: 0°C -40°C; 15 - 95 % relative maximum humidity

Battery: 1,5 / 1,55 V; battery lifetime: approx. 4500 measurements

Battery replacement: impossible for safety reasons

Expected service life: 2 years or 4500 measurements

Reference to standards: EN12470-3: 2000, clinical thermometers

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US-CFR 16 part 1511, requirements for pacifier.

ASTM F963-96a, Specifications on toy safety

This device complies with the requirements of the Medical Device Directive 93/42/EEC. Technical alterations reserved.

GUARANTEE

This product is warranted for a two-year period from the purchase date. Any damage causes by improper handling shall not be covered by the guarantee. Batteries and packaging are also excluded from the guarantee. All other damage claims excluded. Please report any accidents to Dorel at the address given here and ensure that you enclose proof of purchase.

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