53794

AET-R1B1

ES	TERMÓMETRO INFRARROJO .....	1
EN	INFRARED THERMOMETER .....	7
FR	THERMOMÈTRE INFRAROUGE .....	13
IT	TERMOMETRO A INFRAROSSI .....	16
PT	TERMÓMETRO DE INFRAVERMELHOS .....	22

# ES MANUAL DE INSTRUCCIONES

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# EN INSTRUCTION MANUAL

## PRODUCT'S PRESENTATION

This Infrared Thermometer is a non-sterile, reusable, handheld device. It can be used by consumers in homecare environment and doctors in clinics in order to measure the human body temperature of people aged three months and above by detecting infrared heat from the center of the forehead.

- The clinical thermometer is an adjusted mode clinical thermometer.
- The patient can use and maintenance the device as an operator.
- For safe and proper use of this product, be sure to read and fully understand Safety Precautions contained in this user manual.
- Keep this manual with you for a convenient reference.
- Keep this manual in a proper way and avoid any loss.
- If you need other information, please contact the manufacturer.


If a reading below 32°C is recorded a "Lo" sign will be displayed followed by 2 consecutive beeps.

If a reading above 37.8°C is recorded there will be 6 consecutive beeps.


If a reading above 42.2°C is recorded a "Hi" sign will be displayed followed by 2 consecutive beeps.

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1	Probe
2	ON/OFF/MEASURE switch
3	Unit Body
4	LED screen
5	Battery cover


## CORRECT MEASUREMENT METHOD


The correct measurement method to follow is to measure the forehead at a distance from 15mm to 50mm:

## INFORMACIÓN EN PANTALLA

1	Location indication
2	Temperature unit
3	Temperature data
4	Battery symbol
5	Memory symbol


## Instructions for use

1. Turn the thermometer on by pressing the "ON/OFF/MEASURE" button. To choose mode, press "ON/OFF/MEASURE" on the shutoff mode. A beep sound will be heard and the screen will flash with the symbol ---°C, indicating the unit is ready for use.

2. Place the thermometer probe point near the forehead and press "ON/OFF/MEASURE" and release after ~ 1 second. The thermometer will beep and display the accurate temperature. Whenever consecutive readings are required, please wait for the ---°C to flash. The above procedure can be repeated as many times as required.

3. To turn off the device: Press and hold the "ON/OFF/MEASURE" switch for ~ 5 seconds. The thermometer will turn off automatically after ~ 60 seconds.

- Location indication: the icon flashes to indicate the correct measuring position when measuring.

- Temperature data indication: the temperature reading is displayed after measuring.

- Temperature unit indication: unit of measuring temperature shows °C or °F.


- Low battery symbol indication: the icon prompts you to replace the battery when the power is low.

- Memory symbol indication: it shows in memory query mode currently.

## LOADING AND REPLACING BATTERIES

When the battery icon "  " is flashing, replace a new battery immediately:

1. Press and hold the battery cover.
2. Push it upward diagonally to push the cover away.
3. Install the 2\*AAA batteries as shown in the positive and negative pole.

## ⚠ WARNING!

Take out the battery if you do not use the thermometer for a long time.

**BEFORE TAKING A MEASUREMENT**

1. Try the device by yourself first to know how to use it.
2. It must be placed in a general environment of room temperature of 16°C-35°C for at least 30 minutes before use to avoid incorrect measurement results.
3. The person measured should be kept at room temperature for at least 20 minutes to make body temperature balanced.
4. After strenuous exercise, you should rest for at least 30 minutes before taking measurement.
5. Wipe the skin dry and push aside hair strands before commencing.
6. When repeated measurements are made continuously, move away the probe between each measurement and rest for at least 5s before making the next to obtain the most correct result.
7. The normal reading value is about 35.5°C 37.8°C. 8. The user should take more temperature to know his or her own health temperature. Pay attention to physical changes or consult your doctor if there is any abnormality.
8. Check the probe is clean before each measurement to ensure accuracy.

**RECALLING READINGS FROM MEMORY**

1. Ensure the thermometer is switched off.
2. Press the "ON/OFF/MEASURE" button for 4 seconds to put the device on memory mode. The last reading will be displayed.
3. Press and release the "ON/OFF/MEASURE" button again to display the next reading in a sequence from the last reading to the first reading.
4. The thermometer has memory to recall the last 32 readings.
5. The thermometer will automatically turn off if it idles for more than 12 seconds.

**CHANGING TEMPERATURE SCALES**

1. Ensure the thermometer is switched off.
2. Press the "ON/OFF/MEASURE" button for 8 seconds to bring it into temperature mode. The current temperature unit will be displayed.
3. Release and press the ON/OFF/MEASURE button again to toggle and select between °C and °F.
4. The thermometer will switch off automatically if it idles for more than 4 seconds and will display the last reading

**WARNING!**

- Do not use the thermometer under extreme temperatures (below 16°C/60.8°F or over 35°C/95°F) or humidity extremes (below 15%RH or over 80%RH).
- Do not use the thermometer under extreme temperatures (below -20°C/-4°F or over 55°C/131°F) or humidity extremes (below 15%RH or over 93%RH).
- If temperature probe has been shocked and lux, contact your retailer immediately.
- It is dangerous for patients to perform a self-evaluation and self-treatment based on the measuring results. Be sure to follow doctor's instruction.
- Do not touch or blow infrared sensor.
- Clean a polluted infrared sensor with a soft dry cloth in a gentle manner. Cleaning with toilet tissue or paper towel may scratch the infrared sensor, causing inaccuracy.
- Install the battery in the right position according to the

manual.

- Do not immerse the thermometer in liquids, this device is not waterproof.
- Do not use a mobile phone nearby when the thermometer is working.
- If the probe glass is damaged, stop using it and contact your retailer immediately.
- Do not modify the device.
- The device is not repairable and contains no user serviceable parts.
- The device requires no calibration.
- If the device is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.
- The operator shall not touch the thermometer and the patient simultaneously.
- Keep the device out of children's reach.
- For temperature difference between storage area and measurement site, condition the thermometer for about 30 minutes in room temperature.
- Ensure that children do not use the instrument unsupervised. Some parts are small enough to be swallowed.
- No servicing/maintenance white the MW equipment in use.
- The patient can measure data and change battery under normal circumstances,
- The patient can maintain the device according to the user manual.

**ERROR MESSAGES**

Error message	Problem	Solution
Hi	Temperature is higher than 42.2°C (108°F).	Operate the thermometer only between the specified temperature ranges. In the event of a repeated error message, contact your retailer.
Lo	Temperature is lower than 32 °C (89.6 °F).	Operate the thermometer only between the specified temperature ranges. In the event of a repeated error message, contact your retailer.
Err	Sensor not welded.	Contact your retailer.
	Low power.	Replace new batteries.
ErH	The temperature of the measuring environment is too high.	Decrease the temperature of the environment, keep it within 16°C-35°C.
ErL	The temperature of the measuring environment is too low.	Increase the temperature of the environment, keep it within 16°C-35°C.

**ELECTROMAGNETIC COMPATIBILITY**

1. This product needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile radio frequency (RF) communications equipment
2. This machine must not be used adjacent to or stacked with other equipment, and if an adjacent or stacked unit is necessary, this device should be checked to verify normal operation in the configuration to be used.
3. Caution: This unit has been thoroughly tested and inspected for its proper performance and operation.
4. Carbon: This machine must not be used adjacent to or stacked with other equipment, and if an adjacent or stacked unit is necessary, this device should be checked to verify normal operation in the configuration to be used.

Guidance and manufacturer's declaration - electromagnetic immunity - for all EQUIPMENT and SYSTEMS

1	Guidance and manufacturer's declaration-electromagnetic emission		
2	The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such an environment.		
3	Emissions test	Compliance	Electromagnetic environment-guidance
4	RF emissions CISPR11	Group 1	The product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
5	RF emissions CISPR11	Class B	The product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
6	Harmonic emissions IEC 61000-3-2	N/A	
7	Voltage fluctuations / flicker emissions IEC 61000-3-3	N/A	

Guidance and manufacturer's declaration - electromagnetic immunity			
The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient / burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% U <sub>T</sub> (>95% dip in U <sub>T</sub> for 0.5 cycle) 40% U <sub>T</sub> (60% dip in U <sub>T</sub> for 5 cycles) 70 % U <sub>T</sub> (30 % dip in U <sub>T</sub> for 25 cycles) < 5 % U <sub>T</sub> (>95 % dip in U <sub>T</sub> for 5 sec	N/A	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Non-contact infrared thermometer requires continued operation during power mains interruptions, it is recommended that the Non-contact infrared thermometer be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE UT is the a. c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration - electromagnetic immunity - for EQUIPMENT and SYSTEM that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM - for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration - electromagnetic immunity			
The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF  IEC 61000-4-6	3 Vrms  150 kHz to 80 MHz	N/A	<p>Portable and mobile RF communications equipment should be used no closer to any part of the product, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$ <p>150 kHz to 80 MHz</p> $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$ <p>800 MHz to 2.5 GHz</p> <p>where <math>p</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in metres (m).<sup>b</sup></p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>
Radiated RF  IEC 61000-4-3	3 V/m  80 MHz to 2.5 GHz	3 V/m  80 MHz to 2.5 GHz	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Non-contact infrared thermometer is used exceeds the applicable RF compliance level above, the Non-contact infrared thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Non-contact infrared thermometer.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Non-contact infrared thermometer			
The product is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the product can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the product as recommended below, according to the maximum output power of the communications equipment			
Rated maximum output of transmitter  W	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz  $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz  $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz  $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	N/A	0.12	0.23
0.1	N/A	0.38	0.73
1	N/A	1.2	2.3
10	N/A	3.8	7.3
100	N/A	12	23

For transmitters rated at a maximum output power not listed above the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## EXPLANATION OF MARKS OR SYMBOLS

The following symbols may appear on the manual, Infrared Thermometer and its accessories.

	<b>MUST OBSERVE</b> Means Obligatory with detailed items expressed in words or figures within or beside the mark. Left one means General Compulsory.
	Follow instructions for use.
	<b>IMPLICATION OF SYMBOL</b> Type-BF applied part.
	Caution: Consult accompanying documents.
	Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.
	Transport package shall be kept away from rain.
	Transport package shall not be exposed to sunlight.
	Indicates correct upright position of the transport package.
	Contents of the transport package are fragile therefore it shall be handled with care.
	Indicates temperature limits within the transport package shall be stored and handled.
	Non-ionizing electromagnetic radiation.
	Do not roll.
	Lot number.
	Production date.
	The device should not be used after the end of the shown or the day.
<b>IP22</b>	Protected against solid foreign objects of 12.5mm diameter and greater. protected against vertically falling water drops when the device is tilted up to 15°.
	Manufacturer.
	CE Mark: conforms to essential requirements of the Medical Device Directive 93/42/EEC.
	Authorized Representative in the European Community.

## MAINTENANCE, STORAGE AND CALIBRATION

### Maintenance

1. Remove any stains on the body with a soft and dry cloth.
2. Clean the probe tip as follows: very gently wipe the surface with a cotton swab or soft cloth moistened with clinical alcohol, after the alcohol has completely dried out, you can use the thermometer. If the probe tip is damaged, please contact the maintenance center.
3. Don't wash with water or detergent containing abrasive or benzene.
4. Do not immerse in liquids.

### Storage

1. The device must not be stored or used at an excessively high or low temperature or humidity, in sunlight, in association with an electrical current or in dusty locations.
2. Remove the battery if non-use for a long time.

### Calibration

This thermometer is calibrated at the time of manufactured. If the thermometer is used according to the instruction, periodic recalibration is not required. If at any time you question the accuracy of the measurement, contact the retailer immediately. Don't attempt to modify or reassemble the thermometer.

### WARRANTY

1 year warranty from the date of purchase, excluding errors made by the user below:

1. Error produced by unauthorized disassembly and modification.
2. Error produced by a fall during use or transport.
3. Error produced by use far from the instructions described in the manual.

### INCLUDED IN DELIVERY

- Infrared thermometer
- Manual

### STANDARD

This product is made under the ISO80601-2-56.

**TECHNICAL SPECIFICATIONS**

Measurement mode . . . . .	Forehead mode
Measurement units . . . . .	Celsius (°C) and Fahrenheit (°F)
Operating conditions . . . . .	16°C-35°C (60.8°F to 95 °F) with a relative humidity of 15% - 80%
Storage conditions . . . . .	-20°C - +55°C (-4°F to 131°F) with a relative humidity of 15% - 93%
Measurement distance . . . . .	15mm – 50mm
Forehead temperature measurement range and accuracy:	
Measurement range . . . . .	32.0°C- 42.2°C (89.6°F - 108°F)
Measurement accuracy . . . . .	±0.2°C (35°C - 42°C)
Measurement accuracy . . . . .	±0.4°F (95°F – 107.6°F)
Measurement accuracy . . . . .	±0.3°C (±0.5°F) other range
Memories recall . . . . .	32 measurement recalls
Dimensions . . . . .	144x101x36mm
Weight . . . . .	81g (without batteries)
Battery . . . . .	2xAAA batteries (DC 3V) (not included)
Switch off . . . . .	Automatically switches off after 1 minute

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Alicn Medical Shenzhen, Inc.
Shanghai International Holding Corp. GmbH (Europe)