BLACKROLL®

FASCIA TOOLS

Instruction Manual

EN | DE | FR | NL | IT | PL | CN

Product description

Congratulations on your new BLACKROLL® fascia product! These fascia rolls and balls made of EPP (expanded polypropylene) are perfect for self-massage and fascia training. They promote blood circulation, release tension, relieve pain and improve mobility. The products are available in different sizes and surface structures and can be used for different parts of the body as needed.

Sizes and degrees of hardness

Roll length:	6 cm, 15 cm, 30 cm, 45 cm
Roll diameter:	3 cm, 5 cm, 10 cm, 15 cm
Degrees of hardness:	soft (80%), medium (100%), hard (150%)
Ball diameter:	8 cm, 12 cm

Surfaces

Smooth:	all-rounder for regeneration and activation.
Saw blade:	textured surface for increased blood flow.
Grooved:	activating effect with intense increase in blood flow.

Safety instructions

- + Use the fascia products on a non-slip surface to avoid slipping or falling.
- + Pay attention to the controlled application of pressure: Excessive pressure can cause bruising or tissue damage.
- + Consult a doctor if you suffer from health problems such as osteoporosis, blood clotting disorders, herniated discs or varicose veins.
- + Avoid application to acute injuries, redness or inflammation.
- + Store the products in a dry place and not in direct sunlight.

Contraindications

Do not use the BLACKROLL® fascia rolls and balls if you have osteoporosis, blood clotting disorders, fresh injuries or inflamed skin areas without consulting your doctor first. The products should not be used on areas with acute pain or inflammation.

Application

1 Preparation

- + Make sure that the surface on which you are using the fascia products is stable and non-slip.
- + Start with light pressure and increase it if necessary.

2 Use

- + Slowly roll the fascia roll or fascia ball over the affected muscle area.
- + Make sure to adjust the intensity of the massage to your well-being.
- + Smaller rolls and balls are better suited for targeted treatments of small muscle groups (e.g. feet, arms).

3 After use

- + Clean the products with water and a soft cloth.
- + Store them in a dry place to ensure the longevity of the material.

Technical data

Material:	EPP (expanded polypropylene)
Hardness:	Varies by product (soft, medium, hard)
Insoluble in water:	Yes
Recyclable:	Yes

Care instructions

- + Clean the fascia products regularly with water and a soft cloth.
- + Avoid aggressive cleaning agents or chemical substances that could damage the material.
- + Store the products in a dry place and protected from direct sunlight.

Possible uses

- + Self-massage: Use the fascia rolls and balls for targeted massage of muscle groups such as the back, neck, thighs or calves.
- + Pain relief: They help to release muscle tension and improve mobility.
- + Fascia training: With regular use, the products can help to prevent muscle soreness and support recovery.

Warranty and liability

The BLACKROLL® fascia products come with a two-year warranty. Damage caused by improper use or failure to follow the care instructions is excluded from the warranty.

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BLACKROLL AG
WSVK Oederan GmbH