SPIS TREŚCI

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DIGITAL BLOOD PRESSURE MONITOR

MODEL: TMA-10

Thank you for purchasing Digital Blood Pressure Monitor TMA-10.

TECH-MED
TECHNICA MEDICINA has always been with you for over 75 years – we use our experience to continually improve our products, which incorporate all the latest technologies in diagnostics and are constantly subjected to quality control. Your choice is the best possible one. We wish you good health. We are always at your disposal.

INTRODUCTION

The TMA-10 blood pressure monitor is a fully automatic digital blood pressure measuring device with advanced technology. Pressing the button after the cuff is in place causes the device to measure blood pressure and heart rate while inflating the cuff, and then display the results on the digital display. TMA-10 also detects symptoms of arrhythmia. After the symptom repeats several times and the «♥» sign is displayed on the LCD panel, it is advised to consult the doctor. Do not interpret neither the results of blood pressure nor the signs of arrhythmia by yourself. These results should be interpreted by a PHYSICIAN or a qualified medical professional who knows the patient's medical history. Taking regular measurements and recording the results allows the doctor to assess the trends in the patient's blood pressure over a longer period of time.


What you need to know about blood pressure

What is blood pressure?

Blood pressure is the force with which flowing blood acts on the walls of the arteries. Blood pressure changes during the cardiac cycle. The highest pressure in the cycle is called **SYSTOLIC BLOOD PRESSURE**. The lowest is called **DIASTOLIC BLOOD PRESSURE**. Measuring both is essential for the doctor to assess the patient's blood pressure status. Blood pressure is influenced by many factors, such as physical activity, anxiety, or time of day. It changes constantly throughout the day.

It rises rapidly in the early morning and decreases before midday, then rises again in the afternoon and finally falls to a low level at night.

It can also change in short intervals. Consequently, the results of subsequent measurements may vary. The graph below illustrates the changes in pressure over the course of a day with measurements taken every 5 minutes. Thick line represents sleep. Pressure spikes at 16 (A) and 24 (B) correspond to an attack of pain and sexual activity.

Measuring blood pressure in the doctor's surgery only gives an instantaneous value. Regularly repeated measurements at home allow the actual blood pressure to be determined in the patient's everyday environment. In addition, patients often have a different pressure measured at home because they are more relaxed than in the doctor's surgery. Regular home measurements give the doctor valuable information about the patient's blood pressure under normal conditions.

What is an arrhythmia?

The heart acts as a pump that pumps blood through its four chambers. Blood is pumped through muscle contractions occurring in a strictly controlled sequence. This process is controlled by bursts of eells that control the electric activity of the heart. If this sequence is disrupted,

heart rhythm disturbances (arrhythmias) develop. They cause a reduction in the efficiency of blood pumping. Most arrhythmias are transient and benign. Most often, the heart misses one beat every so often or extra beats occur. Such occasional disturbances may be caused by strong emotions or intense physical activity. However, some types of arrhythmias can be life-threatening and require treatment.

Symptoms of arrhythmia

General symptoms of arrhythmia: palpitations or a feeling of a sudden strong heartbeat, a feeling of fatigue or dizziness, loss of consciousness, breathlessness, and chest pain.

Can arrhythmias be treated?

The treatment of an arrhythmia depends on the type of arrhythmia and on the patient's age and physical condition. There are methods to prevent arrhythmias. These involve using relaxation techniques to reduce stress and limiting caffeine, nicotine, alcohol, and stimulants. Many forms of arrhythmia do not require treatment. They are naturally eliminated by the body's immune system.

However, other forms of arrhythmia should be controlled. For this purpose, cardiac medications, implantable automatic defibrillators, or artificial heart pacemakers

Cardiac arrhythmias can cause very serious consequences, e.g. in the United States they cause almost 250 000 deaths a year. Arrhythmias that start in the ventricles are more serious than those that start in the atria.

SAFETY INFORMATION

The signs below might be in the user manual, labeling or other component. They are the requirement of standard and using.

Caution: These notes must be observed to prevent any damage to the device.

Symbol for „THE OPERATION GUIDE MUST BE READ“

Symbol for „TYPEBAPPLIED PARTS“

Symbol for „COMPLIESWITH MDD 93/42/EEC REQUIREMENTS“

Symbol for „ENVIRONMENT PROTECTION -

Electrical waste products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice“

Symbol for „MANUFACTURER“

Symbol for „DIRECT CURRENT“

Symbol for „SERIAL NUMBER“

THE PRODUCT FEATURES

1. The device has the WHO warning strip „■” symbol display function, with the blood pressure value change, the higher blood pressure value „■” symbol will appear relatively in the higher position (refer to WHO definition of high blood pressure levels-reference table).
2. The device has the average of the last three measurements function (press the memory key for the first time to display the value that is the average of the last three measurements).
3. The device has a large screen display and digital legible.
4. The device has kPa and mmHg display switching and measurement functions.
5. The device has double memory lookup function and can store 99 set of measurements for each person to remember the normal status of your blood pressure.
6. Clock function: year, month, day, hour, minute setting.
7. Low power detection: detecting low power under any working state, LCD displays „□” symbol indicates low power.
8. Overpressure protection function: when the pressure is more than 300mmHg, the device will exhaust fast automatically.
9. Automatic shutdown function: the device will shut down automatically after 1 minute of no operation.
10. Heartbeat prompting voice.
11. Measurement completed prompting voice.
12. Incorrect prompting function.

SPECIFICATIONS

Name	Full-automatic arm blood pressure monitor
Display mode	LCD Digital Display
Measurement method	Oscillographic measurement method
Measurement range	Pressure:0~280mmHg (0kPa-37.3kPa), pulse: 40~199 beats /min
Accuracy	Pressure: Within ± 3 mmHg, pulse: Within $\pm 5\%$
Power supply	d.c. 6V (4 AA batteries) or USBC Power Supply – not included

Operating environment	5 ~40°C , 15%RH~80%RH
Storage temperature	-20~55°C , 10%~93%RH
Measuring range	About 220mm~320mm
Weight	About 450g
Dimensions	About 98mm*126mm*60mm
Security classification	Internal power supply type B
Pressure mode	Automatic pressure mode
Exhaust way	Automatic rapid exhaust
Pressure testing	Semiconductor pressure sensor
Pulse testing	Semiconductor pressure sensor
Note: If the device is used outside the specified temperature and humidity ranges, the system may not be able to reach the above-mentioned performance specifications. This device is not intended to be used by children and it applies only to adults. This product has got the product registration certificate, in case of doubt about the effectiveness of the product, see its certificate on the company's website.	

SECURITY CONSIDERATIONS

Please read this manual carefully before using this product.

The warning signs and illustrations shown on the manual are for the purpose of using the product safely and properly as well as to avoid any damage to you and others.

The warning signs, illustrations and related meaning are as below.

Warning signs

NOTE

Do not use the device if an error occurs. Otherwise a personal injury may happen.

Items damage means the damage about houses, possessions and livestock, pets.

ILLUSTRATIONS

The sign represents prohibited.

The specific prohibited contents are in the sign  or marked on the side with words or pictures. The sign on the left represents "General prohibition".

The sign represents enforcement.

The specific prohibited contents are in the sign  or marked on the side with words or pictures. The sign on the left represents "General enforcement".

The sign means type B applied part.

The sign mean class II equipment.

When discarding batteries or the product must comply with local regulations.

Please read instruction manual before using this product.

Electronic information products pollution control flags.

This product environmental protection use period are 10 years, and can be recycled and should not be discarded. Please read instruction manual before using this product.


NOTE

-  It is dangerous for the patients to self-treat and self-judge according to the measurement result, so please follow the doctor's instructions.
-  Self-judgment may accelerate the progression of the disease.
-  If the pressure does not stop, please press the on/off key
-  If do not treat in time, there may be a risk.
-  The person who is suffering from blood circulatory disorders, please use the device under a doctor's guidance.
-  If the user is suffering from an irregular heartbeat, the user needs to consult professional doctors to confirm the measure blood pressure values.


- ❶ When the common arrhythmia such as atrial premature, premature ventricular and atrial fibrillation appears, the measured value may not be accurate or blood pressure cannot be measured.
- ❶ Do not install the wrong battery polarity.
- ❶ After the battery is exhausted, please instal four new batteries.
- ❶ If the device is not used for more than 3 months, please remove the battery as it may cause leakage, overheat, rupture or damage.
- ❶ If you get the device from less than 5°C, please do not use immediately and place in a warm environment at least 1 hour.
- ❶ Otherwise it cannot measure correctly.
- ⊗ This product is only used to measure blood pressure, otherwise accidents may occur.
- ⊗ There are sensitive electronic components in this blood pressure monitor, so please do not use the device under the environment of strong electromagnetic interference, such as the nearby of mobile power, microwave ovens, etc.
- ⊗ If you find the cuff is damaged, please do not measure the blood pressure. Otherwise it may damage skin. The cuff accessories shall not be changed casually.
- ⊗ Please do not use the device in vehicles such as cars, planes, otherwise it can't measure correctly.
- Ⓜ Please do not disassemble or repair the device.
- ⚠ Please use four AA alkaline batteries, do not use other batteries, otherwise it may cause a fire.
- ⚠ Be sure to use the product-specific power adapter.

COMPOSITION OF THE PRODUCT AND COMPONENT NAMES

CONSIDERATIONS FOR MEASURING

Measure the blood pressure at the same time every day, using the same arm and the same posture. Rest for 4-5 minutes before taking the measurement.

1. The correct using method
 - 1 Elbow on the table.
 - 2 Keep the center of the cuff and the heart or nipples at the same height.
 - 3 Do not put the clothes into the cuff.
 - 4 Palms up, keep body relaxed.
 - 5 Sit in a chair, keep your feet flat on the ground.
2. Ideal environment for measuring blood pressure
 - 1 Get up in the morning, when the mood is relaxed.
 - 2 No consciousness to the toilet
 - 3 Room temperature is about 20°C
 - 4 A quiet place, no noise around.
3. Measure the blood pressure at the same time every day. Blood pressure is constantly changing and cannot be judged after a single measurement. Repeated measurements are more reliable. Choose a specific time of a day to measure the blood pressure, for example in the morning after getting up.

Note: be sure to have more than 2-3-minute interval between two measurements. Depending on your individual physiological characteristics, between the two measurements may need to rest for a longer time.

Please keep the wristband around the wrist before you start measuring.

Do not measure your blood pressure on bus.

CONNECT THE POWER SUPPLY

This machine can use an external power adapter or alkaline batteries. Or both can be used simultaneously. Power Supply is not included.

1. Use an external power supply

First connect the product with the DC power adapter and then plug in AC220V power source. After use, please unplug the power source connected AC220V and then unplug the DC connected to the product.

2. Use the batteries.

Please be sure to use alkaline batteries

- 1 Turn host over and remove the battery cover in the direction of arrow.
- 2 After you confirm the battery positive, install four AA5 alkaline batteries.
- 3 Replace the battery cover.

Replace the batteries

Replace the batteries in the following cases:

The screen appears low battery „“symbols during measured.

Voice broadcast is slower.

The screen does not display anything after pressing the on/off key.

Please deal with the waste dry battery after use in accordance with the relevant provisions of city environmental protection.

CUFF STRAPPED METHOD

1. Tie the cuff to the bare arm, the export of rubber hose should be aligned with the Palm side of the little finger. Do not put the clothes into the arm band, otherwise it will affect the measurement accuracy.
2. The cuff can't be tied too tight nor too loose, it is most appropriate to put a finger.
3. Put the elbows on the table, palms up, relax the body, keep the cuff and heart at the same height when measuring.

If the cuff is tied and positioned incorrectly, the device may not measure the blood pressure correctly.

The right arm and left arm blood pressure may be different, so the measurement of the blood pressure may also be different. We recommend that you always use the same arm to measure. If the measurements of two arms vary greatly, please use the device under a doctor's guidance.

The cuff is only used to measure blood pressure when strapped to the arm.

1. Do not bend the cuff and the air pipe of the cuff.
2. When removing the air tube, please hold the head of the air tube to unplug.
3. Do not pressurize before wrapping the cuff around the upper arm.

FUNCTION KEYS INSTRUCTION


1. **ON/OFF key [START/STOP]**
 - 1) In sleep mode, short press to enter the automatic measuring mode.
 - 2) In other mode, short press for the off key.
2. **MEMORY key [PAMIĘĆ]**
 - 1) In sleep mode, short press to enter memory query mode and for the turn key.

- 2) In the time setting mode, short press for the switching key.
- 3) In the unit selection mode, short press to switch between mmHg and kPa. (Usually use mmHg as the unit of blood pressure).
3. **SETTING key [USTAW]**
 - 1) In sleep mode, short press to view the current time and the current user.
 - 2) In sleep mode, long press for 3 seconds to enter the time setting mode and short press for mode switch.
4. **Set the user time and unit**
 - 1 Set the user
In sleep mode, press the set key until the symbol „“ flashes up, press the memory key to switch „“ or „“, if don't choose the system default user „“.
5. **Set the time**
In user mode, press the set key to switch to set the year or screen displays time mode automatically after the power supplied ; Press the memory key to adjust, press key confirm and switch modes.
6. **Set the unit**
In time mode, press the memory key the screen enter into “PA” mode, “OFF” means the blood pressure unit is mmHg, if don't switch the system default the blood pressure unit is mmHg.

MEASURE BLOOD PRESSURE

When measuring, you should relax and sit properly, the room temperature should be appropriate. Do not eat, drink, smoke, exercise and bath within 30 minutes and it is better to rest for 15 minutes before the measurement, do not move and speak during the measurement.

Application: It is used to measure human blood pressure and pulse.

1. Press the on / off key to display the following screen
 - 1 The pressure begins to rise from 0.
 - 2 The pressure began to decline after the rise, and measured. After the pulse is detected, „♥“ heart symbol flashes.

Measurement end

If you have an irregular heartbeat, the screen will display „♥“

3. If it cannot be measured properly, an error symbol will be displayed and measured values cannot be stored.

If it cannot be measured properly:

When an error symbol is displayed, see symbol description. To continue measuring, press the on / off key to temporarily shut down. Be sure to wait for 4-5 minutes before measurement.

If you feel discomfort in the course of using, for example, when the air bag in a state of lasting inflatable, poor blood flow leads to dizziness, you must terminate the measurement of blood pressure at this time, you can press the on/off key and the device will immediately release the air inside the air bag.

4. Irregular heartbeat display

When the screen displays „♥“ at the end of measurement, it indicates Irregular heartbeat.

5. Press the on/off key to shut down the power supply.

The blood pressure monitor will turn off automatically after about 1 minute of non use. What is irregular heartbeat?

When this device measured systolic blood pressure and diastolic blood pressure, detected the Instant heartbeat rhythm more than upper and lower 25% scope of average heart beat rhythm. In the process of measuring, if the irregular heartbeat is detected more than twice, when the measurement results are displayed, the screen will appear irregular heartbeat icon „♥“.

It needs medical examination and diagnosis to judge whether the icon „♥“ flash means irregular heartbeat.

If irregular heartbeat icon „♥“ often appears, it is necessary to go to the hospital to be checked. Self-judgment and self-treatment according to the measurement results are dangerous, please be sure to follow the doctor's advice.

CHECK AND DELETE THE MEMORY VALUE

The device can store 2x99 sets of memory the device can save the measurement data automatically after each measurement ,once the memory capacity is full, the old measurement data will be overwritten by new data. The device can display the average of last 3 measurements.

- A. Check the average
In sleep mode, press the memory key the screen will display the average of last three measurements. The bottom left of the screen will display „AVG“.
- B. Check the result
In average mode, press the memory key, you can view the last measured value (such as the first group of memory). Then press the memory key to view the second group of memory, repeat this operation continue to view the third group ,the fourth group.....and so on.
- C. Delete memory
 - 1. Delete a single set of memory
Press the memory key to select the number of the memory group you want to delete.
Long press set key until the screen display „d EL“.
Press the memory key again remove the current memory number of group.
 - 2. Delete all the memory
In sleep mode, press and hold the memory key for 3 seconds, and then press the set key and hold until the screen displays ALL EE, all the memory will be deleted.

MEASUREMENT Q&A

- Q. Why is the blood pressure value measured at home lower than in hospital?
- A. People often feel nervous when blood pressure measured in hospital. Because mood is stable when measured at home, sometimes the values are lower than the values measured in hospital about 20mmHg~30mmHg (2.7kPa~4.0kPa).
Know your blood pressure at home when calm is very important.
- A. If the cuff is positioned higher than the heart, the blood pressure measurement value is low. Please use correct posture measuring.

Q. Why is the blood pressure value measured at home higher than in hospital?

A. Whether taking blood pressure medicine?

People taking blood pressure medicine, lose its efficacy can cause high blood pressure values. Please accept the guidance of a doctor.

A. Is the cuff wrapped around too loose?

The cuff is wrapped around too loose will make the pressure cannot transfer to arteries, so the measured blood pressure values higher than the actual values. Don't have a gap between arm and cuff, please roll cuff tightly.

A. Is the measuring posture correct?

The measurement of posture such as bend over, sit cross-legged and bent body on the low table or sofa can produce abdominal pressure and so on and elevate blood pressure.

Why is the blood pressure value measured is not the same each time?

If it is measured at the same time

In order to properly manage the blood pressure, please measure at the same time every day.

Q. Which factors can influence the fluctuation of the blood pressure monitor?

A. Even measured at home, the measurement results under the following conditions will change. Within an hour after a meal

- After drinking, drinking coffee, drinking tea
- After smoking After the bath After the campaign
- After the toilet
- Speech during the measurement fretful caused by stress and anxiety Sudden temperature changes

Measurement site or the environment is not the same with the past

When congestion is caused continuous measurement? Can the continuous measurement cause congestion?

Situation due to the arm is oppressed when measuring to make fingertips poor blood circulation and cause congestion.

When congestion, loosen the cuff, hands over your head, repeated about 15 times the palm grip and stretching exercises, you can clear the congestion.

Why do cuff wrapped around the arm can feel pain and numbness?

This is a temporary phenomenon, do not worry about it.


When measuring blood pressure, we need to tighten the cuff to the arterial blood flow is temporarily stopped, so you may feel some pain and numbness in the arm. After removing the cuff, a little rest can be alleviated.

ABNORMAL STATE OF SELF-TEST

Symptom	Processing Method
Measurement values too high or too low	• When measuring, the location of the cuff relative to heart is too high or too low. • Cuff is not properly tied. • Speech or body shaking when measuring. • Keep the cuff and heart at the same height, correct measurement position.
„Err1“	Because the speed of air leakage is too fast, or the pulse signal is too weak. Please check the cuff, tie up and measure again.
„Err2“	Unable to measure the signal of blood pressure due to too much miscellaneous signal. Please remove the signal interference source and measure again.
„Err3“	The measurement is abnormal, please measure again.
„ErrP“	Inflatable fails, please check the wristband, tie appropriately and then measure again.
„HI“	Inflation pressure is greater than 290+/-8mmHg (38,6kPa+/-1kPa), please measure again.
	When the battery is low, please replace four new alkaline batteries at the same time.

WHO BLOOD PRESSURE LEVEL OF THE REFERENCE TABLE

The device comes with hypertension warning strip alerts functions. As shown in the figure, we can learn blood pressure is distributed in the low level, ideal level, normal level (green), high-normal level (yellow), stage1hypertension, stage2hypertension, stage3hypertension(red) from the world health organization (WHO) blood pressure index level. After the measurement is completed „■“ present in the corresponding location according to the blood pressure level.

Systolic blood pressure (high blood pressure) is higher than 140mmHg or diastolic blood pressure (low blood pressure) is higher than 90mmHg, which is called hypertension by the world health organization. Generally believed systolic blood pressure is lower than 90 mmHg or diastolic blood pressure is lower than 60 mmHg, which is called hypotension.

MAINTENANCE AND PRESERVATION

Maintenance methods

1. Please clean the blood pressure monitor regularly.
2. If the body is dirty, please use a dry, soft cloth to wipe.
3. If the body is very dirty, you can use the cloth moistened with water or neutral detergent and then wipe the body.

Maintenance suggestions

1. Please do not put the blood pressure monitor body in the water.
2. Please do not make wristbands wet.
3. Please do not wipe the device with benzene, thinner, gasoline, etc.
4. Do not tightly fold cuff and air tube.
5. Please do not disassemble the device.
6. Do not repair the device yourself. If the products have quality problems or you have any questions about the correct measurement of the blood pressure monitor, please call customer service hotline.

Preservation methods

Please do not put the machine in the following places:

The places are easy to be splashed into the water.

The places are in direct sunlight, high temperature, high humidity, dust, or may contact with the corrosive gas or water.

The places are sloping and can produce vibration, shock. Do not drop or hit the body.

If you do not comply with the above considerations and other proper use methods, the company will not assume responsibility for the quality.

Disposal

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

EMC INFORMATION

Precautions:

This product conforms to the electromagnetic compatibility requirements of YY0505:2012/IEC 60601-1-2:2014 standard;

The user shall use the device according to EMC information provided in this manual. Portable and mobile RF communication equipment may affect the performance of the product. Avoid strong electromagnetic interference, such as close to mobile phones and microwave ovens. The guidance and manufacturer's statement are detailed in the table below

Warnings:

This product should not be used close to or stacked with other devices. If it must be used close to or stacked with other devices, it should be observed to verify that it can operate normally under the configuration used.

Electromagnetic Compatibility Information

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

EMC Information

Guidance and declaration of manufacturer-electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group1	The device uses RF energy only for its internal
RF emissions CISPR 11	Class B	Therefore, its emissions are very low and are not likely to cause any interference in nearby electronic equipment.

Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	N/A	

Warning:

1. The Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.)
2. Portable RF communications equipment(including peripherals such as antenna cables and external antennas)should be used no closer than 30 cm(12 inches) to any part of the [ME EQUIPMENT or ME SYSTEM], including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Guidance and declaration of manufacturer-electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge(ESD) IEC61000-4-2	±8kV contact ±2kV,±4kV t8kV,±15kVair	±8kV Contact 2kV, ±4kV ±8kV, ±15kVair	Floors should be wood,concrete or ceramic tile. If floors are covered with synthetic material,the relative humidity should be at least 30%
Electrostatic transient/burst IEC61000-4-4	±2kVfor power supply lines +1kVforinput/output lines	N/A	
Surge IEC 61000-4-5	+1kV differential mode ±2kV common mode	N/A	

Voltage dips, short interruptions and voltage variations on power supply line-put lines IEC 61000-4-11	<5% UT(>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycle 70% UT(30% dip in UT) for 25 cycle <5% UT(>95% dip in UT) for 5 sec	N/A	
Power Frequency (50/60 Hz) magnetic field IEC 61000-4-8	30A/m; 50 Hz or 60 Hz	30A/m; 50 Hz or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3 RF Wireless Communication Equipment IEC 61000-4-3	3V 0.15 MHz to 80 MHz 6V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80% A _{max} at 1 kHz 10V/m 80 MHz to 2.7 GHz 380 MHz, 27V/m 450 MHz, 28V/m 710 MHz, 745 MHz, 780 MHz 9V/m 810 MHz, 870 MHz, 930 MHz 28V/m 1720 MHz, 1845 MHz, 1970 MHz 28V/m 2450 MHz, 28V/m 5240 MHz, 5500 MHz, 5785 MHz 9V/m	N/A 10 V/m 380 MHz, 27V/m 450 MHz, 28V/m 710 MHz, 745 MHz, 780 MHz 9V/m 810 MHz, 870 MHz, 930 MHz 28V/m 1720 MHz, 1845 MHz, 1970 MHz 28V/m 2450 MHz, 28V/m 5240 MHz, 5500 MHz, 5785 MHz 9V/m	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance 80 MHz to 800 MHz 800 MHz to 2.7 GHz

			where P is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer and d is the recommended separation distance in meters(m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:
--	--	--	---

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