Legally responsible person

Producer
Artsana S.p.A.
Saldarini Catelli 1
Grandate
22070
info@chicco-service.pl
48221100847

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reference. **WARNING:** Before using the digital thermometer, carefully read the user manual. This device is not a toy and must be kept in a safe place out of the reach of children. Never leave a child unattended during temperature measurement. The thermometer contains small parts (batteries etc.) which may be harmful if swallowed. Choking hazard: swallowing the thermometer tip and batteries can be fatal. Do not allow children to use this device without parental supervision. Do not use the thermometer by inserting it into the ear. The device is intended only for readings from the oral cavity, rectum, and under the armpit. Do not place the thermometer battery near sources of strong heat, as it could explode. If the device is not intended to be used for a long period, remove the battery. Using temperature readings for self-diagnosis is dangerous. Consult a physician to interpret the results. Self-diagnosis may cause worsening of existing conditions. Do not take measurements when the thermometer is wet, as the reading may be inaccurate. Do not bite the thermometer. Such action may cause breakage and/or injury. Do not attempt to disassemble or repair the thermometer. Such action may result in inaccurate readings. After each use the thermometer should be disinfected, especially if the device is used by multiple people. Do not insert the thermometer into the rectum by force. Stop insertion and discontinue measurement if pain occurs. Failure to follow this recommendation may cause injury. Do not use the thermometer orally after rectal use. For children aged two years or younger do not use the device orally or vaginally. If the device has been stored at a temperature above 5°C ~ 40°C (41°F ~ 104°F), before use leave it at room temperature for about 15 minutes. Do not clean or disinfect the thermometer during use. Any serious incident related to the device should be reported to the manufacturer and the competent authorities of the Member State in which the user and/or patient is located. **INTENDED**

USE Digital thermometers are intended for measuring human body temperature in normal modes orally, rectally or axillary. The devices are suitable for repeated clinical or home use by persons of all ages, including children under 8 years of age, under adult supervision. **INDICATIONS FOR**

USE The device is intended for monitoring patients' body temperature **READ CAREFULLY**

BEFORE USE This digital thermometer provides a quick and highly accurate reading of a person's body temperature. To better understand its functions and ensure reliable results for many years, read all operating instructions first. The device complies with the following standards: ISO 80601-2-56 Medical electrical equipment - Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement, EN 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Supplementary standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and is compliant with the requirements of EN 60601-1-2 (EMC), IEC/EN 60601-1 (Safety). The manufacturer holds ISO 13485 certification. It complies with EU Regulation 2017/745 on medical devices. **CONTAINS** 1 thermometer, 1 operating manual, 1

case. **PRECAUTIONS** The effectiveness of the device may be impaired if at least one of the

following occurs:❑Using the device outdoors;❑Storing outside the temperature and humidity range specified by the manufacturer;❑Mechanical shock (e.g. test

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