Instructions for Use

Integrated Valve Tissue Expanders (IVTE)

904764 v1 – 02/2026

EUROSILICONE™

 GC Aesthetics®

CONTENTS

EN	Integrated Valve Tissue Expanders (IVTE).....	3
FR	Expanseurs tissulaires à valve intégrée.....	9
DE	Gewebeexpander mit integriertem Ventil	16
IT	Espansori tissutali con valvola integrata	23
PT	Expansores de tecido com válvula integrada	30
ES	Expansores tisulares con válvula integrada	37
DA	Vævsudvidere med integreret ventil	44
NL	Weefselexpanders met geïntegreerd ventiel	51
EL	διατεταμένες ιστών ενσωματωμένης βαλβίδας.....	58
SV	Vävnadsexpandrar med integrerad ventil.....	65
HU	Integrált szelepes szövetszítók	72
LV	Audu espanderi ar integrēto vārstu	79
LT	Audinių plėtikliai su integruotu vožtuvu	85
SK	Tkaníkové expandéry s integrovaným ventilom	92
BG	Тъканни експандери с интегрирана клапа.....	99
RO	Expandere tisulare cu valvă Integrată.....	106
TR	Entegre valfli doku genişleticiler	113
RU	Тканевые экспандеры со встроенным клапаном	119
CZ	Tkářňové expandéry s integrovaným ventilem.....	127
ET	Integreeritud klapiga koeekspanderid	133
HR	Ekspanderi tkiva s integriranim ventilom.....	140
PL	Ekspandery tkankowe ze zintegrowanym dozownikiem	147
SL	Tkivni ekspanderji z integriranim ventilom	154
FI	Kudoslaajentimet integroidulla venttiillä.....	160
NO	Vevsekspandere med integrert ventil	167
IS	Vefjapenjarar með innbyggðum ventli	174

Integrated Valve Tissue Expanders (IVTE)

DESCRIPTION

Integrated Valve Tissue Expander (IVTE) is a sterile, single-use medical device that is surgically implanted designed to gradually stretch the tissue adjacent to the implantation area (normally the skin). Once the tissue has expanded sufficiently, the device is removed and the extra tissue can then be used where it is required, such as preparing for a permanent prosthesis implantation. Integrated Valve Tissue Expanders must not be used as permanent prostheses as they are designed as temporary devices.

Integrated Valve Tissue Expanders consist of a silicone elastomer shell with a magnetically detectable integrated injection site. Integrated Valve Tissue Expanders are of anatomical shape, are available in a range of sizes and are designed for temporary subcutaneous implantation. Please refer to the manufacturer's website (<https://www.gcaesthetics.com/>) and/or promotional literature for full details of product range. Each Integrated Valve Tissue Expander is supplied in a cardboard box with one External Locating Device (ELD), both packed in a double peel pouch. The packaging also contains 8 removable labels indicating the characteristics of the implant (Lot number, Reference number) and a Volume Fill Record and Implant card.

Expansion of the Integrated Valve Tissue Expander is achieved by inflating the shell of the expander with varying volumes of sterile isotonic saline solution at regular intervals, as per the 'Filling Procedure' below. Integrated Valve Tissue Expanders have an integrated injection site with self-sealing technology which helps to facilitate the transfer of saline from the needle to the product.


Figure 1 Schematic of the Integrated Valve Tissue Expander device

INTENDED PURPOSE/USE

The IVTE is a temporary, single use, sterile tissue expander designed for the development/enlargement of breast skin-flap dimensions. Expansion is achieved by inflation of the shell of the expander with varying volumes of sterile isotonic saline solution at regular intervals.

LINK TO E-DOCUMENT (Safety and Performance information relevant to the user)

www.gcaesthetics.com/surgeons/IFU

INDICATIONS

Integrated Valve Tissue Expanders are indicated for the following:

- Reconstruction of the breast following subcutaneous mastectomy and other mastectomy procedures or trauma.
- Breast underdevelopment and combined breast and chest wall abnormalities.
- Scar/defect revision.

CONTRA-INDICATIONS

The use of these devices is contra-indicated in patients who have one or more of the following conditions:

- Use of devices which emit a static magnetic field (ferromagnetic materials, MRI).
- Implanted devices that would be affected by static magnetic fields.
- Radiation Therapy is a relative contraindication for the use of tissue expanders in breast reconstruction.
- Use of Radiation Therapy in conjunction with Integrated Valve Tissue Expanders is at the discretion of the surgeon and radiation Oncologist. The presence of rare earth magnets may affect radiation delivery thereby leading to dosimetric inaccuracies of the target. Should adjuvant radiation therapy take place, implanted devices should not be serially inflated during the expansion phase.
- Insufficient tissue coverage.
- Tight scar tissue or post-radiotherapy fibrotic tissue may affect suitability (may increase risk of capsular contracture).
- Existing local or metastatic carcinoma in the area intended for expansion.
- Deficient vascularization of tissue in local area.
- Any existing pathology in the implantation area; treatments by irradiation, microwave or steroid diathermy.
- Existing painful conditions.
- Known allergy to Silicone or history of implant rejection.
- Patient signs of radiodermatitis: Sclero-atrophic skin, dermatitis, ulceration.
- Surgeon determined that the patient is physiologically / psychologically unsuitable.
- Active infection anywhere in the body.
- A history of compromised wound healing.
- Any other serious medical condition.
- Use of any drug, including those that would interfere with blood clotting or affect tissue viability, which may result in high surgical risk and/or significant postoperative complications.
- Pre-existing infection not resolved before device placement increases the risk of periprosthetic infection.
- Since there is no data concerning the concomitant use of these devices with dermal fillers, such use of dermal fillers in the proximity of the device (e.g., in the décolletage skin area) is contraindicated. Physicians should also advise patients against future use of dermal fillers in proximity to these devices.
- These devices are not intended for the treatment of children or pregnant or nursing women.

Magnetic Resonance Imaging (MRI) is strictly contra-indicated due to the integration of a rare earth magnet within the Integrated Valve Tissue expander.
N.B. It is the surgeons' responsibility to make a medical assessment as to the patient's suitability for implantation.

INTENDED USER

The Integrated Valve Tissue Expander is intended to be used by suitably qualified surgeons.

IMPORTANT: The surgeon is responsible for conducting the medical assessments regarding the patient's eligibility for implantation and for deciding on the most suitable surgical technique for both the patient and the type of implant chosen.

TARGET POPULATION

The Integrated Valve Tissue Expander is intended to be used in adult populations (excluding pregnant and nursing women), with no known contraindications.

INTENDED CLINICAL BENEFITS

The clinical benefits of the Integrated Valve Tissue Expander is as follows:

- Used for reconstructive purposes to increase the tissue area of the breast in preparation for permanent breast implants.
- Current textured, anatomical tissue expanders have the following benefits in patients who have undergone a mastectomy:
 - Highly predictable outcome
 - Generate a customised tissue pocket perfectly shaped for permanent implant placement
 - Increase tissue area of the breast
 - Minimise capsular contracture and implant migration
 - Maintain well-defined inframammary line
 - Fewer nipple deviation
 - Satisfaction with breast aesthetics following mastectomy
 - Improved body image, self-esteem, sexual feeling, and attractiveness

PERFORMANCE CHARACTERISTICS

The Integrated Valve Tissue Expanders are designed and manufactured according to the current EN ISO 14607 standard concerning mammary implants. The performance characteristics of the Integrated Valve Tissue Expanders in terms of:

- Shell integrity
- Implant resistance
- Injection Site competence
- Volume
- Surface texture

Had been validated and comply with the requirements of EN ISO 14607.

LINK TO THE SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The summary of Safety and effectiveness will be available on the Eudamed Website when it is live <https://ec.europa.eu/tools/eudamed/#/screen/home> and the company Website (<https://www.gcaesthetics.com/>)

PATIENT INFORMATION

Any surgical procedure can have complications and risks. The patient should be well counselled on the risk/benefit relationship, by the surgeon. Each of the possible complications and warnings should be discussed with the patient prior to the decision to proceed with surgery.

The possibility of explant surgery taking place at any time during implantation period should also be discussed with the patient. It is the surgeons' responsibility to advise the patient of all potential complications and risks before surgery.

Each patient should receive the patient information booklet during their initial consultation, to allow time for the patient to read and understand the information on risks, follow-up recommendations and benefits associated with breast implants in order to make an informed decision about proceeding with the surgery. The patient information booklet also contains a 'Patient Informed Consent Form' which allows the patient to retain copies of their implant details such as lot number for their records. This document is available in PDF format from our website www.gcaesthetics.com

LIFETIME

Tissue expanders are temporary devices that are not intended for permanent implantation. They should be removed once adequate tissue has been developed. Tissue expanders are typically implanted for up to six months.

Our products must be stored under normal conditions, must be protected from humidity, kept out of direct sunlight, and must be stored with the arrows symbols pointing upwards. Where these conditions are respected. The products have a five-year shelf-life.

IMPLANT CARD

After surgery, the surgeon should provide to the patient:

- The implant card which should be filled out according to the leaflet Instructions for Healthcare Professionals filling out the implant Card.
- The patient leaflet which is available on the Information website for patients referenced on the card.

The patient should be informed that the implant card must be kept on their person at all the times and enables:

- The patient to identify the implanted devices and to get access to other information related to the implanted device (e.g. via EUDAMED, and other websites).
- The patients to identify themselves as persons requiring special care in relevant situations e.g. security checks.

Emergency clinical staff or first responder to be informed about special care/needs for relevant patients in case of emergency situations.

WARNINGS

It is the responsibility of the surgeon to provide the patient with advice of the possible warnings associated with the use of Tissue Expanders prior to surgery.

INTERFERENCE WITH INVESTIGATIVE AND TREATMENT TECHNIQUES

- As Integrated Valve Tissue Expanders contain magnetic elements, Magnetic Resonance Imaging (MRI) examinations are strictly contra-indicated.
- The decision to use Radiation Therapy in conjunction with Integrated Valve Tissue Expanders should be made by the surgeon and the radiation oncologist (Also see contraindications).

MAGNETIC FIELDS

DO NOT use Integrated Valve Tissue Expanders in patients with implanted devices that would be affected by a static magnetic field, such as Pacemakers, Implantable Cardioverter defibrillators (ICD), Artificial Sensing Devices and Drug Infusion Devices.

TEMPORARY DEVICE

A tissue expander is a temporary device and is not intended for permanent implantation. The tissue expander should be removed once adequate tissue has been developed. Tissue expanders are typically implanted for up to 6 months.

DAMAGE OF IMPLANTS - BREAKAGE AT SURGERY

Extreme care should be taken in the use and handling of implants to minimize damage and to ensure aseptic presentation. The structural integrity of the device must be verified before implantation. Damaged products should not be implanted, do not attempt to repair damaged products.

All devices have been produced by established manufacturing techniques and under strict quality control standards, however improper handling may result in some breakage of implants during handling or in surgery, both in initial and any subsequent surgeries. The utmost care must be taken to avoid contact with sharp objects, such as surgical instruments, suture needles or hypodermic needles.

An extra product should be available at the time of surgery in the event of damage to, or contamination of, the device.

SINGLE USE

Tissue Expanders are intended for single use only. Do not reuse explanted products. Do not re-sterilize any product. Explanted products should not be reused because re-cleaning and resterilization procedures may not adequately remove biological residues, such as blood, tissue and other matter, which could retain resistant pathogens.

ALTERATION OF DEVICE

Do not attempt to repair damaged products. Alterations to the original design/fabrication should not be attempted. Unauthorized alteration of devices voids all warranties, express or implied.

USE OF DRUGS

If antibiotic or steroid drug therapy is indicated, it is the individual surgeon's responsibility to ensure that the drugs are compatible and efficacious when in contact with medical grade silicone elastomers and stainless steel. Therapy which relies on the diffusion of drugs through the device wall is not endorsed.

INFLATION OF THE ENVELOPE

Meticulous preparation of the skin over the injection site must be completed each time the injection site is penetrated. The physician must verify the location and proper orientation of the injection site prior to each insertion of the filling needle, which can be achieved by using the External Location Device (ELD) supplied with the Integrated Valve Tissue Expander device. Failure to do so may lead to puncture of the expansion envelope, damage to the injection site and/or subject the patient to undue trauma. The site must be penetrated with a 21-gauge standard bevel hypodermic needle. Other types of needle or larger needles may damage the injection site and prevent resealing. A new needle should be used each time the injection site is penetrated, as contact of the needle tip with the stainless steel needle guard may deform the needle tip and render it unsuitable for re-use. Inflation of the expansion envelope beyond the indicated fill volume (detailed on the Tissue Expander label) is not endorsed. Following each inflation injection the patient should be observed for possible complications. Failure of the expansion envelope to inflate may result from improper placement of the filling needle into the inner lumen of the injection site or from damage to the expansion envelope.

OVER-INFLATION OF THE DEVICE

Over inflation of the device may affect the integrity of the device and cause trauma to the surrounding tissue.

PLACEMENT OF DEVICE

Care must be taken during placement of the device to ensure that the envelope does not fold over the injection site. This may lead to accidental needle penetration of the envelope.

CLEANING PRIOR TO RE-USE OF ELD

The ELD is supplied sterile (upon first use) and shall be used initially to inflate the expander during the expander placement procedure in theatre. Immediately prior to any subsequent use in a NON-STERILE environment, the ELD must first be cleaned following the instructions within this document.

MAXIMUM NUMBER OF USES OF ELD

The ELD can be cleaned for re-use in a non-sterile environment up to a maximum of 15 times. The manufacturer recommends that users utilise the Volume Fill Record supplied with each device in order to track the number of cleaning cycles used.

ELD CONTACT WITH BROKEN SKIN

The ELD should not come into contact with broken skin during subsequent inflations of the expander following initial surgery, as although clean the ELD is no longer sterile.

ELD USED FOR ONE PATIENT ONLY

Each ELD should be used for one patient only. The manufacturer recommends that users utilise the Volume Fill Record supplied with each device in order to track the patient details associated with each ELD. The detachable labels provided with each device may be attached to the Volume Fill Record in order to better track the ELD associated with each patient.

NON-USE OF ALTERNATIVE ELD

The device is only to be used with manufacturer's ELDs, no alternative ELD from a different manufacturer can be used during the implantation process of the device.

COMPLICATIONS

Potential complications associated with all surgical procedures should be discussed with the patient. These include but are not limited to: infection (see below); haematoma (see below); serous fluid accumulation (see below); loss of sensation (see below); reaction to medication; nerve damage; patient intolerance to any foreign implant; and poor wound healing.

BREAST IMPLANT ASSOCIATED ANAPLASTIC LARGE CELL LYMPHOMA (BIA-ALCL)

European safety information, the US FDA and current scientific literature have identified an association between breast implants and the development of a cancer of the immune system called breast implant-associated anaplastic large cell lymphoma (BIA-ALCL). This specific entity is included in the WHO 2016 classification under the terminology "ALCL – AIM". This cancer occurs more commonly in patients with textured breast implants rather than smooth implants, and typically develops many years after breast implant surgery. At this time rates of BIA-ALCL are considered to be low.

In most patients, BIA-ALCL is found within the fluid or scar tissue next to the implant. The main symptoms of BIA-ALCL are persistent swelling aforementioned symptoms, or if BIA-ALCL is suspected, it is advised that patients should be assessed to rule out per-implant ALCL. When testing for BIA-ALCL, collect fresh seroma fluid and representative portions of the capsule and send for pathology tests to rule out ALCL. Diagnostic evaluation should include cytological evaluation of seroma fluid with Wright Giemsa stained smears and cell block immunohistochemistry testing for cluster of differentiation (CD) and Anaplastic Lymphoma Kinase (ALK) markers. Most cases of breast implant-associated ALCL are treated by removal of the implant and the capsule surrounding the implant and some cases have been treated by chemotherapy and radiation. Although treatment is usually successful BIA-ALCL can be fatal if left untreated, hence, early detection and treatment is key.

In the case that your patient is diagnosed with BIA-ALCL, a report should be made to the manufacturer via the e-mail address: productsupport@gcaesthetics.com and an individual treatment program should be initiated by a multidisciplinary team.

BREAST IMPLANT ILLNESS

A small proportion of women who have breast implants or tissue expanders, for both aesthetic and reconstructive purposes, self-identify as having a number of symptoms they believe arise from the presence of their breast implants. Whilst not a medical diagnosis, they refer to their symptoms as 'Breast Implant Illness' (BII). These multiple generic symptoms are very varied and all cases have been self-diagnosed and self-reported. They include, but are not limited to, flu-like symptoms such as extreme tiredness, brain fog, joint aches, immune-related symptoms, sleep disturbance, depression, hormonal issues, headaches, hair loss and chills.

There are a variety of other reasons these symptoms might be found including background illnesses or hormonal changes. In addition, there are a number of scientific studies investigating similar symptoms experienced by women in the general population both with and without breast implants. On average, around 50% of women who self-identify as having 'BII' feel that their symptoms improve after implant removal – sometimes temporarily and sometimes permanent. It therefore appears that removing breast implants does not necessarily improve symptoms in everyone. To date, there is no research demonstrating which symptoms may or may not improve with implant removal.

Currently, there are no tests that can confirm 'BII'. Research continues in this area, specifically around patients who have autoimmune disease or a predisposition to autoimmune disease. Breast implantation is contraindicated in women with a suppressed or compromised immune system - see CONTRA-INDICATIONS section above. Medical advice should be sought if patients suspect they have BII. Their symptoms may not be related to the implants and other medical investigations should not be overlooked or ignored in determining the cause of their symptoms

CAPSULAR CONTRACTURE

Formation of a fibrous tissue capsule surrounding the device is a normal physiological response to a foreign body. Contracture of the fibrous capsular tissue surrounding the device may result in firmness, discomfort or pain, distortion, palpability or displacement of the device.

DEFLATION OF TISSUE EXPANDERS AND RUPTURE

Deflation of the expansion envelope is possible at any time and the containment of saline cannot be guaranteed. If leakage is suspected the device should be removed. The patient should be informed of the deflation potential of the device prior to the decision to proceed to surgery. Wrinkling and/or creasing of the expansion envelope is inherent to the intended use of this device and may result in weakening and deflation of the expansion envelope, especially when left under filled for a period of time. Careful placement would be advised to avoid fold flaws to ensure well-tolerated expansion. Damage of devices during handling or in surgery can cause early or late deflation. Rupture can occur post-operatively (within the implantation period) from damage to the implant during handling or surgery. Envelopes which have been weakened by nicking can rupture at the nick-point postoperatively. The patient should be advised that excessive manipulation or trauma could rupture the expansion envelope and cause deflation. Leakage can also occur through the injection site. It is possible to overlay the lower pole of the expander over the lower edge of the filling port thus care should be taken when placing un-inflated integrated valve expanders in the pocket to avoid mistaking the lower pole as the filling port.

EXTRUSION

Inadequate tissue coverage may result in exposure and extrusion of the device. Skin breakdown can occur with the use of steroid drugs. Tissue sloughing may result from poor vascularization. Signs of skin inflammation, such as tenderness and erythema, may indicate thinning or erosion and must be promptly investigated.

INFECTION

Do not expose the device or filling accessories to contaminants, which increases the risk of infection.

Infection is an inherent risk following any type of invasive surgery. Signs of acute infection include erythema, tenderness, fluid accumulation, pain and fever. Signs of subclinical infection may be difficult to detect. Postoperative infections should be treated aggressively according to standard medical practices to avoid more serious complications. Infection that is unresponsive to treatment or necrotizing infection may require device removal.

Capsular contracture may be related to infection in the area surrounding the device.

HAEMATOMA/SEROMA

Haematoma and serous fluid accumulation are complications associated with any type of invasive surgery. Postoperative haematoma and seroma may contribute to infection and/or capsular contracture. Postoperative haematoma and seroma may be minimized by meticulous attention to haemostasis during surgery, and also possibly by postoperative use of a closed drainage system. Persistent, excessive bleeding must be controlled before the device is implanted. Any postoperative evacuation of haematoma or seroma must be conducted with care to avoid contamination or damage to the device.

BURNING SENSATION

There may be a local skin reaction following implantation and use of tissue expanders. This can cause local temporary redness (erythema) or rashes during the expansion process

ITCHING SENSATION

The risk of temporary or permanent dysesthesia exists following any invasive surgical procedure. Careful surgical technique can minimize, but not preclude, this risk. It is normal for the surrounding skin to experience itching sensation, chest pain, tenderness or tightness during the expanding process. However, persistent itching or burning sensations, swollen lymphs/glands, and/or local reactions (including inflammation and/or fluid discharge) may be a sign of an infection or foreign body reaction. The risk of neurological impairment increases with more extensive surgery.

PAIN

Any unexplainable pain, not associated with any of the complications listed herein, must be promptly investigated.

PHYSIOLOGICAL CONDITIONS

Underlying physiological conditions may contribute to reconstructive failure and may include diabetes, obesity and hypertension.

IMMUNE RESPONSE

The risk of developing an immune response following implantation of tissue expanders, although not biologically plausible cannot be completely ruled out.

NECROSIS

Impaired vascular perfusion may lead to mastectomy flap, skin and/or nipple necrosis. This can cause poor wound healing and device exposure requiring revision surgery or removal of the device.

HYPERTROPHIC SCARRING

Raised scarring may occur as a result of tissue expansion. This thickened scarring may be temporary or permanent.

CHEST WALL DEFORMITY

Chest wall deformity is a potential complication following implantation and expansion of tissue expanders in the chest region. This may lead to temporary or permanent alteration of the chest wall.

DISPLACEMENT/MALPOSITION

Displacement and malposition of the tissue expander is an associated complication. Surgical intervention and device explanation may be required.

SKIN DISCOLOURATION

The skin in or near the surgical site can appear either lighter or darker than surrounding skin. Although uncommon, swelling and skin discoloration may persist for long periods of time and, in rare situations, may be permanent.

INSTRUCTIONS FOR USE

STORAGE AND TRANSIT CONDITIONS

Our products must be stored under normal conditions, must be protected from humidity, kept out of direct sunlight, and must be stored with the arrow symbols pointing upwards. Where these conditions are respected, the products have a five year shelf-life.

SUPPLIED STERILE

Devices are supplied in sterile form, processed by validated strictly controlled sterilization cycles using ethylene oxide. Sterility of the device is maintained only if the package is intact and undamaged. If the package has been opened or damaged the product must not be used. Authorization must be requested to return the product.

SINGLE USE

The tissue expander is supplied sterile and is intended for SINGLE USE ONLY. DO NOT REUSE EXPLANTED PRODUCTS. DO NOT RESTERILIZE ANY PRODUCT. The Tissue Expander should not be re-used because re-cleaning and re-sterilization procedures may damage the device integrity and may not adequately remove biological residues, such as blood, tissue and other matter, which could retain resistant pathogens.

PACKAGING

Sterile product is supplied in a sealed, double primary package, inside a protective carton box. Sterility is not guaranteed if the package has been damaged or opened. DO NOT USE THE DEVICE IF THE PACKAGING IS DAMAGED OR OPENED. The labelling supplied with the device provides the necessary information for traceability and identification of the product. The packaging contains 8 removable labels indicating the characteristics of the implant (lot number and reference number). These labels should be attached to the patients' medical file. The manufacturer recommends that the patient is provided with one of these labels to maintain the traceability of the implanted product.

TO OPEN PACKAGED STERILE PRODUCT

- Peel open outer package under clean, aseptic conditions, over sterile field, allowing sealed inner package to fall gently into the field.
- Attach patient record portion of the inner label to the patients' chart.
- Using aseptic precautions, peel open the inner package.

PRODUCT EXAMINATION AND HANDLING

- Product should be visually examined for any evidence of particulate contamination or damage. The product should be tested for leakage by partially filling with sterile isotonic saline and gently compressing or by filling with air and compressing the device while submerged in a sterile isotonic saline solution. Do not use if there is any evidence of leakage, damage or contamination. The device needs to be deflated following this test and prior to implantation.
- The device should be kept submerged in sterile water or normal saline prior to implantation to prevent contact with airborne and surgical field particulate contaminants.
- Do not immerse the implant in Betadine or Iodine containing solutions. If Betadine or Iodine containing solutions are used in the pocket, ensure that it is rinsed thoroughly so no residual solution remains in the pocket.

SURFACE CONTAMINANTS (TALCS, DUST, LINT, OILS) ON SURFACE OF IMPLANTS CAN CAUSE FOREIGN BODY REACTION, HANDLE WITH CARE WITH SURGICAL GLOVES (RINSED FREE OF TALC) WITH STRICT ASEPTIC TECHNIQUE. DO NOT IMPLANT CONTAMINATED PRODUCT.

Back-up implants must be readily available at the time of surgery for use in the event of contamination.

RESTERILISATION

Tissue Expanders are intended for SINGLE USE ONLY. DO NOT REUSE EXPLANTED PRODUCTS. DO NOT RESTERILIZE TISSUE EXPANDER OR ELD. Explanted products should not be reused because recleaning and resterilization procedures may damage the device integrity and may not adequately remove biological residues, such as blood, tissue and other matter, which could retain resistant pathogens.

SURGICAL PROCEDURE

Correct surgical procedures and techniques are the responsibility of the medical profession. Each surgeon must evaluate the suitability of the procedure based upon current accepted techniques, individual judgement and experience. Proper size and shape of the device must be determined for the individual patient by the surgeon. An incision should be of appropriate length to accommodate the style, size, and profile of the device. This will reduce the potential for creating excessive stress to the device during insertion.

POSITIONING OF INTEGRATED VALVE TISSUE EXPANDER

For optimal expansion, the pocket created should be larger than the base of the device. Typically inflation starts after a period of wound healing. Care must be taken during placement of the device to ensure that the envelope does not fold over the injection site. This may lead to accidental needle penetration of the envelope.

LOCATING THE INTEGRATED INJECTION SITE USING EXTERNAL LOCATING DEVICE

1. Remove the external locating device (ELD) from the sterile packaging.
2. Ensure that the ELD magnetic hourglass can freely pivot on the ELD crosshairs. If the hourglass becomes detached from the holder, it can be re-attached by clipping the hourglass neck through the gap in the holder, as per figure 2. Please observe the hourglass orientation configuration. The magnet should be placed downwards as indicated.
3. Holding the ELD as indicated in figure 3; place the base of the ELD holder on the patient tissue over the implanted tissue expander.
4. Move the ELD in a circular motion over the tissue until the ELD hourglass detects (points towards) the location of the injection port.
5. Carefully adjust the ELD position until the hourglass is pointing straight down. The tissue expander injection site is located when the hourglass is perpendicular to the tissue surface.
6. As indicated in figure 4, with the ELD hourglass centred on the tissue expander injection site, make a mark with an ink pen in each of the four notches around the anterior perimeter of the ELD holder.
7. Remove the ELD and carefully connect the opposing dots with a line, creating an injection 'crosshair'. The point where the two perpendicular lines intersect is the location of the injection site.


Figure 2


Figure 3


Figure 4

FILLING PROCEDURE

Integrated Valve Tissue Expanders are filled using a self-sealing valve system. When inflating, the surgeon should become familiar with the "click" that can be felt due to the needle making contact with the stainless-steel needle guard. This "click" should help assure the user that the needle is in the inner lumen of the injection site.

The maximum fill volume is specified on the product package label.
DO NOT use excessive force during any of the steps in the following procedure.
• At the time of insertion expel all air from the device.

- Identify the centre of the Injection Port using the External Locating Device provided, marking the area with an ink pen.
- Fill the device only with sterile saline for injection and only through the valve as described before. NEVER make injections through the device shell.
- Use a small gauge needle (21-gauge standard bevel hypodermic needle).
- Ensure the needle is introduced into the reservoir vertically.
- Ensure the point of the needle is in contact with the needle guard during the filling procedure.
- DO NOT WITHDRAW THE NEEDLE during the filling procedure.

EXTERNAL LOCATING DEVICE CLEANING INSTRUCTIONS

Initial expansion is performed in the operating theatre and must only be performed with a STERILE ELD as supplied by the manufacturer. The ELD must first be cleaned following the instructions below immediately prior to any subsequent re-uses of the ELD when used in a NON-STERILE environment. The ELD should be washed thoroughly with isopropyl alcohol (IPA). If necessary, a mild, non-oily soap may be used prior to washing in IPA to remove organic or particulate matter. Rinse the ELD copiously with clean, non-pyrogenic water following the IPA wash. The ELD should be dried immediately following the non-pyrogenic water wash. Do not use synthetic-based detergents or oil-based soaps.

The ELD can be cleaned for re-use in a non-sterile environment up to a maximum of 15 times. The ELD should not come into contact with broken skin during subsequent inflations of the expander following initial surgery, as although clean the ELD is no longer sterile. Each ELD should be used for one patient only.

REMOVAL OF DEVICE

When the desired amount of tissue has been created, the device is removed, and the flap is advanced. A noticeable amount of contraction will be observed so that the surface area first generated will be decreased.

MATERIALS DEVICE DESCRIPTION

These medical devices are manufactured from medical grade silicone materials which are suitable for long term implantation. Master Files have been filed with the U.S FDA. The following three tables provide quantitative and qualitative information on the materials and substances which patients can be exposed based on chemical characterization of representative devices. The materials and device have been subjected to biocompatibility testing and evaluation, and risk assessments to demonstrate their biological safety. However, individual responses to chemicals may vary, and all reactions cannot be predicted.

Table 1a. IVTE Device Materials

Device materials	Implant component
Silicone dispersions	Shell
Silicone dispersions	Closure Patch
Silicone dispersions	Base re-inforcement
Silicone dispersions	Integrated valve silicone septum
Silicone adhesive	Shell
Stainless Steel / Neodymium **	Integrated valve needlestop

Table 1b. IVTE Device Materials Quantities*

Component	TE-XF3 (%)	TE-XM3 (%)	TE -XL3 (%)
Dimethyl silicone elastomer	26.2-35.9	25.6-39.8	30.6-40.8
Fluorosilicone elastomer	12.8-17.8	12.6-19.8	15.2-20.3
Masterbatch	7.8-10.2	8.6-10.7	9.27-10.72
HCRA-35M	19.4-28.6	16.4-28.7	15.2-24.2
RTV adhesive	0.3-0.5	16.4-28.7	0.3-0.4

* A range of proportions of silicones and adhesive has been summarized in Table 1b. for each product family of IVTE.

** Please note that stainless steel needlestop and Neodymium magnet encased within pose minimal risk of being in contact with the body as this would need to be released via a rupture within the IVTE. Furthermore, Stainless steel is a biologically inert material and therefore no further reporting of data is deemed appropriate.

Table 2. Chemicals Released from the IVTE

Compounds	Whole device (µg/g or ppm)	Compounds	Whole device (µg/g or ppm)
Volatiles¹ - Chemicals that can be released by breast implants as a gas			
Benzene	Not Detected	Ethylbenzene	Not Detected
Bromobenzene	Not Detected	Hexachlorobutadiene	Not Detected
Bromodichloromethane	Not Detected	Isopropylbenzene	Not Detected
Bromoform	Not Detected	p-Isopropyltoluene	Not Detected
Butylbenzene (N- / tert- / sec-)	Not Detected	N-propylbenzene	Not Detected
Carbon tetrachloride	Not Detected	Naphthalene	Not Detected
Chlorobenzene	Not Detected	Styrene	Not Detected
Chloroform	Not Detected	1,1,1,2-Tetrachloroethane	Not Detected
2,4-Chlorotoluene	Not Detected	Tetrachloroethane	Not Detected
Dibromochloromethane	Not Detected	Trichloroethane	Not Detected
1,2-Dibromo-3-chloropropane	Not Detected	Trichloroethane (1,1,1- / 1,1,2-)	Not Detected
Dichlorobenzene (1,2- / 1,3- / 1,4-)	Not Detected	1,2,3-Trichloropropane	Not Detected
Dichloroethane (1,1- / 1,2-)	Not Detected	Trimethylbenzenes (1,2,4- / 1,3,5-)	Not Detected
Dichloroethene (1,1- / cis-1,2- / trans-1,2-)	Not Detected	Toluene	Not Detected
Dichloropropane (1,2- / 1,3-)	Not Detected	Trichlorobenzenes (1,2,4- / 1,2,3-)	Not Detected
1,1-Dichloropropene	Not Detected	o-Xylene	Not Detected
Ethanol, 2-trimethylsilyl	0.17 max.	m- & p-Xylene	Not Detected
Total volatiles	0.17		
Extractables² - Chemicals that can be released by breast implants following soaking in water			
Cyclic siloxanes (D4, D5, D6...)	Not Detected	Linear siloxanes (L3, L4, L5...)	Not Detected
Total extractables	43.97 max		

1. Not Detected means that the level of the individual volatile was below the quantitation limit of the test method. The quantitation limit was 0.021 µg/g for m- & p-Xylene and 2,4-Chlorotoluene, and 0.42 µg/g for all other compounds.

2. Not Detected means that the level of the individual extractable was below 24.2 µg/g, the quantitation limit of the test method.

Table 3. Heavy Metals Found in the IVTE

Heavy metals	Concentration (ppm)	Heavy metals	Concentration (ppm)
Arsenic	Not Detected	Chromium	Not Detected
Cadmium	Not Detected	Copper	Not Detected
Lead	Not Detected	Tin	Not Detected

Heavy metals	Concentration (ppm)	Heavy metals	Concentration (ppm)
Mercury	10 max.	Nickel	Not Detected
Vanadium	Not Detected	Platinum	Not Detected
Molybdenum	Not Detected	Zinc	0.015 ppm/ 15 ppb max.
Selenium	Not Detected	Magnesium	Not Detected
Cobalt	Not Detected	Manganese	Not Detected
Antimony	Not Detected	Thallium	Not Detected
Barium	Not Detected		

Not Detected means that the level of the individual element was below 0.25 ppm, the quantitation limit of the test method.

DISPOSAL OF DEVICE

There are no special disposal instructions for this device. Normal hospital protocol should be followed to dispose of this device.

REPORTING OF SERIOUS INCIDENT:

Any serious incident (incident that directly or indirectly led, might have led or might lead to any of the following: (a) the death of a patient, user or other person, (b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health, (c) a serious public health threat) should be reported to the manufacturer via the e-mail address: productsupport@gcaesthetics.com and to the competent authority of the Member State in which the user and/or patient is established.

RETURN OF GOODS POLICY

Prior to the return of any product:

1. Authorization must be requested in writing.
2. A decontamination certificate must be supplied with any product returned which has been in contact with body fluid.
3. If it is necessary to puncture the shell of any sealed shell product to assist with safe decontamination sterilization procedures, the area of mechanical interference should be marked with indelible marker on the surface of the product and reference made on the decontamination certificate.

WARRANTY

The manufacturer warrants that reasonable care was used in the manufacture of these products and will replace any product that the manufacturer's investigation shows to be defective at time of shipment. Patient selection, surgical procedures, post-surgical treatment and stresses, and handling of the devices are totally the responsibility of the customer. The manufacturer has no control over the conditions of use and cannot warrant good effect or against ill effect following the use of the device and shall not be responsible for any incidental or consequential loss or damage or expenses directly or indirectly arising from the use of this device. All other warranties (whether implied by law or otherwise) are excluded to the fullest extent permitted by law.

CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

BIBLIOGRAPHIC REFERENCES

Literature references are available upon request from the manufacturer.

LABELLING SYMBOL:

Indicates product dimensions.

Indicates not for use with MRI (magnetic resonance imaging)

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