Instructions for Use

Round Collection™
Round Silicone Gel-filled Breast Implants

904821-v1 - 09/2026

EUROSILICONE™

CONTENTS

EN	Round Silicone Gel-filled Breast Implants.....	3
FR	Implants mammaires ronds remplis de gel de silicone	9
DE	Runde mit Silikongel gefüllte Brustimplantate	16
IT	Impianti mammari rotondi riempiti in gel di silicone	23
ES	Implantes mamarios redondos rellenos de gel de silicona	29
PT	Implantes Mamários Redondos com Enchimento de Gel de Silicone	36
DA	Runde silikongelfyldte brystimplantater	43
NL	Ronde met siliconengel gevulde borstimplantaten	49
EL	Στρογγυλά εμφυτεύματα στήθους γέλης σιλικόνης.....	56
LV	Apalje, ar silikongēlu pildītie krūšu implanti.....	63
SK	Okrúhle prsné implantáty plnené silikónovým gélom	69
RO	Implanturi mamare rotunde cu gel de silicon	76
TR	Yuvarlak Silikon Jel Dolgulu Meme İmplantları.....	82
CS	Kulaté prsní implantáty vyplněné silikonovým gelem	89
ET	Ümarad silikoongeeltäitega rinnaimplantaadid.....	95
HU	Kerek, szilikongéllel töltött mellimplantátumok.....	102
PL	Okragłe implanty piersi wypełnione żelem silikonowym.....	108
SL	Okrogli prsni vsadki, napolnjeni s silikonskim gelom	115
SV	Runda silikongelfyllda bröstimplantat.....	121
BG	Кръгли гръдни импланти с пълнеж от силиконов гел.....	128
HR	Okrugli implantati za grudi punjeni silikonskim gelom	135
FI	Pyöreät silikonigeelitäytteiset rintaimplantit.....	141
NO	Runde silikongelfylte brystimplantater.....	147
LT	Apvalūs silikoniniu geliu užpildyti krūtų implantai.....	154

Only to be implanted in an appropriate medical environment by appropriately trained medical doctors who are qualified or accredited in accordance with national law.

DESCRIPTION

Eurosilicone Round Gel-Filled Mammary Implants are sterile, single-use medical devices that are surgically implanted to replace tissue lost by trauma, surgical intervention or to augment tissue for therapeutic purposes. Eurosilicone Round Gel-Filled Mammary Implants are intended for cosmetic breast augmentation surgery or breast reconstruction surgery. The round silicone elastomer shell may be produced smooth or with the outer surface micro textured to offer a disruptive contact surface for collagen interface and filled with a cohesive or high cohesive gel. Each gel-filled mammary implant is supplied in a cardboard box and packaged within a sealed, double primary pack. The packaging also contains removable labels indicating the characteristics of the implant (Lot number, Reference number...), implant card, patient, and user leaflets.

INTENDED PURPOSE/USE

Eurosilicone Round Gel-Filled Mammary Implants are intended for cosmetic breast augmentation surgery or breast reconstruction surgery.

MEDICAL INDICATIONS

Round Gel-Filled Mammary Implants are indicated for

- Surgical Augmentation and contour correction of congenital anomalies of the breast,
- Surgical Reconstruction of the breast following subcutaneous mastectomy and other suitable mastectomy procedures or trauma,
- Surgical correction of Combined breast and chest wall abnormalities,
- Replacement of breast implants for medical reasons.

COSMETIC INDICATIONS

Round Gel-Filled Mammary Implants are indicated for

- Cosmetic Augmentation Surgery,
- Replacement of breast implant for cosmetic reasons.

CONTRA-INDICATIONS

The use of these implants is contraindicated in patients who have one or more of the following conditions:

- FIBROCYSTIC disease,
- INSUFFICIENT TISSUE covering due to radiation damage on the chest wall, tight thoracic skin grafts or radical resection of the pectoralis major muscle,
- EXISTING LOCAL OR METASTATIC CARCINOMA of the breast,
- Suppressed/compromised immune system,
- A history of sensitivity to foreign materials or atopy,
- PHYSIOLOGICALLY/PSYCHOLOGICALLY unsuitable patient,
- ACTIVE INFECTION anywhere in the body,
- Recent history of breast abscess,
- A history of compromised wound healing or regional swelling,
- KELOID SCARS,
- Any other serious medical condition,
- These devices are not intended for persons who are less than 18 years old or pregnant or nursing women.

INTENDED USER

The Eurosilicone Round Gel-Filled Mammary Implants are intended to be used by suitably qualified surgeons who have received appropriate training on how to use the device safely.

IMPORTANT: The surgeon is responsible for conducting the medical assessments regarding the patient's eligibility for implantation and for deciding on the most suitable surgical technique for both the patient and the type of implant chosen. The user should consider any specific risks that may be applicable to activities of the patient (for example profession, sports or other activities regularly performed by the patient) before the implantation. Eurosilicone recommends the surgeons to perform a post implantation monitoring in order to identify any potential undesirable side-effects.

TARGET POPULATION

Eurosilicone Round Gel-Filled Mammary Implants are intended to be used in adult populations (excluding pregnant and nursing women), with no known contraindications.

INTENDED CLINICAL BENEFITS FOR MEDICAL INDICATIONS

The clinical benefits of the Eurosilicone Round Gel-Filled Mammary Implants are as follows:

- Replace breast tissue that has been removed following mastectomy
- Replace breast tissue that has been removed due to trauma
- Enhance breast tissue affected by tuberous or congenital breast abnormality
- Improve patient's quality of life (based on self-esteem, body image, sexual satisfaction)
- Improve patient's psychology in terms of levels of depression and anxiety

EXPECTED EFFECTS FOR COSMETIC INDICATIONS

The expected effects of the Eurosilicone Round Gel-Filled Mammary Implants are as follows:

- Improve patient's quality of life (based on self-esteem, body image, sexual satisfaction)

PERFORMANCE CHARACTERISTICS

The Eurosilicone Round Gel-Filled Mammary Implants are designed and manufactured according to the current EN ISO 14607 standard concerning mammary implants.

The performance characteristics of the Eurosilicone Round Gel-Filled Mammary Implants in terms of:

- Shell integrity
- Implant resistance
- Compatibility between the filling gel and the shell
- Diffusion from the implant
- Volume
- Surface texture
- Biocompatibility

Had been validated and comply with the requirements of EN ISO 14607.

MATERIALS

Shell and patch: Silicone elastomer. Filling: Silicone gel

LINK TO THE SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The summary of Safety and effectiveness will be available on the EUDAMED Website when it is live <https://ec.europa.eu/tools/eudamed/#/screen/home> and the company Website (<https://www.gcaesthetics.com/>)

For access via Eudamed:

- Click on the link above for the Eudamed home page
- Click on the module called Certificates (Issued or Refused)

- Within the search criteria under Certificate number box insert MDR 765923 and press Search at the bottom of the page
- On the right-hand side of the page under View Detail click on the eye icon.
- On the left-hand side of the page click on SSCP
- Download the SSCP master document by clicking on the link

PATIENT INFORMATION

Any surgical procedure can have complications and risks. Mammary implant surgery is known to provide psychological satisfaction to patients, but like any surgical procedure it can have potential complications and risks. Mammary implantation is an elective procedure, and the patient should be well counselled on the risk/benefit relationship, by the surgeon. Each of the possible complications and warnings should be discussed with the patient prior to the decision to proceed with surgery. The recommendation for the patient has to consider any previous procedures, accidents, conditions, medications, or other simultaneous treatments of the consumer that may affect the procedure (for example skin, diseases, traumas, and auto-immune diseases).

Each patient should receive the patient information leaflet (provided together with the products) and the patient information booklet during their initial consultation, to allow time for the patient to read and understand the information on risks, follow-up recommendations and benefits associated with breast implants in order to make an informed decision about proceeding with the surgery. The patient information booklet also contains a 'Patient Informed Consent Form' which allows the patient to retain copies of their implant details such as lot number for their records. This document is available in PDF format from our website www.gcaesthetics.com/.

LIFETIME

The patient should be informed that this device is not a lifetime device and may require removal or replacement, which could entail another surgical procedure. Patients must be informed that the implant must not be considered definitive due to possible individual and potential physiological reactions, the inherent to the design of silicone implants or the various implant procedures. The implant may be removed or replaced at any time, which requires revision surgery. The manufacturer guarantees no implant duration period for its implants.

Clinical monitoring of Eurosilicone implants have enabled, to date, to establish implant durability of 82.4% after 10 years (risk of a first reoperation = 17.6%, using the Kaplan Meier method).

MRI COMPATIBILITY

Although Eurosilicone implants have not been specifically tested for use in an MRI, please note that Eurosilicone's Gel-filled Mammary Implants are all manufactured from medical implant grade silicone materials which are compatible with MRI Scan. The implantable silicone grade material is the same as that used in other Eurosilicone silicone implants where patients have undergone MRI scans and no compatibility issues have been reported to date.

STEROID USE

The patient should be informed to consult a physician before using steroid drugs in the implant area to avoid extrusion of the implant.

INTERFERENCE

During a mammography, the radiologist must be informed that there is an implant in the breast in order to adapt the mammography pressure accordingly. The surgeon or doctor must be notified about the implant in the event of any surgical intervention in the breast area.

The presence of breast implants could delay the detection of breast cancer by self-exam. For this reason, patients must be informed that they need to consult a surgeon and/or doctor for appropriate medical monitoring, as well as regular breast cancer screening. Patients must also be informed that they need to consult a doctor if they suspect complications particularly in the event of trauma or compression caused, for example, by violently massaging the breasts, a sporting activity or wearing seat belts.

IMPLANT CARD

After surgery, the surgeon should provide to the patient:

- the implant card which should be filled out according to the instructions for Healthcare Professionals filling out the EUROSILICONE International Implant Card.
- the patient leaflet which is available on the Information website for patients referenced on the card.

The patient should be informed that the implant card must be kept on their person at all the times and enables:

- the patient to identify the implanted devices and to get access to other information related to the implanted device (e.g. via EUDAMED, and other websites).
- the patients to identify themselves as persons requiring special care in relevant situations e.g. security checks.
- emergency clinical staff or first responder to be informed about special care/needs for relevant patients in case of emergency situations.

WARNINGS

ALTERATION OF IMPLANTS

No alteration should be made to the mammary implant before implantation. Alteration of implants voids all warranties, express or implied.

DAMAGE OF IMPLANTS - BREAKAGE AT SURGERY

Extreme care should be taken in the use and handling of implants to minimize the potential for breakage of shells. All implants have been produced by established manufacturing techniques and under strict quality control standards, but there may be some breakage of implants during handling or in surgery, both in initial and any subsequent surgeries.

Extreme care must be taken to avoid unintentional damage to the implant during the implantation or explantation surgery.

- Do not contact the implant with sharp surgical instruments or devices such as scalpels, forceps, hemostats, suture needles and hypodermic needles.
- Do not contact the implant with blunt surgical instruments such as clamps, retractors, and dissectors.
- Do not contact the implant with cautery devices.
- Do not use excessive manipulation, force, or stress.

Implants should be carefully inspected for structural integrity prior to use. Damaged products should not be implanted, do not attempt to repair damaged products. An extra product should be available at the time of surgery in the event of damage to, or contamination of, the implant. Recommended procedures for testing, examination and handling of products should be meticulously followed to assure proper use of implants. Patients should be instructed to inform other treating physicians of the presence of implants to minimize the risk of damage.

INTERFERENCE WITH MAMMOGRAPHY

Standard positioning techniques have shown significant limitations when used for imaging augmented breasts. The implant can interfere with the detection of early breast cancer through mammography by obscuring some underlying breast tissue and/or by compressing overlying tissue which can 'hide' suspicious lesions in the breast. Patients should be instructed to request radiologists who are experienced with the most current radiological techniques and equipment for imaging breasts with implants, and to inform their radiologists of the presence, type, and placement of implants.

SINGLE USE

Implants are intended for SINGLE USE ONLY. DO NOT REUSE EXPLANTED PRODUCTS. DO NOT RESTERILIZE ANY PRODUCT.

Explant products should not be reused because recleaning and re-sterilization procedures may not adequately remove biological residues, such as blood, tissue, and other matter, which could retain resistant pathogens and may also affect the performance of the implant.

USE OF DRUGS

Eurosilicone can neither predict, nor warrant, the safety of intraluminal introduction of any drugs, including, but not limited to, aesthetic, steroid, antibiotic, and vitamin solutions. If such use is contemplated, the appropriate drug manufacturer should be consulted.

COMPLICATIONS

Possible complications known to be associated with gel-filled mammary implants are as follows:

Other potential complications associated with all surgical procedures should be discussed with the patient. These include but are not limited to: infection (see below); hematoma (see below); serous fluid accumulation (see below); loss of sensation (see below); reaction to medication; nerve damage; patient intolerance to any foreign implant; and poor wound healing.

ACCIDENTAL RUPTURE OF IMPLANTS IN SURGICAL POCKET

If this occurs, the gel mass may be penetrated with index finger of double gloved hand, whilst exerting pressure on the breast with other hand. The gel can be manipulated out of the cavity, within the hand, and the outer glove can then be pulled over the mass for disposal. Blot the surgical pocket with gauze sponges. Isopropyl alcohol will aid the removal of gel from instrumentation.

ANAESTHESIA

Risks of Anaesthesia

There are three types of anaesthesia used during surgery, all of which carry some level of risk.

LOCAL ANAESTHETIC

The lowest level of risk is a local anaesthetic, which involves minimal I.V. sedation with injection of local anaesthetic in the area beneath the breasts. Some patients may have an allergic reaction to the local anaesthetic or experience a rapid heartbeat due to the epinephrine that is used to reduce bleeding. There may be some discomfort intermittently throughout the procedure with a local anaesthetic.

SEDATION:

A second method of anaesthesia is I.V. sedation (intravenous medications delivered without a tube in the throat). This carries risks of respiratory distress, reactions to the medications or medication overdose. It is recommended that a professional trained in the use of I.V. sedation be available to monitor and administer the medications. The cost for this type of anaesthesia is generally higher due to the cost of the medication and personnel to administer them.

GENERAL ANAESTHESIA

A third method is general anaesthesia in which a patient is asleep during the surgical procedure. The risks of general anaesthesia are the same as those of a general anaesthetic used for other operations and can involve respiratory problems and/or blood clots in the legs. The expense of general anaesthesia is higher because of the anaesthetist needed and the equipment and medications used.

ASIA SYNDROME

Autoimmune/Inflammatory syndrome induced by adjuvants. Patients can present myalgia, myositis, muscle weakness, arthralgia and/or arthritis, chronic fatigue, neurological manifestations, cognitive impairment, pyrexia, dry mouth. Diagnosis is confirmed by tests and markers of immune response. Removal of inciting agent generally induces improvement.

In the case that your patient is diagnosed with ASIA syndrome, a report should be made to Eurosilicone via the e-mail address: productsupport@gcaesthetics.com and an individual treatment should be initiated.

ASYMMETRY

Asymmetry may result from improper initial placement, displacement, or failure to correct pre-existing asymmetry through variation in individual implant size. Pre-existing asymmetry may not be entirely correctable. Asymmetry may also be a symptom of capsular contracture, fluid accumulation, infection, post-operative breast dysplasia, unilateral discrepancy in muscle development, or deflation, requiring further investigation.

BREAST IMPLANT ASSOCIATED ANAPLASTIC LARGE CELL LYMPHOMA (BIA-ALCL) (and other lymphomas)

European safety information, the US FDA and current scientific literature have identified an association between breast implants and the development of a cancer of the immune system called breast implant-associated anaplastic large cell lymphoma (BIA-ALCL). This specific entity is included in the WHO 2016 classification under the terminology "ALCL – AIM". This cancer occurs more commonly in patients with textured breast implants rather than smooth implants, and typically develops many years after breast implant surgery. At this time rates of BIA-ALCL are considered to be low.

In most patients, BIA-ALCL is found within the fluid or scar tissue next to the implant. The main symptoms of BIA-ALCL are persistent swelling or pain, late onset persistent seroma, capsular contracture, or a mass in the vicinity of the breast implant. If women present with any of the aforementioned symptoms, or if BIA-ALCL is suspected, it is advised that patients should be assessed to rule out peri-implant ALCL. When testing for BIA-ALCL, collect fresh seroma fluid and representative portions of the capsule and send for pathology tests to rule out ALCL. Diagnostic evaluation should include cytological evaluation of seroma fluid with Wright Giemsa stained smears and cell block immunohistochemistry testing for cluster of differentiation (CD) and Anaplastic Lymphoma Kinase (ALK) markers. Most cases of breast implant-associated ALCL are treated by removal of the implant and the capsule surrounding the implant and some cases have been treated by chemotherapy and radiation. Although treatment is usually successful BIA-ALCL can be fatal if left untreated, hence, early detection and treatment is key.

WARNING: AS BIA-ALCL OCCURRENCE IS A RESIDUAL RISK OF THE IMPLANT: In the case that your patient is diagnosed with BIA-ALCL, a report should be made to Eurosilicone via the e-mail address: productsupport@gcaesthetics.com and an individual treatment program should be initiated by a multidisciplinary team.

BIA-SCC

Breast implant-associated squamous cell carcinoma is a rare complication of breast implantation that can result in significant morbidity and mortality. BIA-SCC is a rare, epithelial-based tumor distinct from BIA-ALCL. Patients can present with unilateral breast seroma, fluid collection, pain, erythema, and capsular contracture (Baker Grade IV).

Initial treatment typically involved removal of the implant and the capsule surrounding the implant and some cases have been treated by chemotherapy with limited success. Aggressive first-stage surgery seems to have a better prognosis and outcome. BIA-SCC can be fatal, early detection and aggressive treatment is key.

In the case that your patient is diagnosed with BIA-SCC, a report should be made to Eurosilicone via the e-mail address: productsupport@gcaesthetics.com and an individual treatment program should be initiated by a multidisciplinary team.

BREAST IMPLANT ILLNESS

A small proportion of women who have breast implants, for both aesthetic and reconstructive purposes, self-identify as having a number of symptoms

they believe arise from the presence of their breast implants. Whilst not a medical diagnosis, they refer to their symptoms as 'Breast Implant Illness' (BII).

These multiple generic symptoms are varied, and all cases have been self-diagnosed and self-reported. They include, but are not limited to, flu-like symptoms such as extreme tiredness, brain fog, joint aches, immune-related symptoms, sleep disturbance, depression, hormonal issues, headaches, hair loss and chills.

There are a variety of other reasons these symptoms might be experienced, including background illnesses or hormonal changes. In addition, there are a number of scientific studies investigating similar symptoms experienced by women in the general population with and without breast implants. On average, around 50% of women who self-identify as having 'BII' feel that their symptoms improve after removal – sometimes temporarily and sometimes permanently. It therefore appears that removing breast implants does not necessarily improve symptoms in everyone. To date, there is no research demonstrating which symptoms may or may not improve with implant removal.

Currently, there are no tests that can confirm 'BII'. Research continues in this area, specifically around patients who have autoimmune disease or a predisposition to autoimmune disease. Breast implantation is contraindicated in women with a suppressed or compromised immune system – see CONTRAINDICATIONS section above. Medical advice should be sought if patients suspect they have BII. Their symptoms may not be related to the implants and other medical investigations should not be overlooked or ignored in determining the cause of their symptoms.

BREAST FEEDING

The implantation may affect the milk production during nursing, particularly, peri-areolar incision may reduce considerably the possibility of breast feeding.

CALCIUM DEPOSITS

Calcification commonly occurs in mature breast tissue with or without implantation. Calcification is also known to occur after implantation of a foreign body, although the etiology is unknown, and reported cases are rare. Microcalcification after implantation typically occurs on or around the fibrous capsule in thin plaques or accumulations. Extensive microcalcifications may cause breast hardness and discomfort and may necessitate surgical intervention.

CAPSULAR CONTRACTURE

Fibrous capsular contracture is a common complication following breast implant surgery. While formation of a fibrous tissue capsule surrounding the implant is a normal physiological response to a foreign body, not all capsules contract. Contracture of the fibrous capsular tissue surrounding the implant may result in firmness, discomfort or pain in the breast, distortion of the breast, palpability of the implant, or displacement of the implant. The etiology of capsular contracture is unknown but is most likely multifactorial. Contracture develops to varying degrees, unilaterally or bilaterally, and may occur within weeks to years after surgery. Severe cases may require surgical intervention. Capsular contracture may recur after capsulotomy or capsulectomy. Lumps perceived as capsular scarring have delayed the diagnosis of palpable tumor. Questionable lumps must be promptly investigated.

COMPARTMENT SYNDROME

Compartment syndrome is a very rare but serious complication that can be linked to any type of surgery. This causes an increase in pressure in the upper or lower limbs that needs to be released thanks to a surgical intervention consisting of compartment-releasing fasciotomy to reduce or possibly prevent acute and long-term complications.

DELAYED WOUND HEALING

Some patients experience delayed healing, and for others the incision site may not heal well. It may open from injury or infection. If the implant is exposed, further surgery will be required. Tissue necrosis is the development of dead tissue around the implant. It will delay wound healing, may cause wound infection and may require surgical correction and/or implant removal. Tissue necrosis has been reported following the use of steroid drugs, chemotherapy, radiation to breast tissue, and smoking, but in some cases, it may occur without any known cause.

DISSATISFACTION WITH RESULTS

The complications of incorrect size, misplaced scar location, hypertrophic scarring and those listed below are usually related to surgical technique. Careful size selection, creation of an appropriate and adequate size surgical pocket, and use of current accepted surgical procedures are the surgeon's responsibility.

DISPLACEMENT OF IMPLANT AND MALROTATION

Implants may displace or rotate with accompanying discomfort and/or distortion in breast shape. Difficult placement techniques may increase the risk of displacement by reducing pocket size and placement accuracy. Displacement may require surgical intervention.

EXPLANATION OF MAMMARY IMPLANTS

Although there are no definite reports in the medical literature, it is anticipated that a large number of mammary implants will be surgically explanted or replaced. This is not unlike many other plastic surgery procedures which are commonly repeated to maintain patient satisfaction. Due to the wide variety of patients' physical responses to breast implant surgery, differences in surgical techniques and medical treatments, as well as potential complications, patients should be advised that these should not be considered lifetime implants and explant surgery may be indicated at any time. Eurosilicone makes no representations for the term of implantation of the device.

EXTERNAL CAPSULOTOMY

Treatment of capsular contracture by external manual compression may cause the shell to weaken or rupture. There have been reports of rupture and Eurosilicone recommends against the external capsulotomy procedure and is not responsible for the integrity of the implant should manual compression or other external stress techniques be used.

EXTRUSION

Unstable or compromised tissue covering and/or interruption of wound healing may result in exposure and extrusion of the implant. Causes or contributing factors may include infection, wound dehiscence, necrosis with or without infection, capsular contracture, closed capsulotomy, unsuitable skin flap, improper size, and placement of implant, and/or tissue erosion associated with implant folds. The incidence of extrusion has been shown to increase when the prosthesis has been placed in injured areas: scarred, heavily irradiated, or burned tissue or crushed bone areas; where severe surgical reduction of the area has been performed; and where steroids are used in the surgical pocket.

GRANULOMA

Granuloma formation is a common tissue response to the presence of foreign materials. It is possible for a tissue reaction to be caused by the presence of the implant leading to formation of silicone granuloma.

HEMATOMA/ EARLY ONSET SEROMA

Hematoma and serous fluid accumulation are complications associated with any type of invasive surgery. Postoperative hematoma and seroma may contribute to infection and/or capsular contracture. Postoperative hematoma and seroma may be minimized by meticulous attention to hemostasis during surgery, and also possibly by postoperative use of a closed drainage system. Persistent, excessive bleeding must be controlled before the device is implanted. Any postoperative evacuation of hematoma or seroma must be conducted with care to avoid contamination or damage to the implant.

HYPERPIGMENTATION

Scar hyperpigmentation is due to melanocyte absence or excessive presence in the scar tissue. The melanogenesis may not happen properly. This adverse event is diagnosable after the full cicatrization of the wound (six to twelve months) and may be irreversible.

INFECTION

Pre-existing infection not resolved before implant placement increases the risk of periprosthetic infection. Do not expose the implant or filling accessories to contaminants, which increases the risk of infection. Infection is an inherent risk following any type of invasive surgery. Infection around a breast implant may occur within days, weeks, or even years, after surgery. Signs of acute infection reported in association with breast implants include erythema, tenderness, fluid accumulation, pain, and fever. Signs of subclinical infection may be difficult to detect. Postoperative infections should be treated aggressively according to standard medical practices to avoid more serious complications. Infection that is unresponsive to treatment or necrotizing infection may require implant removal. Capsular contracture may be related to infection in the area surrounding the implant.

PAIN

Any unexplainable pain, not associated with any of the complications listed herein, must be promptly investigated.

PHLEBITIS SUPERFICIAL

In rare cases, damage to the vein walls may occur during implantation of the mammary prosthesis. Subsequently, this may result in phlebitis superficial, which is an inflammation of the veins caused by a blood clot below the surface of the skin. This complication is rarely serious and with proper care, can be resolved rapidly.

PTOSIS

Ptosis, also known as sagging of the breast, is a deformity which can occur pre- and post- breast surgery. Pre-breast surgery, due to the breast having no bony or cartilaginous support, time and gravity can affect the natural breast including glandular and skin stretching. Post-breast surgery, the size, weight, and position of breast implants can have an effect on the degree of ptosis. If significant, this complication can often lead to re-intervention.

RUPTURE OF GEL IMPLANTS

Gel rupture can occur post-operatively from damage to the implant during handling or surgery. Rupture of the shell can also occur from contracture, trauma, or excessive manipulation. Ruptures of unknown etiology have also been reported. Because of the cohesiveness of the gel material, some gel ruptures can go undetected unless surgery occurs for another reason (e.g. size exchange). Despite the cohesive properties of the gel, extravasation out of the surgical pocket can occur under pressure, in which case additional surgery may be required to retrieve the gel. Inflammation and formation of silicone granulomas have been reported. If shell rupture is suspected, the implant should be removed. Some silicone gel may diffuse or "bleed" through the shell of an intact implant. Implant shell and silicone gel improvements have reduced this bleeding to a minimum.

SENSATION

The risk of temporary or permanent dysesthesia exists following any invasive surgical procedure. Careful surgical technique can minimize, but not preclude this risk. Dysesthesia of the nipple/areolar complex, and less frequently, the breast area in general, has been reported following implantation and may be temporary or permanent. The risk of neurological impairment increases with more extensive surgery. Surgical implantation of a mammary prosthesis may interfere with the ability to breast feed. However, it should be noted that previous breast reconstruction surgery, may be the initial cause of this interference.

SEROMA (late)

A seroma is an accumulation of serous fluid around the implant. It is categorized late one year after the surgery. A seroma can usually be felt as a mass under the skin. The area is swollen, it may cause erythema (red skin) and pain. If you have a late seroma, please immediately refer to your surgeon.

SILICONE ALLERGY

Although silicone allergy is extremely rare, this risk cannot be completely ruled out.

SUPERFICIAL WOUND

Superficial wounds result from incisions made on the skin that are required for the insertion of implants during surgery. Treatment of superficial wounds can include cleaning the area with saline solution (salt water) or mild soapy water. After cleaning, gently pat the wound dry using a clean, soft cloth or piece of gauze.

Superficial wounds are minor injuries, but if not properly managed, they can potentially lead to complications such as infection due to an increased risk of bacterial contamination, delayed healing, or undesirable scarring.

SILICONOMA

Siliconoma is a foreign body granulomatous inflammation with dense infiltrate consisting of plasma cells, lymphocytes and histiocytes. Siliconomas can occur locally or present at a distant site due to its high fat solubility.

UPPER LIMB LYMPHEDEMA

Lymphedema is a chronic condition characterized by swelling of the arm, hand, or breast, which can result from breast cancer treatment, due to the accumulation of lymphatic fluid in the interstitial tissues. This condition is known to significantly compromise quality of life due to body image changes, alterations in arm function, and increase the risk of other complications including infection and cellulitis.

Recent studies suggest immediate breast reconstruction may be associated with a reduced risk of post-mastectomy lymphedema. Treatment options include conservative management or surgical treatment using a super-microsurgical technique. However, no complete cure for lymphedema is known.

WRINKLING AND FOLDS

Thin or inadequate overlying tissue, patients with little or no subcutaneous fat, implants which are too large for the surgical pocket, or the anatomical structure of the patient, and subcutaneous placement may contribute to palpable, or visible, wrinkles and folds. Folds may result in thinning and erosion of adjacent tissue and erosion of the implant. Signs of skin inflammation, such as tenderness and erythema, may indicate thinning or erosion and must be promptly investigated. Palpable wrinkling / and or folds may be confused with palpable tumor, and questionable cases must be promptly investigated.

INSTRUCTIONS FOR USE

STORAGE CONDITIONS

Eurosilicone products must be stored in normal conditions, must be protected from humidity and direct sunlight and must be stored with the arrows pointing up. In these conditions, Eurosilicone products have a shelf life of five years.

SUPPLIED STERILE

Implants are supplied in sterile form (Ethylene Oxide Sterilization or Dry Heat Sterilization, corresponding symbol on product label), processed by validated strictly controlled sterilization cycles with ethylene oxide or dry heat. Sterility is verified in accordance with standards. Sterility of the implant is maintained only if the package is intact and undamaged.

If the sterile packaging has been damaged or unintentionally opened before use do not use.

SINGLE USE

Implants are intended for SINGLE USE ONLY. DO NOT REUSE EXPLANTED PRODUCTS. DO NOT RESTERILIZE ANY PRODUCT. Explanted products should not be reused because recleaning and re-sterilization procedures may not adequately remove biological residues, such as blood, tissue, and other matter, which could retain resistant pathogens and may also affect the performance of the implant.

PACKAGING

Sterile product is supplied in a sealed, double primary package. Sterility is not guaranteed if the package has been damaged or opened. Tear-Off patient record labels are attached to the primary package. These labels are to be attached to the patients' records. The implant card label intended to be placed on the implant card is also provided on the primary package.

TO OPEN PACKAGED STERILE PRODUCT

- Peel open outer package under clean, aseptic conditions, over sterile field, allowing sealed inner package to fall gently into the field.
- Attach patient record portion of the inner label to patient's chart.
- Using aseptic precautions, peel open the inner package.

PRODUCT EXAMINATION AND HANDLING

- Product should be visually examined for any evidence of particulate contamination, damage, or leakage.
- The implant should be kept submerged in sterile water or normal saline prior to implantation to prevent contact with airborne and surgical field particulate contaminants.
- Do not immerse the implant in Betadine or Iodine containing solutions. If Betadine or Iodine containing solutions are used in the pocket, ensure that it is rinsed thoroughly so no residual solution remains in the pocket.
- Gel-filled mammary implants may contain air bubbles. This occurs normally during sterilization and has no effect on product integrity or performance.
- Use disposable, non-shedding wipes only to ensure there are no particles left on the surface of the implant.

CONTAMINATION OF IMPLANTS

Care must be taken to prevent surface contaminants such as talc, dust, and skin oils from coming into contact with the implant. Products should be inspected for contamination prior to insertion. Contamination at the time of surgery increases the risk of periprosthetic infection, and possibly capsular contracture. Textured silicone surfaces are potentially more susceptible to contamination than smooth silicone surfaces. Extra care should be taken in the handling of textured surface implants. Eurosilicone accepts no liability for products contaminated by other substances after the product leaves our possession.

SURFACE CONTAMINANTS (TALCS, DUST, LINT, OILS) ON SURFACE OF IMPLANTS CAN CAUSE FOREIGN BODY REACTION, HANDLE WITH CARE WITH SURGICAL GLOVES (RINSED FREE OF TALC) WITH STRICT ASEPTIC TECHNIQUE. DO NOT IMPLANT CONTAMINATED PRODUCT.

Back-up implants must be readily available at the time of surgery for use in the event of contamination.

SURGICAL PROCEDURE

Correct surgical procedures and techniques are the responsibility of the medical profession. Each surgeon must evaluate the suitability of the procedure based upon current accepted techniques, individual judgement, and experience. Proper size and shape of implants must be determined for the individual patient by the surgeon. An incision should be of appropriate length to accommodate the style, size, and profile of the implant. This will reduce the potential for creating excessive stress to the implant during insertion.

SPECIFIC GUIDELINES FOR THE USE OF ROUND EUROSILICONE MAMMARY IMPLANTS

- Size selection should be decided by base dimensions, the patient's dimensions, and individual breast characteristics rather than preferred volume or projection.
- The preformed shapes of the Eurosilicone implants are firmer than other cohesive gel filled implants, and the surgical incisions should be longer than those required for compressible implants.
- Where appropriate the new inframammary fold may be required to be approximately 1cm below the natural fold.
- The recommended approach is via the inframammary fold. It is the responsibility of the surgeon to decide if axillary or trans areolar placements may be possible. We caution that insertion and orientation may be more difficult due to the increased firmness of the gel.
- Sub-muscular placement (with possible lower ligament release) is preferred where possible. As the gel is highly cohesive it has a more defined radius on the upper apex margin. Sub-muscular placement reduces the probability of this apex radius being palpable.
- DO NOT use force during placement of Eurosilicone implants. Over manipulation of Eurosilicone implants may result in permanent distortion of the cohesive gel and loss of the implant's shape
- Maintenance of correct alignment and avoidance of rotation in the post-operative phase is aided by:
- Achieving a surgical pocket closely aligned to the implant.
- Achieving a "dry" surgical pocket.
- Practicing the procedure of implant stabilization by use of post-operative external immobilization.

MATERIALS DEVICE DESCRIPTION

These medical devices are manufactured from medical grade silicone materials which are suitable for long term implantation. Master Files have been filed with the U.S FDA.

The following tables provide quantitative and qualitative information on the materials and substances which patients can be exposed based on chemical characterization of representative devices. The materials and device have been subjected to biocompatibility testing and evaluation, and risk assessments to demonstrate their biological safety. However, individual responses to chemicals may vary, and all reactions cannot be predicted.

Table 1. Round Breast Implant Device Materials

Device materials	Implant component
Silicone dispersions	Shell
Silicone adhesive	Shell
Silicone gel	Gel fill
Silicone dispersions	Closure Patch
Platinum catalyst	Shell and Gel fill
Organotin catalyst	Shell (silicone adhesive)

Table 2. Chemicals Released by Round Breast Implants

Compounds	Whole device (µg/g or ppm)	
	Round Collection Textured Cristalline Paragel	Round Collection Smooth Paragel
Volatiles¹ - Chemicals that can be released by round breast implants as a gas		
Benzene	Not Detected	Not Detected
Bromobenzene	Not Detected	Not Detected
Bromodichloromethane	Not Detected	Not Detected
Bromoform	Not Detected	Not Detected
Butylbenzene (N- / tert- / sec-)	Not Detected	Not Detected
Carbon tetrachloride	Not Detected	Not Detected
Chlorobenzene	Not Detected	Not Detected
Chloroform	Not Detected	Not Detected

Compounds	Whole device (µg/g or ppm)	
2,4-Chlorotoluene	Not Detected	Not Detected
Dibromochloromethane	Not Detected	Not Detected
1,2-Dibromo-3-chloropropane	Not Detected	Not Detected
Dichlorobenzene (1,2- / 1,3- / 1,4-)	Not Detected	Not Detected
Dichloroethane (1,1- / 1,2-)	Not Detected	Not Detected
Dichloroethene (1,1- / cis-1,2- / trans-1,2-)	Not Detected	Not Detected
Dichloropropane (1,2- / 1,3-)	Not Detected	Not Detected
1,1-Dichloropropene	Not Detected	Not Detected
Ethanol, 2-trimethylsilyl	1.25 max.	0.75 max.
Ethylbenzene	Not Detected	Not Detected
Hexachlorobutadiene	Not Detected	Not Detected
Isopropylbenzene	Not Detected	Not Detected
p-Isopropyltoluene	Not Detected	Not Detected
N-propylbenzene	Not Detected	Not Detected
Naphthalene	Not Detected	Not Detected
Styrene	Not Detected	Not Detected
1,1,1,2-Tetrachloroethane	Not Detected	Not Detected
Tetrachloroethene	Not Detected	Not Detected
Trichloroethene	Not Detected	Not Detected
Trichloroethane (1,1,1- / 1,1,2-)	Not Detected	Not Detected
1,2,3-Trichloropropane	Not Detected	Not Detected
Trimethylbenzenes (1,2,4-/1,3,5-)	Not Detected	Not Detected
Toluene	Not Detected	Not Detected
Trichlorobenzenes (1,2,4- / 1,2,3-)	Not Detected	Not Detected
o-Xylene	Not Detected	Not Detected
m- & p-Xylene	Not Detected	Not Detected
Total volatiles	1.25	0.75

Compounds (for Round Collection – Cristalline Paragel)	Whole device (µg/g or ppm)		Compounds (for Round Collection Smooth – Paragel)
	Round Collection Textured Cristalline Paragel	Round Collection Smooth Paragel	
Extractables ² - Chemicals that can be released by round breast implants following soaking in water			
Cyclic siloxanes (D4, D5, D6...)	Not Detected	Not Detected	Semi-volatile organic compound
Linear siloxanes (L3, L4, L5...)	Not Detected	0.91	Non-volatile organic compound
Total extractables	< 5.0	< 7.2	
¹ Not Detected means that the level of the individual volatile was below the quantitation limit of the test method. The quantitation limit was 0.50 µg/g for m- & p-Xylene and 2,4-Chlorotoluene, and 0.26 µg/g for all other compounds.		¹ Not Detected means that the level of the individual volatile was below the quantitation limit of the test method. The quantitation limit was 0.72 µg/g for m- & p-Xylene and 2,4-Chlorotoluene, and 0.37 µg/g for all other compounds.	
² Not Detected means that the level of the individual extractable was below 5.0 µg/g, the highest quantitation limit of the test method.		² Not Detected means that the level of the individual semi-volatile organic compound was below 7.2 µg/g (LOQ of the GC-MS method).	

Table 3. Heavy Metals Found in Round Breast Implants

Heavy metals	Concentration (ppm)	
	Round Collection Textured Cristalline Paragel	Round Collection Smooth Paragel
Arsenic	Not Detected	Not Detected
Cadmium	Not Detected	Not Detected
Lead	Not Detected	Not Detected
Mercury	10 max.	10 max.
Vanadium	Not Detected	Not Detected
Molybdenum	Not Detected	Not Detected
Selenium	Not Detected	Not Detected
Cobalt	Not Detected	Not Detected
Antimony	Not Detected	Not Detected
Barium	Not Detected	Not Detected
Chromium	Not Detected	Not Detected
Copper	Not Detected	Not Detected
Tin	Not Detected	Not Detected
Nickel	Not Detected	Not Detected
Platinum	Not Detected	Not Detected
Zinc	0.05 max.	
Magnesium	Not Detected	Not Detected
Manganese	Not Detected	Not Detected
Thallium	Not Detected	Not Detected
Not Detected means that the level of the individual element was below 0.05 ppm, the quantitation limit of the test method.		
Not Detected means that the level of the individual element was below 0.07 ppm, the quantitation limit of the test method.		

DISPOSAL OF DEVICE

All mammary implants must be eliminated according to requirements related to medical wastes involving an infectious risk.

REPORTING OF SERIOUS INCIDENT:

Any serious incident (incident that directly or indirectly led, might have led or might lead to any of the following: (a) the death of a patient, user or other person, (b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health, (c) a serious public health threat) should be reported to Eurosilicone via the e-mail address: productsupport@gcaesthetics.com and to the competent authority of the Member State in which the user and/or patient is established.

RETURN OF GOODS POLICY

Prior to the return of any product:

- Authorization must be requested in writing.
- Decontamination certificate must be supplied with any product returned which has been removed from the packaging.
- If it is necessary to puncture the shell of any sealed shell product to assist with safe decontamination sterilization procedures, the area of mechanical interference should be marked with indelible marker on the surface of the product and reference made on the decontamination certificate.
- Products File Notes and returned products must be sent to Eurosilicone S.A.S. by your distributor only.

WARRANTY

The manufacturer warrants that reasonable care was used in the manufacture of these products and will replace any product that the manufacturer's investigation shows to be defective at time of shipment. Patient selection, surgical procedures, post-surgical treatment and stresses, and handling of the devices are totally the responsibility of the customer. Eurosilicone has no control over the conditions of use and cannot warrant good effect or against ill effect following the use of the device and shall not be responsible for any incidental or consequential loss or damage or expenses directly or indirectly arising from the use of this device. All other warranties (whether implied by law or otherwise) are excluded to the fullest extent permitted by law.

CAUTION:

Federal (USA) law restricts this device to sale by or on the order of a physician.

HARMONISED STANDARD AND COMMON SPECIFICATIONS APPLIED

The following common specification and relevant harmonized standards are applied for all mammary implants: Common specification 2022/2346; EN ISO 13485; EN ISO 14971; EN ISO 62366-1; EN ISO 14630; EN ISO 14607; EN ISO 10993 series; EN ISO 15223-1; EN ISO 20417; EN ISO 11607; EN ISO 11135 or EN ISO 20857 (depending sterilization method).

BIBLIOGRAPHIC REFERENCES

Literature references are available upon request from Eurosilicone.

LABELLING SYMBOLS:

Indicates the products dimensions.

„Round Collection is a silicone gel pre-filled mammary implant intended to be used for both non-medical (cosmetic) and medical purposes. It is not intended to be used in persons who are less than 18 years old. Only to be implanted in an appropriate medical environment by appropriately trained medical doctors who are qualified or accredited in accordance with national law”.

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- gaminių failo pastabas ir gražinamus gaminius „Eurosilicone S.A.S.“ gali išsiųsti tik jūsų platintojas.
Eurosilicone SAS