


Solus™ Laryngeal mask airways

Solus™

Due to growing concerns regarding the risk of cross infection in hospitals, Intersurgical offer the Solus single use laryngeal mask airway which is ideal for use both in anaesthesia and emergency medicine. The Solus compliments the range of Intersurgical products and allows us to provide the option of a complete respiratory and airway management package.



The Solus™

The Solus is now available in seven different product sizes from adult to neonate.

Flexibility and strength

The materials used in the manufacture of the Solus have been specifically chosen to provide the product with the quality, flexibility and strength required to work effectively within a range of demanding clinical applications.

Convenience and safety

The Solus range has been designed with convenience and safety in mind. Because Solus is a single use product it requires only two simple pre-use checks: Observation of the airway tube to inspect for occlusion or foreign objects and cuff inflation to check for patency. The Solus is supplied sterile so it is a convenient ready-to-use product.



The Solus is designed to be disposed of immediately after use so this removes the need for any cleaning processes.







Latex Free.

Why single use?

Some reusable laryngeal mask airways can go through an estimated 800 process handlings during their recommended 40-use lifespan; Solus undergoes just 8 process handlings, including disposal. When cross infection is a major concern⁽¹⁾ this offers a distinct advantage: typically Solus will be handled by only three individuals during its use, whilst reusable devices can be handled by any number of individuals.






Adult laryngeal mask airways



 <p>Laryngeal mask airway - Size 5</p> <p>Code: 8005 Box Quantity: 20</p>	 <p>>70kg <40ml <60cm H₂O</p>
 <p>Laryngeal mask airway - Size 4</p> <p>Code: 8004 Box Quantity: 20</p>	 <p>50-70kg <30ml <60cm H₂O</p>




SOLUS™ Laryngeal mask airway



Adult laryngeal mask airways (continued)




			30-50kg	
Laryngeal mask airway - Size 3				
Code: 8003	Box Quantity: 20			

Paediatric laryngeal mask airways

			20-30kg	
Laryngeal mask airway - Size 2 1/2				
Code: 8025	Box Quantity: 20			

			10-20kg	
Laryngeal mask airway - Size 2				
Code: 8002	Box Quantity: 20			

			5-10kg	
Laryngeal mask airway - Size 1 1/2				
Code: 8015	Box Quantity: 20			

			<5kg	
Laryngeal mask airway - Size 1				
Code: 8001	Box Quantity: 20			

References:

(1) **Biological contaminants identified from laryngeal mask airways** - LUJR Professor and head of Biological physics, UMIST Manchester. *ISSM Journal Sept 2004 Vol 9 N°1*

IS6.1 KW Issue 4 06.05

Solus™

Laryngeal mask airway

A complete range of single use devices
for adults through to neonates

