

### M3001A A01

Parameters: 3 or 5 Lead ECG, Philips Fast SpO<sub>2</sub>, NiBP (Non-Invasive Blood Pressure)



### M3001A A02 (Nellcor Oximax)

Parameters: 3 or 5 Lead ECG, Nellcor Oximax SpO<sub>2</sub>, NiBP (Non-Invasive Blood Pressure)



### M3001A A03 (Masimo)

Parameters: 3 or 5 Lead ECG, Masimo SpO<sub>2</sub>, NiBP (Non-Invasive Blood Pressure)



### M3001A A04 (Nellcor Oximax)

Parameters: 3 or 5 Lead ECG, Nellcor Oximax SpO<sub>2</sub>, NiBP (Non-Invasive Blood Pressure)



### M3001A A01C06

Parameters: 3 or 5 Lead ECG, Philips Fast SpO<sub>2</sub>, NiBP (Non-Invasive Blood Pressure), Temperature, iBP (Invasive Blood Pressure)



### M3001A A02C06 (Nellcor Oximax)

Parameters: 3 or 5 Lead ECG, Nellcor Oximax SpO<sub>2</sub>, NiBP (Non-Invasive Blood Pressure), Temperature, iBP (Invasive Blood Pressure)



### M3001A A03C06 (Masimo)

Parameters: 3 or 5 Lead ECG, Masimo SpO<sub>2</sub>, NiBP (Non-Invasive Blood Pressure), Temperature, iBP (Invasive Blood Pressure)



### M3001A A04C06 (Nellcor Oximax)

Parameters: 3 or 5 Lead ECG, Nellcor Oximax SpO<sub>2</sub>, NiBP (Non-Invasive Blood Pressure), Temperature, iBP (Invasive Blood Pressure)



### M3001A A01C12

Parameters: 3, 5 or 12 Lead ECG, Philips Fast SpO<sub>2</sub>, NiBP (Non-Invasive Blood Pressure)



### M3001A A02C12 (Nellcor Oximax)

Parameters: 3, 5 or 12 Lead ECG, Nellcor Oximax SpO<sub>2</sub>, NiBP (Non-Invasive Blood Pressure)



### M3001A A03C12 (Masimo)

Parameters: 3, 5 or 12 Lead ECG, Masimo SpO<sub>2</sub>, NiBP (Non-Invasive Blood Pressure)



### M3001A A04C12 (Nellcor Oximax)

Parameters: 3, 5 or 12 Lead ECG, Nellcor Oximax SpO<sub>2</sub>, NiBP (Non-Invasive Blood Pressure)



### M3001A A01C18

Parameters: 3, 5 or 12 Lead ECG, Philips Fast SpO<sub>2</sub>, NiBP (Non-Invasive Blood Pressure), Temperature, iBP (Invasive Blood Pressure)



### M3001A A02C18 (Nellcor Oximax)

Parameters: 3, 5 or 12 Lead ECG, Nellcor Oximax SpO<sub>2</sub>, NiBP (Non-Invasive Blood Pressure), Temperature, iBP (Invasive Blood Pressure)



### M3001A A03C18 (Masimo)

Parameters: 3, 5 or 12 Lead ECG, Masimo SpO<sub>2</sub>, NiBP (Non-Invasive Blood Pressure), Temperature, iBP (Invasive Blood Pressure)



### M3001A A04C18 (Nellcor Oximax)

Parameters: 3, 5 or 12 Lead ECG, Nellcor Oximax SpO<sub>2</sub>, NiBP (Non-Invasive Blood Pressure), Temperature, iBP (Invasive Blood Pressure)

