

- placed over the head set tubing
- ☐ Cushions are installed during the assembly
- ☐ The *TLC* cushions move easily back and forth on the head set tubes
- ☐ Unique, patented, anatomically correct face piece design
- ☐ *TLCannulas* are available in a wide selection of popular Salter-Style cannulas
- ☐ Now available in a dual port demand type cannula

- of nonscheduled home visits
- ☐ No assembly is required, means time savings for the care provider and the cushions will not fall off
- ☐ Easy for patient to reposition for the most comfortable and secure fit
- ☐ Provides a more "natural" fit, and reduces patient anxiety
- ☐ Easy to find a comfortable fit and reduces unnecessay time consuming service visits
- ☐ Provides maximum comfort for the LTOT patient using either a dual port, portable or stationary conserving device



TLCannula[™] is now available in the following popular models of Salter-Style[®] Cannulas

1600TLC...A unique, patented design which sets the standard for comfort. All *Salter-Style* cannulas are anatomically shaped with soft, curved nasal prongs.

1606TLC...All the comfort features of the 1600 model, connecting tubing is not included.

1600QTLC...All the comfort features of the 1600 model. However, the curved face piece and nasal prongs are larger to accommodate higher flows (up to 8 LPM). It is considerably quieter and provides a more soothing environment for the patient.

1616TLC...Micro... all the comfort features of the 1600 model. This version uses a smaller molded face piece, nasal prongs and cannula tubes. It is lighter in weight, less conspicuous and is intended for low flow (up to 3 LPM) applications.

4804TLC and 4807TLC...Demand TLCannulas for use with dual port oxygen conservation delivery systems. A permanent barrier in the cannula's face piece allows one prong to sense inhalation/exhalation while oxygen is delivered through the other prong.

TLCannula

Gray colored, soft foam ear cushions are pre-assembled over the cannula headset tubing for comfort and convenience.



All *TLCannulas* incorporate the patented features common to all *Salter-Style* cannulas. Over the years our cannulas have earned customer acclaim, as an industry standard for long term "wearability comfort". The gray colored, closed-cell foam cushions on *TLC* models are permanently installed on the cannula tubing. The foam cushions slide easily on the head set tubing to achieve the most comfortable fit and to diminish chafing and tissue injury over the ear. An added benefit of *TLCannulas* is a reduction of unnecessary home visits to your more demanding customers.

	Ordering Information Salter-Style® TLCannula™	Units (per case)	Part#
	Adult, TLCannula, clear with foam cushions and 7' supply tube	25	1600TLC*
	Adult, TLCannula, clear with foam cushions and 1' supply tube	25	1600TLC-1*
	Adult, TLCannula, clear with foam cushions and 4' supply tube	25	1600TLC-4*
	Adult, TLCannula, clear with foam cushions and 14' supply tube	25	1600TLC-14*
	Adult, TLCannula, clear with foam cushions and 15' supply tube	25	1600TLC-15*
	Adult, TLCannula, clear with foam cushions and 25' supply tube	25	1600TLC-25*
	Adult, Quiet <i>TLCannula</i> , clear large face piece, foam cushions and 7' supply tube (for oxygen flows up to 8 LPM)	25	1600QTLC*
	Small Adult/Pediatric, TLCannula, clear with foam cushions and 7' supply tube	25	1602TLC
	Adult, TLCannula, clear with foam cushions, no supply tube	25	1606TLC
	Adult, Micro <i>TLCannula</i> , clear, foam cushions, and 7' supply tub (for oxygen flows up to 3 LPM)	e 25	1616TLC
	Adult, Micro <i>TLCannula</i> , clear, foam cushions, and 4' supply tub (for oxygen flows up to 3 LPM)	e 25	1616TLC-4
	Adult, TLCannula, clear with female thread grip connector and 7' supply tube	25	1699TLC
E W	Adult , Demand <i>TLCannula</i> , clear with foam cushions and 4' smooth bore supply tube	10 or 25	4804TLC**
EW	Adult, Demand <i>TLCannula</i> , clear with foam cushions and 7' smooth bore supply tube	10 or 25	4807TLC**
	Adult , Conventional Style <i>TLCannula</i> , clear with foam cushions and 7' supply tube	25	1053TLC*
	Adult , Conventional Style <i>TLCannula</i> , clear with foam cushions and no supply tube	25	1057TLC

The closed cell foam cushions are available as separate items for self-installation on nasal cannulas.

(see SLC-71 Salter Accessories E-Z Wrap^m P/N 1016).



 $Printed in U.S.A @ Copyright 2004, Salter Labs^*, *U.S. Patent No. D369,409 ** U.S. Patent No. 5,137,017. Other U.S. and worldwide patents pending. Revised April, 2004 SLC-00$