



	Avante Wildcat 750 Patient Warmer	3M Bair Hugger Patient Warmer	Stryker Mistral Air Patient Warmer
Dimensions, Height & Weight	H 11.0in x D 8.7in x W 8.7in 11.9lb	H 13in x D 15in x W 15.5in 16lb	H 15.7in x D 6in x W 13.4in 11.5lb
Calibration	Digital Calibration - No need to calibrate	Analog Calibration - Must calibrate every 6 months for accurate temperature	Calibration Device needed to Calibrate
Noise Level	High Level (<52 dB) Low Level (<40dB)	High Level (<53 dB) Low Level(<48dB)	51 dB
Filters	HEPA, 99.99%, 0.2 µm high efficiency filter, 2,000 hours	Non-HEPA, 99.7%, 0.3 µm high efficiency filter, 500 hours	HEPA, 99.99%, 0.2 µm high efficiency filter, 2,000 hours
Storage Cart	Height Adjustable Storage Cart	Fixed Height Storage Cart	Storage Cart connected to pole clamp
Sterilized Blankets	Yes	Only Cardiac	Only Cardiac and Surgical Access
Blanket Inlets	Up to 4 inlets	1 hose inlet	Up to 2 inlets
Machines needed to Heat Gowns	1	2	1