



|                                      | Avante Wildcat 750   | 3M Bair Hugger   | Stryker Mistral Air  |
|--------------------------------------|--|--|--|
| <b>Dimensions &amp; Weight</b>       | H 11.0in x D 8.7in x W 8.7in<br>11.9lb                         | H 13in x D 15in x W 15.5in<br>16lb   | H 15.7in x D 6in x W 13.4in<br>11.5lb                          |
| <b>Calibration</b>                   | Digital Calibration:<br>No need to calibrate                   | Analog Calibration:<br>Must calibrate every 6 months<br>for accurate temperature | Calibration Device<br>needed to Calibrate                      |
| <b>Noise Level</b>                   | High Level (<52 dB)<br>Low Level (<40dB)                       | High Level (<53 dB)<br>Low Level(<48dB)  | 51 dB  |
| <b>Filters</b>                       | HEPA, 99.99%, 0.2 µm<br>high efficiency filter,<br>2,000 hours | Non-HEPA, 99.7%, 0.3<br>µm high efficiency filter,<br>500 hours                  | HEPA, 99.99%, 0.2 µm<br>high efficiency filter,<br>2,000 hours |
| <b>Storage Cart</b>                  | Height Adjustable<br>Storage Cart                              | Fixed Height<br>Storage Cart   | Storage Cart connected<br>to pole clamp                        |
| <b>Sterilized Blankets</b>           | Yes  | Only Cardiac   | Only Cardiac and<br>Surgical Access                            |
| <b>Blanket Inlets</b>                | Up to 4 inlets   | 1 hose inlet   | Up to 2 inlets   |
| <b>Machines Needed to Heat Gowns</b> | 1  | 2  | 1  |