

Designed for adult, pediatric, and neonatal patients.

The refurbished Siemens Servo 300A Ventilator offers a wide range of ventilation modes with an Automode feature for greater treatment possibilities, shorter weaning times and increased patient comfort.

Automode Feature: After two consecutive patient triggers, the ventilator shifts from controlled to supported ventilation. It remains in the support mode as long as the patient continues triggering. If the patient stops breathing, the ventilator shifts back to the control mode after twelve seconds for adult patients, eight seconds for pediatric patients, and five seconds for neonatal patients.



Technical Specifications

Ventilation mode	s Controlled • Supported • Combined
Patient range	Adult/Pediatric/Neonate
Classification	IEC60601-1/EN60 601-1 Type B
Method of trigger	ring Flow and Pressure
Dimensions	
Patient unit	9.5" W x 14.6" D x 9.4" H
	(242 mm W x 370 mm D x 240 mm H)
Interface	17" W x 5.9" D x 12.8" H
	(431 mm W x 150 mm D x 325 mm H)
Weight	Approximately 52.9 lbs. (24 kg)

A Quick Look

- » Automode
- » Neonate to adult patient range
- » Wide range of ventilation modes
- » Easy to use
- » O₂ air gas module filters bacteria

