


DETECTING AND REGULATING PFSA CONTAMINATION IN WATER: ADVANCEMENTS IN LC-MS/MS METHODOLOGY

Florentina Laura Chiriac, Florinela Pirvu, Iuliana Paun, Antonia Ioana Cimpean, Vasile Ion Iancu


National Research and Development Institute for Industrial Ecology-ECOIND, Drumul Podu, Dambovitei Street 57-73, 070762 Bucharest, Romania

Background

- PFAS Overview:** Synthetic chemicals linked to environmental and health issues.
- Regulatory Focus:** EPA's timeline for drinking water regulations targeting PFOS and PFOA.
- Analytical Demand:** Increased need for effective detection methods for PFAS.
- Detection Techniques:** Essential use of LC-MS/MS for accurate PFAS measurement.
- Public Health Impact:** Need for advanced methods to meet standards and reduce contamination risks.
- Aim of the Work:** Development and validation of a new LC-MS/MS method for PFAS detection in drinking water.



Methodology




LC conditions

- ❖ Chromatographic column: Zorbax eclipse XDB-C18
- ❖ Colum temperature: 30°C
- ❖ Injection volume: 20 µl
- ❖ Mobile phase: 5 mM ammonium acetate/ MeOH
- ❖ Flow: 0.4 mL/min
- ❖ Sample solvent: water
- ❖ Elution mode: gradient program
- ❖ Chromatographic run-time: 20 min


Time (min)	A (%)	B (%)
0.00	85	15
1.00	75	25
1.50	45	55
5.50	30	70
9.00	0	100
10.0	0	100
10.1	85	15

Chromatogram



Calibration standard with PFAS and labelled IS standards, at 500 ng/L (in sample concentration after 500x sample pre-concentration during sample preparation)

Method performance



Performance parameters for PFAS analysis in drinking water using LC-MS/MS

No.	Analyte	LOD (µg/L)	LOQ (µg/L)	Linearity	Concentration coefficient	Matrix effects (%)	Intra-day precision (%)	Inter-day precision (%)
1	PFBA	0.01	0.03	0.05 - 0.50	500	102 ± 6.3	3.2	9.4
2	PFMPA					100 ± 4.8	3.7	9.2
3	4:2 FTSA					98 ± 7.4	4.1	10.6
4	PFEESA					103 ± 4.2	4.8	10.2
5	PFHxA					108 ± 6.3	6.9	12.3
6	PFBS					102 ± 6.3	5.7	11.7
7	NFDHA					96 ± 6.4	6.3	10.4
8	HFPO-DA					94 ± 7.8	8.1	13.8
9	PFMBA					99 ± 7.0	4.9	10.2
10	PFPeA					103 ± 7.1	5.7	10.9
11	PFHxS					96 ± 8.0	6.4	11.7
12	DONA					98 ± 6.6	8.3	12.4
13	PFHpA					104 ± 8.7	5.6	11.7
14	PFHpS					107 ± 8.5	7.1	12.8
15	6:2 FTSA					93 ± 6.0	5.8	11.4
16	PFOA					90 ± 8.3	6.9	12.3
17	PFNS					105 ± 6.1	5.8	12.2
18	9Cl-PF ₃ ONS					96 ± 5.4	7.4	13.5
19	8:2 FTSA					94 ± 4.5	7.6	13.8
20	PFDA					97 ± 8.9	5.3	10.7
21	PFOS					105 ± 5.0	6.4	11.2
22	PFNA					102 ± 6.6	6.8	12.4
23	11Cl-PF ₃ OUdS					106 ± 7.6	7.2	13.5
24	PFDODA					94 ± 7.8	6.9	12.6
25	PFDS					97 ± 9.1	8.2	13.9
26	PFUnDA					92 ± 8.4	6.5	11.8
27	PFOSA					103 ± 6.4	7.3	12.7

Conclusions

- **Effective Detection:** LC-MS/MS achieved LODs as low as 0.01 µg/L for 27 PFAS, meeting EPA and Ord. 7/2023 requirements.
- **Reliable Performance:** The method demonstrated strong linearity and excellent precision, ensuring accurate quantification for routine monitoring.
- **Minimal Interference:** The analysis exhibited low matrix effects (90-108%), confirming robustness in complex water samples.
- **Comprehensive Analysis:** Detection of 27 PFAS highlights the method's applicability for regulatory compliance and risk assessment.
- **Public Health Impact:** This analytical approach supports effective monitoring, enhancing protections against PFAS contamination in drinking water.

Acknowledgments

This work was carried out through the “Nucleu” Program within the National Research Development and Innovation Plan 2022–2027 with the support of Romanian Ministry of Research, Innovation, and Digitalization, contract no. 3N/2022, Project code PN 23 22 01 01.

Disemination