

Welcome to the Supportline Catalogue





TÜRKİYE • AFRICA • AMERICA • ASIA • EUROPE

WE ARE ALL AROUND THE WORLD

As Ersamed, our pursuit of excellence has been continuing since the first day we started production in 2005;

We test ourselves every day in research and development processes with our innovative working model. With our expert staff and state-of-the-art machinery, exceeding our limits every day and being better than the previous day has become a life model for us.

We would like to thank our partners and customers who stand with us during this process, supported us in our progress and believed in us and our products.

Thanks to your feedback and constructive criticism, you have been our biggest supporters on this path we set out to be the best.



Index

Neck Collars	04-09
Shoulder & Arm	10-19
Hand & Wrist Orthosis	20-29
Corsets	30-41
Knee & Hip	42-51
Ankle Orthosis	52-59
Universal Group	60-67
Foot Health	69-77
Knitted Group	78-85
Elastic Group	86-89
Kids	90-99
Compression Stockings	100-105
Pillows	106-109

Neck Collars







SL-05D Neck Collar (Foam Backed)



S	M	L
28 cm 33 cm	34 cm 39 cm	40 cm 45 cm

Features

This is a sponge collar which is designed according to the anatomical shape of the neck.

Remind the patient to avoid excessive neck movements Provides feedback and keeps the neck area warm

Indications

- Cervical Hernia
- Muscle to be strained
- Osteoarthritis
- Muscel Spasms



SL-05 Neck Collar (Sponge)



S	M	L
28 cm 33 cm	34 cm 39 cm	40 cm 45 cm

Features

This is a sponge collar which is designed according to the anatomical shape of the neck.

Remind the patient to avoid excessive neck movements. Provides feedback and keeps the neck area warm

Indications

- Cervical Hernia
- Muscle to be strained
- Osteoarthritis
- Muscel Spasms



ERS-110 Nelson Collar



XS	S	M	L	XL
28 cm 31 cm	32 cm 35 cm	36 cm 39 cm	40 cm 43 cm	43 cm 46 cm

Features

It's made of 12 mm perforated plastozot.

Indications

- Cervical Hernia
- Straines
- Used for surgery, post-traumatic rest and protection.

ERS 110K Nelson Collar With Fabric



Features

It's made of 12 mm perforated plastozot. It is covered with a cotton fabric that does not sweat.

Indications

- Cervical Hernia
- Straines
- Used for surgery, post-traumatic rest and protection.

XS	S	M	L	XL
28 cm 31 cm	32 cm 35 cm	36 cm 39 cm	40 cm 43 cm	43 cm 46 cm

ERS-111 Philedelphia Collar



Features


Made of 12 mm plastosote material and supports semi-rigid plastic on the front and back.

Indications

- Cervical Hernia
- Osteorthritis
- Used for protetion after surgery

S	M	L
33 cm 37 cm	38 cm 42 cm	42 cm 46 cm

ERS-112 Philadelphia Tracheotomy Collar



Features

Made of 12 mm plastosote material and supports semi-rigid plastic on the front and back.


Indications

- Cervical Hernia
- Osteorthritis
- Used for protetion after surgery

S	M	L
33 cm 37 cm	38 cm 42 cm	42 cm 46 cm



ERS-108 Cervical Collar W/O Chin Support



S	M	L	XL	XXL
32 cm 35 cm	36 cm 39 cm	40 cm 43 cm	43 cm 46 cm	46 cm 50 cm

Features

Designed according to the anatomical shape of the neck. It's made of plastic collar that do not support the jaw, Custom height and length adjustments can be made with Velcro.

Indications

- Cervical Hernia
- Osteoarthritis
- As a protection in case of injuries.



ERS-109 Cervical Collar With Chin Support



S	M	L	XL	XXL
32 cm 35 cm	36 cm 39 cm	40 cm 43 cm	43 cm 46 cm	46 cm 50 cm

Features

Designed according to the anatomical shape of the neck. It's made of plastic collar which support the jaw, Custom height and length adjustments can be made with Velcro.

Indications

- Cervical Hernia
- Strains
- Used for surgery, post-traumatic rest and protection.



ERS-114 Wire Frame Collar



Features


Designed according to the anatomical shape of the neck. It's a collar that can be adjusted in bending and stretching mode and its outer surface is semi rigid.

Indications

- Cervical Hernia
- Osteoarthritis
- Used for protection after surgery

SHORT-REGULAR

ERS-116 Minerva Orthosis



Features

It is a rigid orthosis used for stabilization of cervical and upper thoracic area. It's made of polyethylene material at outer side and terry cloth at inner side which providing maximum comfort in the area, and flexion, extension lateral flexion and also rotation movements of cervical area will be limited.

Indications

- C3-T2 stable fractures
- Post-trauma, muscular /ligamentous injury
- Degenerative joint disease

Chest Size

S	M	L
70 cm	90 cm	110 cm
90 cm	110 cm	130 cm

Neck Size

S	M	L
33 cm	38 cm	42 cm
37 cm	42 cm	46 cm

Shoulder & Arm







SL-01 Arm Sling (Lux)



Features

Inner and outer surface made of sponge material which is provide extra protection supports.

Indications

- Used for full or half dislocations in shoulder
- Used for post operative prorection and reheal
- Used to prevent frctures formation of the distall edema

S	M	L	XL	XXL
30 cm 33 cm	33 cm 36 cm	36 cm 39 cm	39 cm 42 cm	42 cm 45 cm



SL-01 F Arm Sling (Mesh)



Features

Inner and outer surface made of perforated fabric which is air permeability and antiperspirant material.

Indications

- Used for full or half dislocations in shoulder
- Used for post operative prorection and reheal
- Used to prevent frctures formation of the distall edema

S	M	L	XL	XXL
30 cm 33 cm	33 cm 36 cm	36 cm 39 cm	39 cm 42 cm	42 cm 45 cm



SL-01 B Arm Sling (Simple)



Features

The inner surface of the towel and upper part surface is made of spongefabric.

Indications

- Used for full or half dislocations in shoulder

UNIVERSAL



SL-13 Acromio-Clavicular Bandage

Features

It is made of breathable, non-sweaty fabric. The shoulder and armpit pads made of sponge provide comfortable use.

Indications

Acromio-clavicular dislocations and semi-dislocations

S	M	L
30 cm 33 cm	33 cm 36 cm	36 cm 39 cm



SL-14 Arm Shoulder Stabilizer



UNIVERSAL

Features

Arm sling made of sponge fabric, used to immobilize the shoulder and upper arm. With its functional design, the front velcro can be opened when necessary, allowing easy treatment and control.

Indications

- Used for full or half dislocations in shoulder.
- Used for post operative protection and reheel.
- To provide non-surgical stabilization.
- To prevent edema formation in case of distal fracture.
- During conservative and post-operative therapy of scapula fractures.



SL-02 Velpeau Bandage



S	M	L	XL	XXL
30 cm	33 cm	36 cm	39 cm	42 cm
33 cm	36 cm	39 cm	42 cm	45 cm

Features

It's made of sponge cloth abduction for shoulder and upper arm and allows immobilization, also the bandage around the chest increase the stability

Indications

- Used for dislocations in shoulder
- Used after surgery for protection and reheel
- Are used for Upper of the scapula and humerus fractures rotator cuff injury



SL-03 Arm Sling (30°)



UNIVERSAL

Features

The inner surface of the towel the outer surface of the sponge cloth produced. Provides 30° abduction immobilization shoulder. Prevent post-operative adhesions. Mobilization and rehabilitation help.

Indications

- The impingement syndrome of post-op
- Rotator cuff tear after repair,
- Frozen shoulder under anesthesia after manipulation.
- Subacromial bursitis and subcoracoid.
- Used after shoulder endoprosthesis.



SL-03 D Arm Sling (45°- 60°)



Features

The inner surface of the towel the outer surface of the sponge cloth produced. Provides 45-60° abduction immobilization shoulder. Prevent post-operative adhesions. Mobilization and rehabilitation help.

Indications

- The impingement syndrome of post-op
- Rotator cuff tear after repair,
- Frozen shoulder under anesthesia after manipulation.
- Subacromial bursitis and subcoracoid.
- Used after shoulder endoprosthesis.

UNIVERSAL



SL-04 Clavicle Bandage



Features

The inner part of the product is made of towel kind of material which is adjustable and provides the shoulder pulled back to the desired tension and ensures proper posture. Its soft and antiperspirant the content of the material provides comfort while use.

Indications

- Used for clavicle fracture
- Clavicle and sternoclavicular Joint Disorders

XS	S	M	L	XL	XXL
18 cm 24 cm	25 cm 30 cm	31 cm 36 cm	37 cm 42 cm	43 cm 48 cm	49 cm 54 cm



REF-620 Shoulder Support



Thigh Support

Made of neoprene material, and with compression effect prevents edema. Reduces pain by reducing muscle tone. Lower-end of the strap at the velcro for easy use, there are under arm pad.

Indications

- Semi-dislocations injuries
- Tendinitis osteoarthritis symptom
- Used after prosthesis shoulder
- Used after surgery for protection

S	M	L	XL	XXL
22 cm 24 cm	24 cm 27 cm	27 cm 30 cm	30 cm 34 cm	34 cm 38 cm



REF-621 Humerus Brace Splint



S	M	L	XL	XXL
30 cm 33 cm	33 cm 36 cm	36 cm 39 cm	39 cm 42 cm	42 cm 45 cm

Features

It is a shoulder orthosis made of hard sponge material and used in all situations that require immobilization of the shoulder joint in internal rotation.

Indications

- Proximal humerus fractures
- Following proximal humeral shaft fractures
- Non-surgical, Post operative



REF-300 Neoprene Elbow Support



S	M	L	XL	XXL
30 cm 33 cm	33 cm 36 cm	36 cm 39 cm	39 cm 42 cm	42 cm 45 cm

Features

Its anatomically designed from Neoprene material.

Indications

- Protect against Trauma
- Using for compression and heat on injured or pain area.
- Protect the area during and after sport activities



REF-301 Neoprene Elbow Support - With Epicondylitis Band



S	M	L	XL	XXL
30 cm 33 cm	33 cm 36 cm	36 cm 39 cm	39 cm 42 cm	42 cm 45 cm

Features

It's made of neoprene material. This anatomical fabric is designed to create compression on the muscle groups attached to the lateral epicondyle and reduce tension and pain.

Indications

- Epicondylitis Lateral issues
- Problem with increamentation of two joint fluid.
- Using for compression on injured or pain area.

REF-302 Neoprene Elbow Support - With Pad



S	M	L	XL	XXL
22 cm 24 cm	24 cm 26 cm	26 cm 28 cm	28 cm 30 cm	30 cm 32 cm


Features

Made from Neoprene material and anatomically designed. Elbow area has a special design by eva material and provides better support and protection.

Indications

- Protect against Trauma
- Olecranon Bursitis
- Protect the area during and after sport activities

ERSA-217 Cubital Tunnel Elbow Support



UNIVERSAL

Features

It is made of neoprene material. There is a curved steel support on the inside of the elbow.

By preventing the pole joint from bending, it prevents the ulnar nerve from being stretched and relieves the symptoms.

Indications

- Cubital tunnel syndrome
- Ulnar nerve problems

REF-303 Epicondylitis Band



UNIVERSAL

Thigh Support

Made of neoprene material ,and with compresion effect prevents edema.Reduces pain by reducing muscle tone. Lower-end of the strap at the velcro for easy use, there are under arm pad.

Indications

- Semi-dislocations injuries
- Tendinitis ostorartrite symptom
- Used after prosthesis shoulder
- Used after surgery for protection



REF-304 Epicondylitis Band



Features

It is made of neoprene material and anatomically designed. Has a special pad with three-point silicone, which provides compression and reduce the tension in tendon. It is used for medial or lateral epicondylitis. Easy to use.

Indications

- Semi-dislocations injuries
- Tendinitis ostarartrite symptom
- Used after prosthesis shoulder
- Used after surgery for protection

UNIVERSAL



SL-906 Arm Orthosis (Adjustable Angle Support)



Features

This is an adjustable elbow angle orthosis manufactured anatomically in accordance with elbow shape. Its made of sponge fabric which inside is made of Terry cloths. Angle adjustment can be fulfilled easily by a pressing button. The orthosis is controlling elbow flexion and extension. Angle flexion adjustment is possible in degrees. According to the arm length, adjustment on the orthosis can be performed.

Indications

- Helping rehabilitation in elbow contracture
- After ligament Injury
- When require to controlled the limited motion of the elbow
- Conservative treatment after elbow dislocation

S	M	L
22 cm 24 cm	24 cm 26 cm	26 cm 28 cm

Shoulder & Arm



Hand & Wrist Orthosis



supportline™







SL-15 Dequervain Splint



Features

It's made of air permeability cotton fabric.
It's restricting the movement with adjustable aluminium Support.

Indications

- De quervain syndrome treatment
- Joint irritations purpose
- Used to protect thumb ligaments laxity

S	M	L	XL	XXL
15 cm 17 cm	17 cm 19 cm	19 cm 21 cm	21 cm 23 cm	23 cm 25 cm



REF-19 Hand & Wrist Splint - Both Sided



Features

Manufactured with respirable perforated fabric.
On the volar side of the wrist there is an aluminium support which is anatomically designed and provides stabilization for the lace system around the wrist.

Indications

- The aluminium rod inside constrains the movements of the wrist.
- Provides better stabilization
- It has 19 cm length

S	M	L	XL
15 cm 17 cm	17 cm 19 cm	19 cm 21 cm	21 cm 23 cm



SL-19 Hand & Wrist Splint - Single Sided



Features

It's made of cotton material and air permeability,
the aluminium part support the wrist and providing neutral position and avoid extra movement.
It has a bandage around the wrist which provide protection and stabilization.

Indications

- Carpal tunnel syndrome
- Tendonitis tenosynovitis, tenosynovitis
- Reduced pain in hand wrist
- Used for protection and support after injury or surgery

S	M	L	XL	XXL
15 cm 17 cm	17 cm 19 cm	19 cm 21 cm	21 cm 23 cm	23 cm 25 cm



SL-21 Hand & Wrist Splint - Single Sided (With Thumb Support)



Features

It's made of air permeability cotton fabric. It's restricting the movement of wrist and thumb with adjustable aluminium support

Indications

- In case of carpal tunnel and de quervain
- De quervain tenosynovitis
- Join degeneration
- Reheel the pain in hand wrist and thumb

S	M	L	XL	XXL
15 cm 17 cm	17 cm 19 cm	19 cm 21 cm	21 cm 23 cm	23 cm 25 cm



REF-601 Hand & Wrist Splint - Single Sided



Features

this brace wrist made of neoprene material,with the anatomic shaped prevent extra movement
Increase support and stabilization around the wrist .

Indications

- Carpal tunnel syndrome .
- Tendonitis, tenosynovitis, tenovajinit issues.
- Used to reduce the hand wrist pain.
- Used for protection and reheal after surgery and injury.

S	M	L	XL	XXL
15 cm 17 cm	17 cm 19 cm	19 cm 21 cm	21 cm 23 cm	23 cm 25 cm



REF-602 Hand & Wrist Splint - Single Sided (With Thumb Support)



Features

Thumb fixation splint made of neoprene material, used in all indications that require thumb fixation and immobilization of the wrist. Available for right and left hand.

Indications

- Distortion of the wrist and/or thumb.
- Postoperative immobilization of the wrist and thumb.
- Lesions of collateral ligaments of the basal thumb joint.
- Tendinitis, tenosynovitis of the wrist and/or thumb.
- Inflammatory phases of rheumatic diseases.
- Rhizarthritis.

S	M	L	XL	XXL
15 cm 17 cm	17 cm 19 cm	19 cm 21 cm	21 cm 23 cm	23 cm 25 cm



REF-603

Lace Hand Wrist Splint (Both Hand)

S	M	L	XL
15 cm 17 cm	17 cm 19 cm	19 cm 21 cm	21 cm 23 cm

Feature

Manufactured with, air permeability perforated fabric. There is an aluminium support which is anatomically formed and a bandage which provides stabilization for the lace system around the wrist. Easy to use and apply

Indications

- Designed with special elastic cotton fabric.
- The aluminium rod which is placed inside restrict the wrist movements.
- Provides better stabilization.



REF-604 Hand & Wrist Splint - Single Sided (With Thumb Support)



S	M	L	XL
15 cm 17 cm	17 cm 19 cm	19 cm 21 cm	21 cm 23 cm

Features


Manufactured from respirable perforated fabric. Aluminum support is anatomically formed to fit the volar side of the wrist. The bandage provides stabilization for the lace system around the wrist. Easy to use and apply.

Indications

Carpal Tunnel Syndrome.
Tendinitis, tenosynovitis, tenovaginitis.
Immobilization for wrist pain recovery
Post-op and protection purposes.
In fractures after encasement in plaster.



ERSA-210 Wrist - Forearm Splint



Features

It is a type of long splint that provides immobilization of the wrist in neutral position. It is designed to start from the forearm and grasp the wrist. There is an aluminum support under the arm inside. It is possible to use both hands by changing the location of the aluminum support.

Indications

- Carpal Tunnel Syndrome.
- Tendon and soft tissue damage.
- Rehabilitation after surgery

UNIVERSAL



ERSA-220 Boxer's Fracture Splint



Features

It provides stabilization in the wrist and in the little and ring fingers. It can be adjusted to the desired tightness with the velcro straps on it. There are two aluminum splints inside for maximum stabilization and immobilization.

Indications

- Metacarpal joint fractures.
- Rehabilitation after surgery

UNIVERSAL



ERSA-225 Middle Finger Fracture Splint



Features

Provides stabilization in the wrist and middle finger. It can be adjusted to the desired tightness with the velcro straps on it. There is an aluminum splint inside for maximum stabilization and immobilization.


Indications

- Finger fractures.
- Tendon and nerve injuries
- Finger injuries and sprains
- Rehabilitation after surgery

UNIVERSAL



SL-901 Thermoplastic Cleinert Splint



Features

Manufactured from thermoplastic material with an inner terry cloth coating. Allows passive flexion and active extensions. In flexor tendon injuries provides early controlled motion after tendon repair and prevents the formation of adhesions.

Indications

- Flexor tendon injuries

M	L
14 cm 18 cm	18 cm 22 cm



SL-902 Thermoplastic Hand Splint



Thigh Support

Hand and wrist orthosis manufactured from thermoplastic material with an inner terry cloth coating. Keeps the fingers and the wrist in recovery position. Enables immobilization of wrist, used for resting position following plaster treatment of fractures and painful positions of the wrist.

Indications

- In fractures after encasement in plaster
- Sprain and injury treatment
- Hemiplegia
- Post-op recovery
- Spasticity deformity protection

M	L
14 cm 18 cm	18 cm 22 cm

SL-904 Thermoplastic Dynamic Splint



M	L
14 cm 18 cm	18 cm 22 cm

Features

Hand and wrist orthosis manufactured from thermoplastic material with an inner terry cloth coating. Keeps the fingers and the wrist in passive extension. Active finger flexion is allowed. In extensor tendon injuries provides early controlled motion after tendon repair and prevents the formation of adhesions.

Indications

- Extension tendon injuries

SL-905 Thermoplastic Anti-Spasticity Splint



M	L
14 cm 18 cm	18 cm 22 cm


Features

Hand and wrist orthosis manufactured from thermoplastic material with an inner plastizode coating. Due to the grooved finger design and materials resistance to tension, splint generally used for hand spasticity.

Indications

- Spasticity
- Wrist problems

SL-602 The Baseball Splint



S	M	L
4,5 cm 5 cm	5 cm 6 cm	6 cm 7 cm

Features


Manufactured from special aluminum with inner plastizode coating. Provides finger immobilization in extension.

Indications

- Finger fractures
- Tendon and nerve injuries
- Finger injuries, sprains and trauma



SL-603 Stabilization Splint



S	M	L
4,5 cm 5 cm	5 cm 6 cm	6 cm 7 cm

Features


Manufactured from special aluminum with inner plastizode coating. Provides finger immobilization in extension.

Indications

- Finger fractures
- Tendon and nerve injuries
- Finger injuries, sprains and trauma



SL-604 Flat Finger Splint



S	M	L
4,5 cm 5 cm	5 cm 6 cm	6 cm 7 cm

Features


Flat Finger Splint.

Indications

- Finger fractures
- Tendon and nerve injuries
- Finger injuries, sprains and trauma



SL-601 Mallet Finger



1	2	3	4	5	6	7
4 cm 4,5 cm	4,5 cm 5 cm	5 cm 5,5 cm	5,5 cm 6 cm	6 cm 6,5 cm	6,5 cm 7 cm	7 cm 8 cm

Thigh Support

Finger splint manufactured from special plastic material, keeping the DIP joint in extension. The fingers are mobile from PIP joint on.

Indications

- Mallet finger deformity
- PIP joint fractures, traumas and sprains

Hand & Wrist Orthosis



Corsets







SL-11 Woolen Corset



Features

It is made of wool fabric and it keeps the waist warm with compression.

Indications

- Abdomen surgery and post-natal period.
- Using to cover the excess abdomen area

S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm



SL-210 Abdominal Corset



Features

It is used for mild back pain, but also to support and shape the abdomen

Indications

- Abdomen surgery and post-natal period.
- Using to cover the excess abdomen area

S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm



SL-215 Abdominal Corset - With Hook



Features

It is used for mild back pain, but also to support and shape the abdomen.

Indications

- Abdomen surgery and post-natal period.
- Using to cover the excess abdomen area.

S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm



SL-220 Work Corset



Features

It is made of neoprene fabric that supports the waist and back area at its height. There are 2 plastic supports in the back area. Two cross-banded belts provide stabilization from the back area to the waist area. blood circulation by maintaining body temperature; thus speeding up the healing process


Indications

- Lumbar spine hernias and calcifications
- In sacralization and lumbalization
- Sacroiliac joint problems
- Sudden and long-term low back pain
- Rehabilitation after surgery

UNIVERSAL



SL-240 Abdominal Corset - Elastic



Features

It made from elstic material , air permeability and antiperspirant fabric also provide support for abdomen.

Indications

- Abdomen surgery and post-natal period.
- Inperitoneum injuries.
- Using to cover the excess abdomen area.

S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm



SL-230 Uni-Thorax Corset



Features

It is a corset that supports the sternum area. It has two plastic supports on the back. It has a shape pads located on two aluminium underwires in the sternum area. Thanks to the bands on the side, it fits the body perfectly. With two plastic buckles in the sternum area, it allows the patient to perform health care without removing the corset during the use of the patient.

Indications

- Rib fractures
- In rehabilitation after surgery of the rib cage

UNIVERSAL



SL-235 Pectus Carinatum



UNIVERSAL

Features

The front support rest to sternum and the back support rest to spine. It consists of soft bands that support the rib cage from the sides. It can be adjustable with velcro. The presence of adjustable shoulder straps stabilizes the orthosis better on the body. The pad that supports sternum can be adjusted according to the location of the deformity.

Indications

- Pectus carinatum deformity



SL-241 Chest Corset



Features

It made from elastic material, air permeability and antiperspirant fabric, also the sternal area is supported with towel fabric in this chest corset. With supplying compression chest cage, support sternum

Indications

- For supporting the sternum after coroner surgery
- For supporting costa and sternum fracture part.
- It is used after chest surgery and med trauma

S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm



SL-242 Colostomy Corset



Features

It made from elastic cotton fabric and antiperspirant material. It used for protection of colostomy bag.

S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm

SL-243 Umbilical Hernia Corset




Features

It made from elastic cotton material , air permeability and antiperspirant fabric with compression in herinated area. Also has the special designed pad for protecting the umbilical hernia.

S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm

SL-244 Maternity Corset



Features

It made from elastic cotton material , air permeability and unsweatability fabric,it has abdominal support brace for pre-natal and post natal purpose. This post-partum corset is helping to recover the wide range of sags for front part of the abdomen. And gives people a feeling of confidence and comfort for movement in everyday life.

S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm

SL-260K Lumbosacral Corset (26 cm)



Thigh Support

It made from elastic material , air permeability and antiperspirant fabric,has steel and elastical underwire with 26 cm height. Also helping to reduce the load of the discs

Indications

- It used for lumbar degenerative disc disease
- Spodilolizis spondilolistezis
- Redusing the acute or cronical pains
- It is used to avoid extra movement after surgical operation

S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm



SL-261 Lumbosacral Corset With Abdominal Support



S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm

Features

Manufactured from flexible, respirable, sweat-preventing fabric. At the back, it has 4 pieces of anatomically designed steel underwire. Steel underwires can be manually adjusted to fit lumbar lordosis. The front side is composed of two pieces. The lower and wider piece supports to the lower abdomen. The narrow piece on the front side is fastened in a transverse manner. Relieves the pain on the disks at the lumbar region, by taking the pressure off.

Indications

- Lumbago, degenerative disc lesions
- Spondylolysis, spondylolisthesis
- Acute or chronic lumbar pain
- Pre-and Post-Surgical Stabilization
- Degenerative Spinal pathologies on Lumbar area



SL-320K Lumbosacral Corset (32 cm)



S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm

Features

It made from elastic material, air permeability and antiperspirant fabric, also has 4 pieces of steel underwire with 32 cm height which is helping to reduce the load of the discs to extra crosswise

Indications

- It used for lumbar degenerative disc disease
- Spodilolizis spondilolistezis
- Redusing the acute or cronical pains
- It is used to avoid extra movement after surgical operation



SL-270 Lumbosacral Corset With Belt (26 cm)



S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm

Features

It made from elastic material, air permeability and antiperspirant fabric, also has 4 pieces of steel underwire with a fabric which is support the lumbr at the back.

Indications

- It used for lumbar degenerative disc disease
- Spodilolizis spondilolistezis
- Redusing the acute or cronical pains
- It is used to avoid extra movement after surgical operation



SL-271 Lumbosacral Corset - Airtex




Features

It made from elastic material, air permeability and antiperspirant fabric, also has 4 pieces of steel underwire which is support the lumbr at the back.

Indications

- It used for lumbar degenerative disc disease
- Spodilolizis spondilolistezis
- Redusing the acute or cronical pains
- It is used to avoid extra movement after surgical operation

S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm



SL-161 Posturex Corset



Features

It made from elastic material, air permeability and unsweatability fabric. It correct the posture and the sinal correction is evident after use

Indications

- Posture disorders
- It is used for weak back and shoulder muscles.

S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm



SL-162 Posturex Corset



Features

It made from elastic material, air permeability and unsweatability fabric. It correct the posture and the sinal correction is evident after use

Indications

- Posture disorders
- It is used for weak back and shoulder muscles.

UNIVERSAL



SL-330K Dorsolumbar Corset



S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm

Features

It made from elastic material , air permeability and antiperspirant fabric,and its covering the whole back with metal sprint on both sides and with an extra belt support.

Indications

- It used for Vertebral trauma
- It used for Osteoporsis disease
- Light kyphos is used in



SL-911 TLSO - Corset



S	M	L
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm

Features

Body orthosis manufactured from sponge laminated terry cloth with inner side. 3 mm Plastic inner under wires provide a rigid support to the abdominal and lumbar parts of the corset as an addition to the 2 mm plastic bearers on the sternum. Plastic bearers on the sternum in particular has a significant difference restricting the flexion movement of the body when compared to dorsolumbar corset. Restricts the motion of the vertebrae to a certain point, secures accurate anatomical posture and enhances stability thus accelerating recuperation process.

Indications

- Vertebral fractures
- Post-op recovery and preservation purposes
- Compression fractures
- Post-operative fusion
- Kyphosis
- Postural deformities



SL-912 Hyperextension Corset (Jweet)




S	L
75 cm 85 cm	95 cm 105 cm

Features

Body orthosis manufactured from sponge laminated terry cloth with inner side. 3 mm Plastic inner under wires provide a rigid support to the abdominal and lumbar parts of the corset as an addition to the 2 mm plastic bearers on the sternum. Plastic bearers on the sternum in particular has a significant difference restricting the flexion movement of the body when compared to dorsolumbar corset. Restricts the motion of the vertebrae to a certain point, secures accurate anatomical posture and enhances stability thus accelerating recuperation process.

Indications

- Vertebral fractures and trauma
- Compression fractures in osteoporosis Kyphosis
- Post-op protection and rehabilitation purposes
- Stable fractures in the lower thoracic and lumbar spine
- Compression fractures
- Osteoporosis



SL-350 Spinal Support

S	M	L
30 cm 33 cm	33 cm 36 cm	36 cm 39 cm

Features


It is a corset consisting of 3 millimeter aluminium with eva support inside. It is made of a sponge cloth with inner surface. It is orthosis can be used in vertebral body fractures. It can be also used in the flattening of the back. It stabilizes the patient with the shoulder straps and waist belt.

Indications

Osteoporotic bone collapse in the thoracic and lumbar spine
Juvenil scheuermann disease
Chronic back pain hyperkiphosis



SL-915 Hyperextension Corset



Features

Body orthosis manufactured from sponge laminated terry cloth with inner side. 3 mm Plastic inner under wires provide a rigid support to the abdominal and lumbar parts of the corset as an addition to the 2 mm plastic bearers on the sternum. Plastic bearers on the sternum in particular has a significant difference restricting the flexion movement of the body when compared to dorsolumbar corset. Restricts the motion of the vertebrae to a certain point, secures accurate anatomical posture and enhances stability thus accelerating recuperation process.

Indications

- Vertebral fractures and trauma
- Compression fractures in osteoporosis Kyphosis
- Post-op protection and rehabilitation purposes
- Stable fractures in the lower thoracic and lumbar spine
- Compression fractures
- Osteoporosis

S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm



SL-920 Taylor Corset



Features

Adapting the treatment according to the needs of each stage can improve stabilisation or provide greater freedom of movement as the patient progresses through the rehabilitation process.

Indications

- Padded shoulder straps that can be adjusted for a comfortable fit.
- Central containment reinforcement with 2 adjustable steel stays.
- Micro velcro for secure fastening in any position.

S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm




SL-513 Inguinal Hernia Belt - Single Sided



Thigh Support

Made of elastic cotton fabric. Used for post operative herina to protect growth and prevent the outbreak of the sewing. pre-and post sagittal used in herina and wash in warm water with soap but without rubbing

RIGHT/LEFT



SL-513B Inguinal Hernia Belt - Both Sided




UNIVERSAL

Features

Made of elastic cotton fabric. Used for post operative herina to protect growth and prevent the outbreak of the sewing. pre-and post sagittal used in herina and wash in warm water with soap but without rubbing.



SL-924 Supporter Of Suspensory



Features

Non-irritating to the part made of special waist bag to hold the cotton fabric is added to the soft elastic band

5	6	7	8	9	10
50 cm 60 cm	60 cm 70 cm	70 cm 80 cm	80 cm 90 cm	90 cm 100 cm	100 cm 110 cm

Knee & Hip




supportline™





SL-10 Woolen Knee Support



S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm

Features

It's made of woolen fabric. Causes heat and compression. Used for rheumatism and knee pain.

Warning

- Can be use by doctor recomendation
- Can be use as shown in photo image
- Only Hand wash with soap and water
- Can be dry in room remprature on a towel

**Color options available.*



REF-100 Neoprene Knee Brace



S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm

Features

It's made of neoprene material and the anatomical designed.

Warning

- Can be use by doctor recomendation
- Can be use as shown in photo image
- Only Hand wash with soap and water
- Can be dry in room remprature on a towel



REF-101 Neoprene Knee Brace - Pattela Support



S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm

Features

Its made of Neoprene material and anatomically designed. The patella area is open and supported by a pad.

Indications

- Osteoathritis
- Problem with increamentation of two joint fluid
- Decrease the pain in knee cap
- Protection and rehabilitation after surgery
- Protect during the sport activities
- It is used for post operative protection
- Used for compression and provide heat on the area



REF-102 Neoprene Knee Brace - Pattela - Ligament Support



S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm

Features

Its made of Neoprene material and anatomically designed. The patella area is open and supported by a pad. Flexible side wires increase stability and reduce the load on the knee.

Indications

- Arthritis, Bursitis and Tendinitis
- In strain and tensile
- Meniscus problems
- It is used for post operative protection
- Ligaments instability
- Used after sport injuries



REF-103 Neoprene Knee Brace - Pattela - Ligament - Band Support



S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm

Features

Its made of Neoprene material and anatomically designed. The patella area is open and supported by a pad. Flexible side wires increase stability and reduce the load on the knee.

Indications

- Arthritis, Bursitis and Tendinitis
- In strain and tensile
- Meniscus problems
- It is used for post operative protection
- Ligaments instability
- Used after sport injuries



REF-104 Neoprene Knee Brace - Steel Hinged



S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm

Features

It's made of neoprene material and anatomically designed. Steel joints on both sides provide maximum support and control the movements. Hyperextension directry blocks.

Indications

- Ligaments injuries and instability
- Meniscus Lesions
- First or second degree strain or tensile
- It is used for protection and avoid extra movement after surgery



REF-105 Neoprene Knee Brace - Closed Padded



S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm

Features


Manufactured from neoprene material and anatomically designed. Patella and surrounding material protect and designed with special eva pad.

Indications

- Used for protection during sport activities
- Patella is surrounded with soft tissue to protect against injury.
- Used to protect the knee bursitis



REF-106 Neoprene Knee Brace - Post-Arthroscopy



S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm

Features

Manufactured from neoprene material and anatomically designed. Patella and surrounding material protect and designed with special eva pad.

Indications

- Treatment of post arthroscopy edema and pain
- Any condition where hot or cold treatment is required



REF-107 Neoprene Knee Brace - Patella & Ligament Crossed Band Support



S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm

Features

It has neoprene structure and anatomical design. Supports patella and surrounding. There are flexible wires on both sides. Also three-point original strips are available and stabilizes the anterior cruciate ligaments and patella.

Indications

- Using in case of cruciate ligament and knee injuries and instability
- For patella stabilization
- For meniscus problems
- It is used for protection and stabilization of postoperative



REF-110 Patellar Tendinitis Band



Features

It's made of neoprene material and is anatomically designed. It has a special silicone pad that support the patellar tendon. Does not interfere with joint movements. The effect is decrease the pressure to reduce the tension and pain in the tendon.


Indications

- Osgood-schlatter disease,
- It is used in patellar tendon disorders.

S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm



REF-111 Patellar Tendon Brace



Features

Made of Neoprene material and anatomically designed .There is a special pad supporting the patellar tendon ,and has compression effect to reduce pain.

Indications

- Osgood-schlatter disease,
- It is used in patellar tendon disorders.

S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm



REF-500 Neoprene Thigh Support



Features

It is made of neoprene material and anatomically designed .support the thighs and provid heat and compression effect increases circulation and reduces pain.


Indications

- Muscle stress relifer
- In painful conditions of the upper calf
- To protect during and after sport activities
- It is used to provide heat and compression

S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm



REF-510 Neoprene Calf Support



S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm

Features


It is made of neoprene material and is anatomically designed. It supports the lower calf, the neoprene material protects the body temperature, increases blood circulation, accelerates the metabolism, the compression effect it applies prevents edema and reduces the feeling of pain in the muscles.

Indications

- Lower calf muscle strains
- Painful conditions of the lower calf
- Protection during and after sports activities
- For temperature and compression purposes



SL-12 Knee Immobilizer



45 - 55 - 65 CM

Features


The inside of sponge cloths towel is made of immobilization extention provides arrectory. It's anatomically designed from strong posterior, medial and lateral stabilization which providing support is made of steel with velcro. Available in standard bodies 45 cm, 55 cm, 65 cm long.

Indications

- Used after the operation for immobilization
- For rest after two injury-bonds
- To post operative meniscal injuries



SL-09 Angle Adjustable Knee Brace



ADJUSTABLE

Features

Adjustable-degree knee orthosis, manufactured from aluminum bar and sponge laminated cloth with a terry inner side. Flexion (0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°, 120°) and extension (-15°-0°, 15°, 30°, 45°, 60°, 75°, 90°) angle adjustment possible. Provides exceptional support and stability following knee injury. Use for protected activity during post-op rehab following ligamentous surgical repairs and reconstruction.

Indications

- Joint mobilization after ligament injury operations
- Stable femoral fractures
- High tibial osteotomies
- Gradual mobilization after operation



SL-09 N

Neoprene Angle Adjustable Knee Brace

Feature

Adjustable-degree knee orthosis, manufactured from aluminum bar and neoprene with a terry inner side. Flexion (0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°, 120°) and extension (-15°-0°, 15°, 30°, 45°, 60°, 75°, 90°) angle adjustment possible. Provides exceptional support and stability following knee injury. Use for protected activity during post-op rehab following ligamentous surgical repairs and reconstruction.

Indications

- Joint mobilization after ligament injury operations
- Stable femoral fractures
- High tibial osteotomies
- Gradual mobilization after operation



SL-09 A Knee Orthosis - 4 Point Principle Adjustment



S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm

Features

Manufactured from light aluminum, functional angle adjusted hinged rigid brace allowing limitations on extension and flexion, anatomically designed to fit the knee shape. Flexion (0°, 15°, 30°, 60°, 90°, 120°) ve extension (0°, 15°, 30°) angle adjustment is possible. Provides more convenient stabilization with the 4-point principle and prevents tibial displacement and knee hyperextension.

Indications

- Isolated rupture of ant. cruciate ligament or combined ant. and post. Cruciate ligament ruptures
- Instability caused by collateral ligament rupture
- Conservative treatment of combined instabilities



SL-900 Neoprene Hip Support



UNIVERSAL

Features

The circular compression on the hip and thigh ensures the required compression of the operation site and supports the direct postoperative positioning of the treated hip joint. At the same time, the elastic tension on the thigh counteracts the usual external rotation of the leg.

The compact design of the hip rotation support limits involuntary and uncontrolled movements, and offers the patient the additional sense of security that is necessary for rapid and wide-ranging mobilisation.

Indications

- High wearing comfort
- Skin-friendly material
- Simple application



SL-908 Hip Orthosis - Adjustable



UNIVERSAL

Features

It is a hip abduction orthosis that providing the stabilization of the hip joint and produced thermoplastic material designed anatomically. The inner surface that contacting to the skin is covered with terry cloth. Flexion adjustment can be made up to 0-125 degrees with 15 degrees steps. 45° abduction 30° abduction 45° hyperextension aspect adjustable.

Indications

- Stabilization of total hip prosthesis or revision
- Femur head resection
- Femoral or acetabular osteotomy
- Recurrent hip dislocation



Ankle Orthosis







REF-400 Eight Bandage - Elastic



UNIVERSAL

Features

It's made of soft cotton material with air permeability and it's easily used in shoes and provides flexibility and support to the ankle

Indications

- To support ligament weakness
- Used for ankle pain



REF-401 Neoprene Ankle Support - Cross Strap



Features

It has neoprene structure and anatomical design. The elastic yellow tape increase the stabilization and compression


Indications

- Used for 01 degree of twisting and stretching
- Used after minor fractures
- Used after injury compression for edema purpose
- Used for stabilization and protection

S	M	L	XL	XXL
21 cm 23 cm	23 cm 25 cm	25 cm 27 cm	27 cm 29 cm	29 cm 31 cm



REF-402 Neoprene Ankle Support - With Malleolar Pad



Features

The internal and external malleolus is supported by anatomically designed pad which is made of neoprene material. So will reduce the pressure applies on ankle and reduce the edema pain.

Indications

- Mild ligament injuries and instability
- Arthritis and osteoarthritis
- Used for support and protection after injuries or support

S	M	L	XL	XXL
21 cm 23 cm	23 cm 25 cm	25 cm 27 cm	27 cm 29 cm	29 cm 31 cm

REF-403 Neoprene Ankle Support - With Ligament



S	M	L	XL	XXL
21 cm 23 cm	23 cm 25 cm	25 cm 27 cm	27 cm 29 cm	29 cm 31 cm

Features

This anatomically designed product is made of neoprene material and covered with plastic sheets to increase stability. The yellow band inside the elastic provides additional support for eight diagonals. So with this support will reduce the pressure applied on ankle.

Indications

- In ligament instability and rupture
- Protective and anti -edema
- Used for protection or support after surgery and injuries

SL-916 Ankle Dorsi-Flexion Bandage



ADJUSTABLE

Features

Anatomically manufactured from neoprene material. Used for foot-drop deformity, the product keeps the ankle within the desired stretch in the dorsal flexion. Provides visible improvement in patient's gait.

Indications

- It is used for drop foot problems due to the various cases.
- Practical application.
- Can be used in any daily shoe.

SL-08 K Ankle Support - Padded



UNIVERSAL

Features

External surface is designed in accordance with anatomical structure, and inner surface is covered with sponge pad for maximum stabilization on ankle.

Indications

- Used for soft tissue damages and during post-operative period.
- Fits well to the foot thanks to the special design soft sponge pad.



SL-08 Ankle Support - Air Cast



UNIVERSAL

Features

External surface is designed in accordance with anatomical structure, and inner surface is covered with plastic airbag for maximum stabilization on ankle.

Indications

- Used for soft tissue damages and during post-operative period.
- The necessary compression can be adjusted thanks to the special air pad which covers well.
- The air pad can be blown up via the mini pump
- The adjustable air system provides great accordance with the foot.



SL-08J Ankle Support - Jelly



UNIVERSAL

Features


External surface is designed in accordance with anatomical structure, and inner surface is covered with gel pad for maximum stabilization on ankle.

Indications

- Used for soft tissue damages and during post-operative period.
- Gel pads provides being used with all kinds of shoes.



SL-903 Peroneal Orthosis



Features

Ankle orthosis manufactured from plastic material. Provides static dorsiflexion assistance and lateral stability for entire ankle/foot area. Provides easy usage with velcro fasteners going around the ankle. Use a scissor to trim insole for exact shoe fit. May be modified with a heat gun.

Warnings

- It is used for drop foot problems due to the various cases.
- Can be used with any kinds of shoes thanks to the anatomic design.
- Produced separately for right and left foot.

S	M	L
37 cm 38 cm	39 cm 40 cm	41 cm 42 cm

ERS-01 Foot Guard



S	M	L
35 cm 37 cm	38 cm 40 cm	40 cm 44 cm

Features

It is made of thermoplastic material the orthosis ankle band covered in a terycloth .the tape will prevent the movement and stabilize the ankle at 90 degrees, is preferred for hospitalized patients.

Indications

- Gives to the ankle 90-degree stability.
- Used in order to keep the ankle at 90 degrees.
- For the inpatients; thanks to the support with the form 'T', it prevents the ankle to turn.

ERS-01 D Foot Guard - Degree-Adjustable



S	M	L
35 cm 37 cm	38 cm 40 cm	40 cm 44 cm

Features

It is a T-band ankle orthosis made of thermoplastic material, covered with terry cloth. It provides stabilization of the ankle at different degrees and prevents rotation with the T band, is preferred for inpatients.

Indications

- Provides stability at different degrees.
- For the inpatients; thanks to the support with the form 'T', it prevents the ankle to turn.

SL-510 Rom Walker - Degree-Adjustable



S	M	L
35 cm 37 cm	38 cm 40 cm	40 cm 44 cm

Features

Helps to early mobilization and adjustment on dorsi and plantar flexions. Plantar and dorsi flexion angle adjustment possible at: 0°-7.5°-15°-22.5°-30°-37.5°-45. Rocker sole design encourages the natural roll of the foot during the walking motion. Available 2 sizes, long and short.

Indications

- Chronic ligament instability
- Post-op stabilization of tendon, ligament and soft tissue injuries
- Achilles tendon rupture
- Post-op or trauma for protection and rehabilitation purposes



SL-511 Achilles Boot



Features

It is an ankle orthosis with a base made of eva and an outer part made of hard plastic material and covered with a towel pad to assist early mobilization. Plantar flexion height at the ankle is adjusted by means of a 3-layer pad on the heel. Available 2 sizes, long and short.

Indications

- Chronic ligament instability
- Post-op stabilization of tendon, ligament and soft tissue injuries
- Achilles tendon rupture
- Post-op or trauma for protection and rehabilitation purposes



SL-521 Air Walker - Short



Features

Ankle orthosis designed to help maximum immobilization of the ankle. Its hard, thermoplastic structure provides maximum protection. Its practical pump system helps reducing edema and pain by applying desired amount of pressure around the ankle. The inner surface is manufactured from skin friendly, antibacterial auraterry fabric. Soft pads in the heel region provides maximum comfort. Extra hook and loop fasteners increase stabilization. The buckle system, which allows wearing the product without changing the adjustment of the bandages, provides ease of use. Nonskid sole with rocker system (convex design) supports walking by moving harmoniously with foot's natural oscillation. Special bamboo fiber socks given as a supplementary product provides thermostatic control and speeds up the healing process by balancing blood flow.

Indications

- For post-operative stabilization in cases of tendon, liga ment and soft tissue injuries
- Chronic ligament instability
- Metatarsus stress fracture and tarsal bone fracture treatment
- Simple navicular and cuneiform fractures
- Ankle sprain



SL-522 Air Walker - Long



Features

Ankle orthosis designed to help maximum immobilization of the ankle during treatment and healing process. Its hard, thermoplastic structure provides maximum protection. Its practical pump system helps reducing edema and increasing stabilization by applying desired amount of pressure. The inner surface is manufactured from skin friendly, antibacterial auraterry fabric. Soft pads in the heel region provides maximum comfort. Long part supporting the calf region and extra hook and loop fasteners increase stabilization. The buckle system, which allows wearing the product without changing the adjustment of the bandages, provides ease of use. Nonskid sole with rocker system (convex design) supports walking by moving harmoniously with foot's natural oscillation. Special bamboo fiber socks given as a supplementary product provides thermostatic control and speeds up the healing process by balancing blood flow.

Indications

- For post-operative stabilization in cases of tendon, ligament and soft tissueinjuries
- Chronic ligament instability
- Stable distal tibia and fibula fracture or distal tibia and fibula fracture where internal fixation is used
- Ankle sprain
- Stable foot and ankle fracture
- Achilles tendon tearing
- To prevent edema

Ankle Orthosis










UNIVERSAL

SL-02
Velpeau Bandage



UNIVERSAL

REF-620
Shoulder Support



UNIVERSAL

SL-04
Clavicle bandage



UNIVERSAL

ERSA-210
Wrist - Forearm Splint



UNIVERSAL

SL-13
Acromio-Clavicular Bandage



UNIVERSAL

ERSA-211
Epicondylitis Band



UNIVERSAL

ERSA-217

Cubital Tunnel Elbow Support



UNIVERSAL

ERSA-206

Hand & Wrist Splint - Single Sided



UNIVERSAL

ERSA-222

Elbow Support



UNIVERSAL

ERSA-207

Hand & Wrist Splint - Single Sided (With Thumb Support)



UNIVERSAL

ERSA-205


Dequervain Splint



UNIVERSAL

ERSA-208

Wrist Bandage



UNIVERSAL

ERSA-209

Wrist support - Simple



UNIVERSAL

SL-240

Abdominal Corset - Elastic



UNIVERSAL

ERSA-220


Boxer's Fracture Splint



UNIVERSAL

SL-241

Chest Corset



UNIVERSAL

SL-162

Posturex Corsert



UNIVERSAL

SL-260K


Lumbosacral Corset (26 cm)



SL-320K
Lumbosacral Corset (32 cm)




SL-330K
Dorsolumbar Corset



ERSA-215
Neoprene Abdominal Belt



ERSA-200
Neoprene Knee Brace - Closed



ERSA-216
Neoprene Lumbosacral Corset



ERSA-201 Neoprene Knee
Brace - Patella Support



ERSA-202

Neoprene Knee Brace - Pattela & Ligament Support



ERSA-219

Neoprene Calf Support



ERSA-203

Neoprene Knee Brace - Steel Hinged



ERSA-204

Neoprene Thigh Support



ERSA-218

Neoprene Thigh Support



ERSA-212 Neoprene Ankle Support - Cross Strap



ERSA-213

Neoprene Ankle Support - With
Malleolar Pad



ERSA-214

Neoprene Ankle Support - With
Ligament

Foot Health








SL-500 Silicone Insoles

Features

It is manufactured from silicone materials, allowing the body to a balanced distribution of the load on foot joints reduces the pressure and stress. Absorbs ground energy forces.

Indications

- Heel spina treatment
- Plantar fasciitis treatment
- Reduce the pain on heels
- Used for metatarsal pains



S	M	L	XL	XXL
34 cm 36 cm	36 cm 38 cm	38 cm 41 cm	41 cm 44 cm	44 cm 47 cm




SL-502 Silicone Insoles - Arch High

Features

It is manufactured from silicone materials, allowing the body to a balanced distribution of the load on foot joints reduces the pressure and stress. Absorbs ground energy forces.

Indications

- Heel spina treatment
- Plantar fasciitis treatment
- Reduce the pain on heels
- Used for metatarsal pains



S	M	L	XL	XXL
34 cm 36 cm	36 cm 38 cm	38 cm 41 cm	41 cm 44 cm	44 cm 47 cm




SL-500 Y Silicone Insoles - Yoyo

Features

It is manufactured from silicone materials, allowing the body to a balanced distribution of the load on foot joints reduces the pressure and stress. Absorbs ground energy forces.


Indications

- Heel spina treatment
- Plantar fasciitis treatment
- Reduce the pain on heels
- Used for metatarsal pains



S	M	L	XL	XXL
34 cm 36 cm	36 cm 38 cm	38 cm 41 cm	41 cm 44 cm	44 cm 47 cm

SL-504 Steel Insoles



34 / 35 / 36 / 37 / 38 / 39 / 40 / 41 / 42 / 43 / 44


Features

Insoles with longitudinal and metatarsal arch support. It's made of leather and supported by a stainless steel bridge

Indications

- Used for Heel spina
- Used for Dome's foot disorders
- Used to support the longitudinal arch

SL-506 Cork Insoles



34 / 35 / 36 / 37 / 38 / 39 / 40 / 41 / 42 / 43 / 44


Features

Arc support manufactured of special thermoplastic material which has leather covered with longitudinal and metatarsal support. Distributes the standing body weight equally and reduces load and stress on the joints.

Indications

- Make sure to use the correct size of product that matches your foot size.
- Please check the left and right orientation of the product before using.
- Insole should be placed at the sole part of the shoe.
- This medical tool is designed for only one patient.
- If you want to perform some little changes on the product for comfort or easy usage, please consult your physician or your dealer for technical support.
- Do not expose the product to extreme temperatures or flame.
- This product cannot be used for other purposes or combined with other product.
- It must be used only with original indications.
- Use only if prescribed by your physician.

SL-501 Silicone Heel Cup



S	M	L
35 cm 38 cm	38 cm 41 cm	42 cm 44 cm

Features

Made of silicone material. It relieves the load on the heel that causes pain.

Indications

- Reduce the load on heels
- Heel spina treatment
- Plantar fasciitis treatment
- Kallaneusun avascular necrosis
- Used to reduce the heel's pain




SL-501 Y Silicone Heel Cup - Yoyo

Features

Made of silicone material. It relieves the load on the heel that causes pain.

Indications

- Reduce the load on heels
- Heel spina treatment
- Plantar fasciitis treatment
- Kallaneusun avascular necrosis
- Used to reduce the heel's pain




S	M	L
35 cm 38 cm	38 cm 41 cm	42 cm 44 cm



SL-519 Silicone Medial and Lateral Wedge

Features

It is made of silicone coated fabric. According to the problem, it can be placed the inner or outer side of the shoe. Adhesive feature can be secured in the shoe. It is used in Calcaneal varus, calcaneal valgus, O-legs, X-Leg, meniscal degeneration.



S	M	L	XL
34 cm 36 cm	36 cm 38 cm	38 cm 41 cm	41 cm 44 cm




SL-505 Epin Heel Cup

Features

It is made of polyurethane, leather covering material. It relieves the load on the heel that causes pain.

Warnings


- Reduce the load on heels
- Heel spina treatment
- Plantar fasciitis treatment
- Kallaneusun avascular necrosis
- Used to reduce the heel's pain



1	2
34 cm 38 cm	38 cm 42 cm



SL-512 Silicone Arch Support - With Textile




Features
It is manufactured from silicone materials, allowing the body to a balanced distribution of the load on foot joints reduces the pressure and stress. Absorbs ground energy forces.

- Indications**
- Pes Cavus
 - Pes Planus

S	M	L	XL
34 cm 36 cm	36 cm 38 cm	38 cm 41 cm	41 cm 44 cm



SL-503 Silicone Toe Separator




Features
It is made of silicone material. It is suitable for day time use. Help to correct the Hallux valgus thumb disorder.

S	M	L
34 cm 37 cm	38 cm 41 cm	42 cm 45 cm



SL-516 Silicone Toe Ring




S	M	L
34 cm 37 cm	38 cm 41 cm	42 cm 45 cm



SL-509 Silicone Metatarsal Support

Features

Protection on the metatarsal area prevents the growth of hard areas and callouses. Reduces the pain and the irritation in this area. It is suitable for sports and daily use.




STANDART




SL-507 Silicone Bunion Protector

Features

It is made of silicone material to prevent friction on the bunion against the trauma. It is used to relieve sensitivity on the bunion.




STANDART



SL-514 Silicone Hammer Toe Pad

Features

It is made of silicone material. It is used for hammer toe.



STANDART



SL-517 Silicone Stress Ball

Features

Silicone hand therapy balls are produced in 3 colors.

BLUE - GREEN - RED



SL-518 Silicone Breast Prosthesis

Features

It is specially designed for breasts, with thick bra straps. It is specially prepared for women with right mastectomy, left mastectomy or both. There are pockets inside to slide into breast prostheses easily. Sizes between 80-130.

80 / 85 / 90 / 95 / 100 / 105 / 110 / 115 / 120 / 125 / 130



SL-925 Thermo Gel

Features

It is reusable. Compress Gel is the natural way to cure minor injuries, pains and pangs without using drugs. A large number of people including the athletes who have muscle strain, the patients who suffers inflammation, people who has painful and rigid joints can use the compress gel effectively.


It is very practical to get rid of the muscle and skin pains after a stressful and tiring day. It is also effective against insect bites and burns.



HV-10 Hallux Valgus Night Splint

Features

The inside is made of sponge fabric with towels material. The thumb prevents the lateral deviation and deformities. Hallux valgus are particularly suitable for night use.




RIGHT / LEFT



HV-11 Hallux Valgus Splint

Features

It's made of elastic material. This product is used to prevent or correct the deformities.




UNIVERSAL



SL-508 Plaster Shoe

Features

Shoe specially designed in accordance with the anatomical form of the foot for the purpose of providing ease of walking when the foot is put in a cast. It is manufactured from skin-friendly, antibacterial materials. Nonskid sole with rocker system (convex design) supports walking by moving harmoniously with foot's natural oscillation. The upper surface is totally made of breathable fabric. Specially designed hook and loop fastener allows it to be used on both feet.




S	M	L
33 cm 35 cm	36 cm 40 cm	41 cm 45 cm



HV-12 Hallux Valgus Shoe

Features

Post-operative shoe designed in accordance with the anatomical form of the foot for the purpose of supporting mobilisation after Hallux Valgus surgery. It is manufactured from skin friendly, antibacterial materials. Its nonskid, platform sole speeds up the healing process preventing pressure on the operated region.



The foot contact surface is soft. It does not put stress on metatarsus and phalanges, normalizes the extra stress on the heel as well. The upper surface is totally made of breathable fabric. Specially designed hook and loop fastener allows it to be used on both feet. It reduces possible post-operative edema on foot.



REF-410 Sand Weight

Features

Paralyzed and disabled patients with hand-arm-foot-leg muscle development in sport and exercise are used in physical therapy and rehabilitations.

500 gr / 1 kg / 1,5 kg / 2 kg / 2,5 kg / 3 kg / 3,5 kg

Knitted Group



supportline™







REF-619 Knitted Shoulder Support



S	M	L	XL	XXL
22 cm 24 cm	24 cm 27 cm	27 cm 30 cm	30 cm 34 cm	34 cm 38 cm

Features

The knitted shoulder support limits shoulder and arm movements and is individually adjustable with two different buckle positions. This allows the shoulder support to be ideally adapted to the user's needs. Thanks to its soft knit fabric design, the compression support is suitable to be worn under clothing. The smart system prevents the tapes in the armpit and chest area from narrowing. This increases the comfort of use under daily wear.

Indications

- Nonsurgical with anterior shoulder dislocation and chronic anterior subluxation
- Painful omarthrosis
- Postoperative after shoulder stabilisations.



REF-709 Knitted Elbow Support



S	M	L	XL	XXL
22 cm 24 cm	24 cm 26 cm	26 cm 28 cm	28 cm 30 cm	30 cm 32 cm

Features

Manufactured from respirable, flexible woven fabric, anatomically suitable for the shape of the elbow bandage is anatomically contoured for maximum fit, comfort, and effectiveness help to support injuries, improve blood circulation, and maintain body warmth hence accelerate the healing process. Provides perfect fit and continuous support.

Indications

- Arthritis
- Intra articular liquid increase
- Post-op or trauma treatment
- Provides protection and support for the elbow
- Support gives pain and stress relief for everyday use



REF-710 Knitted Elbow Support - With Epicondylitis Band



S	M	L	XL	XXL
22 cm 24 cm	24 cm 26 cm	26 cm 28 cm	28 cm 30 cm	30 cm 32 cm

Features

Knitted Elbow Support - With Epicondylitis. Band is an elastic elbow support with pads. It includes a breathable, elastic and light material, and anatomically-shaped knit. The support has a anatomically-shaped silicon pads for effective pressure relief and stabilisation.

Indications

- Epicondylitis (e.g. tennis or golfer's elbow)
- Osteoarthritis and arthritis

REF-722 Knitted Wrist & Hand Splint



S	M	L	XL	XXL
15 cm 17 cm	17 cm 19 cm	19 cm 21 cm	21 cm 23 cm	23 cm 25 cm

Features

Knitted Wrist & Hand Splint is a support with a compressive action that stabilises the wrist. Uniting compressive knitted fabric, a hand splint and an additional strap, the support provides stabilisation, relief and can reduce pain. The supple border prevents constriction. One of the benefits of the Knitted Wrist & Hand Splint is its high elasticity in the thumb area- this optimised comfort zone enables it to mould itself to the body.

Indications

- Carpal tunnel syndrome.
- Tendonitis, tenosynovitis, tenosynovitis issues.
- Used to reduce the hand wrist pain.
- Used for protection and rehaul after surgery and injury.

REF-730 Knitted Lumbosacral Corset



S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm

Features

Knitted Lumbosacral Corset is a lumbar support that can relieve pain in the lumbar spine through stabilisation. The integrated inner hand loop on the left side and the outer loop on the right make it easier to put the support on. The skin-friendly towelling material on the inner surface of the fastener makes it very comfortable to wear. In addition, special knit and fastener shapes in the groin area of the hip ensures a non-slip, correct fit of the Knitted Lumbosacral Corset.

Indications

- Lower back pain (lumbago)
- Muscular insufficiency
- Painful conditions of the lumbosacral junction

REF-731 Knitted Lumbosacral Corset - With Pad



S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm

Features

Knitted Lumbosacral Corset is a lumbar support that can relieve pain in the lumbar spine through stabilisation. The integrated inner hand loop on the left side and the outer loop on the right make it easier to put the support on. The skin-friendly towelling material on the inner surface of the fastener makes it very comfortable to wear. In addition, special knit and fastener shapes in the groin area of the hip ensures a non-slip, correct fit of the Knitted Lumbosacral Corset.

Indications

- Lower back pain (lumbago)
- Muscular insufficiency
- Painful conditions of the lumbosacral junction



REF-699 Knitted Knee Brace



S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm

Features

It is produced from breathable and non-sweating knitted fabric. Thanks to the ingenious knitting technology, it consists of a firm, compressive material. This gives the joint the required stability, yet is still easy to put on. The soft Comfort Zone with its rippled structure and diamond centre in the hollow of the knee provides the necessary comfort in wear, also when moving around.

Indications

- Slight instability of the knee joint.
- Osteoarthritis of the knee.
- Arthritis.
- Protection during sport activities.



REF-700 Knitted Knee Brace - With Patella & Ligament Support



S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm

Features


It is produced from breathable and non-sweating knitted fabric. Thanks to the ingenious knitting technology, it consists of a firm, compressive material. This gives the joint the required stability, yet is still easy to put on. The soft Comfort Zone with its rippled structure and diamond centre in the hollow of the knee provides the necessary comfort in wear, also when moving around. Medial (inside) and lateral (outside) spiral stays offer overall knee joint stability. Thanks to the open border of the patella silicone ring, the pressure on the kneecap is noticeably lower.

Indications

- Arthritis, bursitis
- Intra articular liquid increase
- Pain induced by knee cap
- Protection during and after the sportive activities



REF-701 Knitted Knee Brace - With Patella Support



S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm

Features


It is produced from breathable and non-sweating knitted fabric. Thanks to the ingenious knitting technology, it consists of a firm, compressive material. This gives the joint the required stability, yet is still easy to put on. The soft Comfort Zone with its rippled structure and diamond centre in the hollow of the knee provides the necessary comfort in wear, also when moving around. Thanks to the open border of the patella silicone ring, the pressure on the kneecap is noticeably lower.

Indications

- Slight instability of the knee joint.
- Osteoarthritis of the knee.
- Arthritis.
- Protection during sport activities



REF-702 Knitted Knee Brace - With Ligament Support



Features

It is produced from breathable and non-sweating knitted fabric. Thanks to the ingenious knitting technology, it consists of a firm, compressive material. This gives the joint the required stability, yet is still easy to put on. Medial (inside) and lateral (outside) spiral stays offer overall knee joint stability.

Indications

- Arthritis, bursitis
- Intra articular liquid increase
- Protection during and after the sportive activities

S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm



REF-703

Knitted Knee Brace - Steel Hinged

S	M	L	XL	XXL
31 cm 34 cm	34 cm 37 cm	37 cm 40 cm	40 cm 44 cm	44 cm 48 cm

Feature

It is produced from breathable and non-sweating knitted fabric. Thanks to the ingenious knitting technology, it consists of a firm, compressive material. This gives the joint the required stability, yet is still easy to put on. Medial (inside) and lateral (outside) two steel hinge stays offer overall knee joint stability.

Indications

- Arthritis, bursitis
- Intra articular liquid increase
- Protection during and after the sportive activities

REF-719 Knitted Ankle Support



S	M	L	XL	XXL
19 cm 21 cm	21 cm 23 cm	23 cm 25 cm	25 cm 27 cm	27 cm 29 cm

Features

Knitted Ankle Support is a support for compression of the soft tissue in the ankle joint. It includes a breathable, elastic and light material, and anatomically-shaped knit. The anatomical shape of Knitted Ankle Support enable an effective pressure relief and stabilisation. The support has an increased comfort by using softer material in the flexion area.

Indications

- Mild Ligament injuries
- Arthritis and osteoarthritis
- Post-operative recovery purpose

REF-720 Knitted Ankle Support - With Malleolar Pad



S	M	L	XL	XXL
19 cm 21 cm	21 cm 23 cm	23 cm 25 cm	25 cm 27 cm	27 cm 29 cm

Features

Knitted Ankle Support - With Malleolar Pad is a support for compression of the soft tissue in the ankle joint with pads. It includes a breathable, elastic and light material, and anatomically-shaped knit. The anatomically-shaped silicon pads enable an effective pressure relief and stabilisation. The support has an increased comfort by using softer material in the flexion area.

Indications

- Slight instability of the ankle joint
- Sprains and distortions
- Joint effusions and swellings due to arthritis and osteoarthritis

REF-721 Knitted Ankle Support - With Achilles Tendon Support



S	M	L	XL	XXL
19 cm 21 cm	21 cm 23 cm	23 cm 25 cm	25 cm 27 cm	27 cm 29 cm

Features

It is a sweat breathable mesh fabric, flexible and anatomically designed. Both sides of Achilles tendon are made from silicone as well as the heel support in order to prevent the imbalance between both legs. It is also used to reduce the tension on the tendon and to protect the heel and Achilles tendon.

Indications

- Achilles tendon, acute and chronic inflammation,
- Achilles tendon pain,
- For protection after surgery or injury.

Elastic Group


Elastic Group







ELS-01 Knitted Elastic Knee Support



S	M	L	XL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm

Features


Manufactured with flexible and antiperspirant woven fabric. Knee support provides resistance to a lateral blow when the knee is in full extension. Prevention of tendinous and articular accidents in sport activities and rheumatic knee diseases prophylaxis and calf injuries.

Supportline

Elastic Knee Brace are designed in order to obtain maximum comfort and benefit with regional support to the knee and it is available in 4 sizes.



ELS-02 Knitted Elastic Elbow Support



S	M	L	XL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm

Features

Manufactured with flexible and antiperspirant woven fabric. Gives full compression and support to the elbow area. Used for arthritic elbow conditions, minor sprains and strains. Provides support and relieves pain due to stiffness, swollen or injured elbows. Also protects elbow against further injury.

Supportline

Elastic Elbow Brace are designed in order to obtain maximum comfort and benefit with regional support to the elbow and it is available in 4 sizes.



ELS-03 Knitted Elastic Ankle Support



S	M	L	XL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm

Features


Manufactured with flexible and antiperspirant woven fabric, preventing of acute ankle sprains, reduces the risk of incurring of future injuries ,distuortion and bursitis.

Supportline

Elastic Ankle Brace are designed in order to obtain maximum comfort and benefit with regional support to the ankle and it is available in 4 sizes.



ELS-04 Knitted Elastic Wrist & Hand Support



S	M	L	XL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm

Features

Manufactured with flexible and antiperspirant woven fabric, preventing of acute hand sprains, reduces the risk of incurring of future injuries, distuortion and bursitis.

Supportline

Elastic Wrist Brace are designed in order to obtain maximum comfort and benefit with regional support to the wrist and it is available in 4 sizes.

Kids







XS
28 cm
31 cm

SL-05-KDS
Neck Collar (Sponge)



XS
28 cm
31 cm

ERS-112-KDS
Philadelphia Tracheotomy Collar



XS
28 cm
31 cm

ERS-110-KDS
Nelson Collar



S	M	L
75 cm	85 cm	95 cm
85 cm	95 cm	105 cm

SL-01-KDS
Arm Sling (Lux)



XS
28 cm
31 cm

ERS-111-KDS
Philadelphia Collar



S	M	L
75 cm	85 cm	95 cm
85 cm	95 cm	105 cm

SL-01F-KDS
Arm Sling (Mesh)



UNIVERSAL

SL-02-KDS
Velpeau Bandage



UNIVERSAL

ERSA-222-KDS
Elbow Support



UNIVERSAL

SL-04-KDS
Clavicle bandage



UNIVERSAL

REF-300-KDS
Neoprene Elbow Support



UNIVERSAL

REF-620-KDS
Shoulder Support



UNIVERSAL

ERSA-205-KDS
Dequervain Splint



UNIVERSAL

ERSA-206-KDS

Hand & Wrist Splint - Single Sided



UNIVERSAL

ERSA-209-KDS

Wrist support - Simple



UNIVERSAL

ERSA-207-KDS

Hand & Wrist Splint - Single Sided
(With Thumb Support)



UNIVERSAL

ERSA-211-KDS

Epicondylitis Band



UNIVERSAL

ERSA-208-KDS

Wrist Bandage



UNIVERSAL

SL-260-KDS

Lumbosacral Corset (26 cm)



UNIVERSAL

SL-161-KDS
Posturex Corsert



UNIVERSAL

SL-241-KDS
Chest Corset



UNIVERSAL

SL-162-KDS
Posturex Corsert



UNIVERSAL

SL-330-KDS
Dorsolumbar Corset



UNIVERSAL

SL-240-KDS
Abdominal Corset - Elastic



S	M	L
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm

SL-245
Umbilical Hernia Corset



UNIVERSAL

SL-513-KDS

Inguinal Hernia Belt



UNIVERSAL

ERSA-200-KDS

Neoprene Knee Brace - Closed



S	M	L
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm

SL-910

Pavlic Bandage



UNIVERSAL

ERSA-201-KDS

Neoprene Knee Brace - Patella Support



S	M	L	XL	XXL
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm	105 cm 115 cm	115 cm 125 cm

SL-909

Frejka Pillow



UNIVERSAL

ERSA-202-KDS

Neoprene Knee Brace - Patella & Ligament Support



UNIVERSAL

ERSA-218-KDS

Neoprene Thigh Support



UNIVERSAL

ERSA-213-KDS

Neoprene Ankle Support - With Malleolar Pad



UNIVERSAL

ERSA-219-KDS

Neoprene Calf Support



UNIVERSAL

SL-12-KDS

Immobilizer



UNIVERSAL

ERSA-212-KDS

Neoprene Ankle Support - Cross Strap



UNIVERSAL

REF-400

Eight Bandage - Elastic



S	M	L
75 cm 85 cm	85 cm 95 cm	95 cm 105 cm

ERS-01-KDS
Foot Guard



SHORT / LONG

SL-950
Dennis Browne Splint



12-13-14 28

SL-952
Dennis Browne Boot



SHORT / LONG

SL-951
Dennis Browne Splint - Moving

KiDS



Compression Stockings



supportline™





The diagram illustrates a human leg model with various segments and measurement points. The segments are labeled as follows: cT (thigh), cH (hip), K (knee), cG (groin), cF (femur), cE (upper shank), cD (lower shank), cC (calf), cB1 (heel), cB (ball), cY (arch), cX (heel), and cA (toe). The measurement points are labeled as follows: IGT (total leg length), IK2T (thigh length), IK1T (hip length), K1 (knee height), cG (groin height), cF (femur height), cE (upper shank height), cD (lower shank height), cC (calf height), cB1 (heel height), cB (ball height), cY (arch height), cX (heel height), and cA (toe height). The diagram also shows the segments of the foot: IX (heel), IA (arch), and IZ (toe).

Compression Stockings



UNDER KNEE
Varicose Stockings



OVER KNEE
Varicose Stockings



MATERNITY
Varicose Stockings



PANTYHOUSE
Varicose Stockings

ANTI - EMBOLISM STOCKINGS

Clinically designed to promote venous return for prevention of thrombosis before, during and after surgery and for medical conditions requiring prolonged bed rest and before, during and after child birth.

HOW TO CHOOSE SIZE

Two measurements are sufficient to choose the exact size. The ankle circumference (cb) and height of the leg (cg) then compare the measurements with the table. These measurements confirm the selected stocking size.

Do not use the stocking if the measurement does not match with the table.



Intelligent Topband

The specially developed intelligent silicone topband ensures a non-slip fit and provides information on the correct fit at a glance.



Colour Coded Heel

The colour-coded heel with the velour inner surface reduces the danger of pressure sores, simplifies sorting, and shows that the heel is properly fitted.



Inspection Opening

The unobtrusive inspection opening on the sole of the foot allows staff to check the circulation and reduces the risk of falling.



Hygienic Toe

The antibacterial and antimycotic toe reduces normal hospital micro-organisms—even after repeated washing.

Knee (AD)

Size	Leg measurements* (circumference in cm)				Colour code	
		Ankle cB	Calf cC	Below knee cD	Thigh cG	
XS	(1)	17 - 19	30 - 35	28 - 34	-	natural
S	(2)	19 - 21	31 - 37	29 - 36	-	red
M	(3)	21 - 23	32 - 40	30 - 39	-	blue
L	(4)	23 - 25	33 - 42	31 - 41	-	yellow
XL	(5)	25 - 27	34 - 44	32 - 43	-	orange
XXL	(6)	27 - 29	35 - 46	33 - 45	-	black
XXXL	(7)	29 - 31	36 - 48	34 - 47	-	brown

Thigh (AG)

XS	(1)	17 - 19	30 - 35	28 - 34	42 - 54	natural	natural
S	(2)	19 - 21	31 - 37	29 - 36	46 - 58	red	red
M	(3)	21 - 23	32 - 40	30 - 39	50 - 62	blue	blue
L	(4)	23 - 25	33 - 42	31 - 41	54 - 66	yellow	yellow
XL	(5)	25 - 27	34 - 44	32 - 43	58 - 70	orange	orange
XXL	(6)	27 - 29	35 - 46	33 - 45	62 - 74	black	black
XXXL	(7)	29 - 31	36 - 48	34 - 47	66 - 78	brown	brown

Compression Stockings



ANTI-EMBOLI UNDER KNEE
Compression Stockings



ANTI-EMBOLI ARTICULATED ARCH
Compression Stockings - Pair



ANTI-EMBOLI OVER KNEE
Compression Stockings

Pillows



supportline™








S - M - L

SL-22

Sleeping Pillow



S - M - L

SL-22J

Sleeping Pillow - With Gel



STANDARD

SL-922

Reflux Pillow



STANDARD

SL-923

Visco Snoring Pillow



STANDARD

SL-06V

Visco Lumbar Cushion



STANDARD

SL-06Y

Anatomic Lumbar Cushion



STANDARD

SL-07V
Visco Seat Cushion



STANDARD

SL-07A
Seat Cushion



STANDARD

SL-07Y
Anatomic Seat Cushion



STANDARD

SL-07K
Square Seat Cushion



STANDARD

SL-23
Travelling Neck Pillow



STANDARD

SL-921
Patient Rotation Pillow



Supportline brand is known as best quality in orthopedic medical product during these years of manufacturing and had the advantage of having diversity in mass production as per customer request. Our goal is to keep on expanding and developing our brand quality in accordance with the latest technology and international standards and trust to our loyal customers.

Our unique Supportline liposuction products present below benefits as per customer reviews;

- ✔ Made from the latest and greatest elastic fabrics, with versatile, breathable and excellent fabric returning.
- ✔ Accelerates the healing process and reduces the risk of complications.
- ✔ By applying equal pressure across the treatment area, it prevents excessive fluid accumulation and helps your body absorb the accumulated fluids.
- ✔ Helps to soften, flatten and reduce the visibility of scars.
- ✔ All these features provide a perfect fit and comfort throughout the recovery period.

Smoother, firmer, desirable results.



LIPOSUCTION

GARMENTS



Healthy Beauty

Please ask for detailed Liposuction catalogue!

Notes

[illegible]