

Leinzett



Single-Use Appendixscope

Application: Observation, diagnosis, and treatment of appendiceal endoscopic surgery

Single-Use Appendixscope



Features

- Real-time scale identification for rapid intraoperative measurement and positioning
- No need for a host machine - flexible processor
- No power required - lightweight and portable
- Compatible with tablets, smartphones, laptops, and hospital monitors
- Leinzett Software (App) for video recording and screenshots - easy installation on tablets and laptops



→ Specifications

Single-Use Appendixscope

Model	Working Length	Viewing Angle	Diameter	Working Channel
EA10	2300mm	120°	2.6 mm/7.8Fr	≤1.2 mm/3.6Fr
EA14	2300mm	120°	3.5 mm/10.5Fr	≤1.6 mm/4.8Fr
EA18	2300mm	120°	3.5 mm/10.5Fr	≤2.0 mm/6.0Fr

Accessories

Basket

Model	Diameter	Length	Basket opening width	Basket length	Working channel	Basket Features
SE1027B/L	1.0mm	2700mm	12mm	16/20mm	≥1.2mm	Semi-threaded / fully-threaded
SE1527B/L	1.5mm	2700mm	12mm	16/20mm	≥1.6mm	Semi-threaded / fully-threaded

Biopsy forceps

Model	Diameter	Length	Opening width	Working Channel	Head Characteristics
BF1027T	1.0mm	2700mm	4.0	≥1.2	With Spike
BF1529	1.5mm	2900mm	4.6	≥1.6	Without Spike

Leinzett Medical Europe GmbH

Rheinfeldweg 1,40667 Meerbusch,Germany

www.leinzett.com

info@leinzett.com



Leinzett