

# Redefine mobility.

The first mechanical standing wheelchair that allows you to move while standing.

Get up and go!

Contact us

 (Incubase) De Hems 10, 7522 NL  
Enschede, The Netherlands

 +31655577682

 gabriel@j58.eu

 @j58\_thenextwheelchair

 www.j58.eu

## The problem with other **standing** wheelchairs.



### Expensive

electric

> €15 - 30K+ and often not reimbursed by insurance.

### Very Heavy

electric

70kg+ and often not dismountable. Need for an adapted vehicle for transportation.

### Battery dependent

electric

Usage of motors, among other complex components.

### Slow transition Sit-stand

mechanical

Slow transitioning between sitting and standing.

electric


### Can't move while standing

mechanical

Mechanical wheelchairs don't have the function to move while standing.



# THE EZZER





# The **EZER** standing wheelchair is different.

## These are the reasons why

Instant Standing & Sitting

Move While Standing

Mechanical Design

Additional Safety



### Smooth transitioning

Allows users to go Effortlessly from sitting to a standing position.



### Move While Standing.

Allows users to move around while in a standing position.



### Safety

The front wheels extend outward, expanding the wheelbase to ensure safety, stability, and comfort while standing.



### Fully Mechanical

A standing wheelchair that's as easy to repair as a bicycle.



# Experience eye level conversations.

For individuals with paraplegia or other mobility challenges, engaging in **eye-level** conversations can be important to maintain a sense of **inclusion** in everyday interactions. Incorporating the use of a wheelchair during rehabilitation supports a smoother **adjustment** to their evolving lifestyle.







# Faster adaptability. More mobility.

Ezer's ability to move while standing helps individuals **transition** more easily from rehabilitation activities to daily activities, promoting greater **autonomy** and **motivation** during rehabilitation sessions.

## Enhanced physical therapy.

Moving while standing activates various muscle groups, contributing not only to improved **core strength** but also to better **bone health**, while reducing the risk of secondary injuries like pressure sores.







## Versatile and cost effective.

A mechanical wheelchair with standing and moving capabilities **reduces the need** for multiple specialized devices, **saving space** and **costs** for rehabilitation centers while offering versatile therapy options for diverse patient needs.



# We Offer Plans Tailored to Your Needs

- **Book a demo:** Test the wheelchair in your facility.
- **Lease Option:** Access an affordable and flexible leasing plan.
- **Pre-Order:** Secure your standing wheelchair(s) and place an order today.



Trusted by

**NovelT**



**umec**

**MEDTECH  
TWENTE**

**UNIVERSITY  
OF TWENTE.**



# Get in touch with us



Gabriel Costa

CEO &  
CO-Founder

✉ [gabriel@j58.eu](mailto:gabriel@j58.eu)



Rivelino Montenegro

COO &  
CO-Founder

✉ [rivelino@j58.eu](mailto:rivelino@j58.eu)



Mauricio Goncalves Tikovsky

VP of Business  
Dev & Design

✉ [sales@j58.eu](mailto:sales@j58.eu)



[www.j58.eu](http://www.j58.eu)