

Warm Touch 6000

Warming Unit

A convenient, cost-effective solution to help maintain normothermia

The WarmTouch convective warming system provides veterinary clinicians with an easy solution to help maintain normothermia. The WarmTouch convective warming unit has an intuitive screen showing the temperature and status of each of the 5 temperature settings. The soft but strong WarmTouch blanket range is tailored for most surgical procedures, offering tear resistance, even airflow and patient comfort.

FEATURES

- Easy to use: New intuitive LCD screen and menu function provide all information at a glance.
- Easy to handle: The low weight makes it easy to transport and attach to a cart, bed rail or an IV pole.
- Easy to control: 5 different settings, including ambient and boost, for more flexibility and precision.
- Easy to service: The onscreen filter expiration warning makes it easy to maintain and service.



2608-2022-06-14

SPECIFICATIONS



Weight with Power Cord:
11.5 lbs



Power Requirements:
100 to 240 V AC

Max. Current at 100 V:
8 A

Max. Current at 240 V:
5 A

Input Frequency:
50/60 Hz \pm 1 Hz



Maximum Contact
Surface Temperature:
44.1°C (111.4°F)

Average Time for Contact
Surface Temperature to Rise:
6 minutes

Average Time for
Temperature of Air
Exiting the Hose to Rise:
< 1 minute

Accuracy of Displayed
Temperature:
 \pm 1.0°C (air entering hose)

Automatic Temperature
Stepdown:
After 45 minutes of continuous
use, blower will step down
from the Boost to High setting