



## B0521A 4 PLY LINT FREE MEDICAL BLACK FACE MASK (Earloop)

4 Ply Lint Free Medical Facemask with Earloop is now available in a sleek and sophisticated black color. This facemask stands out not only for its elegant appearance but also for its premium quality material. It is made from Lint Free Fabric which is designed for extended wear, providing reliable protection without sacrificing comfort.



### Specialized For :



F & B



INDUSTRIAL




HEALTHCARE



### PHYSICAL DIMENSION

LAYER	FEATURE
Outer Layer	Fluid Resistance Layer
Third Layer	Polypropylene Spunbond Non-Woven
Second Layer	Meltblown Filtration Layer
Inner Layer	Lint-free and Water Absorbent

	Size
Adult	17.5cm X 9.5cm
Kid	14.5cm X 9cm



### Product Specification & Instruction Use

#### PRODUCT SPECIFICATIONS:

**Classification:**  
ASTM Level 3  
MDA No. : GA 43088421-58036  
Cytotoxicity Test Report  
Fluid Resistance Test Report  
Dermal Sensitization Test Report  
Animal Irritation Test  
BFE & PFE Final Report

**Product Material:**  
Lint Free

**Packaging Mode:**  
30pcs/box

**Physical Features**

- Hypoallergenic
- Dual Band Adjustment
- Quad Fold
- Single Use Only

#### STORAGE:

- Store in dry, clean conditions in original packaging.
- Store away from direct sunlight, source of high temperature and solvent vapors.
- Store with the temperature range 0°C to 39°C.
- Shelf life is 60 months from the date of manufacture when stored as stated above.
- Replace if damaged or heavily contaminated or in accordance with local practice.
- Dispose of contaminated face mask with care and in accordance with national regulations.

#### ASTM F2100 MASK PERFORMANCE

BFE(%)	≥98%
Delta P, (mm H <sub>2</sub> O/cm <sup>2</sup> )	<6.0
PFE, @0.1μ	≥98%
Fluid Resistance (mmHg)	160
Flame Spread	Class 1

LEVEL  
3

#### DISCLAIMER

- Do not reuse
- Stay away from fire
- Use with caution if allergy with non-woven
- Discard the mask directly into the bin to avoid cross contamination
- Wear the face mask in right way and do not touch the inner layer with hand