

# VIBRANT

Vibrant 4 Ply Lint Free Medical Facemask is the perfect fusion of safety and style medical face mask. It consists of five vibrant colours such as Neon Green, Army Green, Purple, Neon Pink, and Yellow. Now, you can express yourself by changing the different vibrant colours facemask while staying protected. Vibrant is made from hypoallergenic, skin-friendly, and irritation-reducing lint-free fabric. Individually packed for maximum hygiene, it's your go-to choice for hospitals, clinics, labs, and public settings where social distancing isn't always possible. Spread Good Vibes and Enjoy Life with every breath.



## Specialized For :



F & B



INDUSTRIAL




HEALTHCARE



## PHYSICAL DIMENSION

LAYER	FEATURE
Outer Layer	Fluid Resistance Layer
Third Layer	Polypropylene Spunbond Non-Woven
Second Layer	Meltblown Filtration Layer
Inner Layer	Lint-free and Water Absorbent

	Size
Adult	17.5cm X 9.5cm



## COLOUR



PURPLE



NEON PINK



YELLOW



NEON GREEN



ARMY GREEN

## Product Specification & Instruction Use

### PRODUCT SPECIFICATIONS:

#### Specification:

MDA  
ASTM LEVEL 3 Test Standard: ASTM F2100  
EU Medical Device Regulation : 2017/745  
BFE & PFE Report Tested

#### Product Material:

Lint-Free

#### Packaging Mode:

30pcs/box (Individual Packaging)

#### Physical Features

- Available in 5 vibrant colours
- Elastic and stretchable round earloop
- Hypoallergenic
- High filtration capacity (BFE  $\geq$  98%)
- Biocompatibility Tested

### STORAGE:

- Store in dry, clean conditions in original packaging.
- Store away from direct sunlight, source of high temperature and solvent vapours.
- Store with the temperature range 0°C to 39°C.
- Shelf life is 60 months from the date of manufacture when stored as stated above.
- Replace if damaged, or heavily contaminated, or in accordance with local practice.
- Dispose of contaminated face mask with care and in accordance with national regulations.

### ASTM-F2100 Mask Performance

BFE (%)

$\geq$  98%

$\Delta$ P (mmH<sub>2</sub>O/cm<sup>2</sup>)

< 6.0

PFE,@0.1 $\mu$ (%)

$\geq$  98%

Fluid Resistance (mmHg)

160

Flame Spread

Class I



### Disclaimer

- Do not reuse
- Stay away from fire
- Use with caution if you are allergic to non-woven
- Discard the mask directly into the bin to avoid cross-contamination
- Wear the face mask in the right way and avoid touch the inner layer with your hand

