








<p>9</p>	<ul style="list-style-type: none"> ● Ultrasound therapy ● Transcutaneous Electrical Nerve Stimulation (TENS) ● Combination therapy (Medium frequency Interferential therapy + Ultrasound therapy) 		<p>The device must be Enraf Nonius, Sonopuls 492 or equivalent. Must be able to apply ultrasound therapy, electrotherapy and combination therapy. Must be able to adjust frequencies between 1 MHz and 3 MHz ultrasound, continuous and pulsed ELECTROTHERAPY Current channels 2 (100% independent) COMBINATION THERAPY Pre-programmed treatment suggestions 2 - Evidence Based Programmable positions 1000+ Must have a touch screen operating panel. Must be able to select the correct parameters with the integrated evidence-based protocols. Must be able to use it as a mobile device and must have the battery and with the mains power supply</p>		<p>3</p>
----------	--	---	--	--	----------






10	<p>Electrical stimulator</p>		<p>M 4-channel unit for electrotherapy Modifiable protocols with pathology indication. Memorized programs with the possibility of sequential, simultaneous or progressive work. Cyclic sequential programs. It involves the patient in the treatment, giving him the possibility of regulating the intensity, types of medium and low frequency currents. Technical characteristics: TENS symmetric and asymmetric two-phase pulsed current: Asymmetric phase duration (pulse width): 10-200 μs. Symmetric phase duration (pulse width): 10-1000 μs ..Pulse rate: 1-200 Hz. Modulation frequency: 1-180 Hz. Modulation programs: 1 / 1s, 6 / 6s, 12/12, 1 / 30s. Burst frequency: 0-9 Hz.</p>		2
----	-------------------------------------	---	--	--	---



11	<p>Transportation cart for the electrotherapy machines</p>		<p>Transport cart for medical devices 1 drawer 2 tray 4 castors with 2 castors brakes Working height, both in standing and in sitting position, so that the parameters on the display of the electrotherapy-units are clearly visible at any time. All wiring and cables can be hidden neatly in the hollow frame.</p>	5
12	<p>Xray Viewer</p>		<p>X-Ray Double Film Viewer, medical X-ray illuminator (view box) and edge lit. Suitable for viewing MRI/CT and X-ray film. Illuminator with a viewing area of 28" x 17" suitable for viewing two films. Illuminator must have LED lamps and has inbuilt universal ac power supply and the entire panel brightness must be controlled with a rotary dimmer, the maximum brightness for this model is > 10,000 lux.</p>	1



13	Disarticulated Human Bone Set		<p>Bilateral Disarticulated Human Skeleton with Skull Made of High Quality PVC Plastic Life Size 170cm tall With Impressions</p>	1
14	Hospital bed screen		<p>Hospital Bed Side Screen, four-Fold Frame work of bed side screen constructed with ms steel tubes. Overall size of bed side screen: 244cm L x 165 cm H (opened) : 60 cm L 165cm (closed) Removable Feet Pre-treated and epoxy powder coated finishing</p>	6
15	Dumbbells set with wall mounting rack		<p>Vinyl dumbbells – 100% iron covered with vinyl – resistant to sweat and moisture. 0.5, 1 kg, 2kg, 3kg, 4kg-5kg (1 pair each) Must have the wall mountable rack for the dumb bells.</p>	1

