

STEENBEEK FOOT ABDUCTION BRACE (SFAB)

TYPE II PATTERNS




WHAT ARE THE TYPE II PATTERNS

Most babies treated for clubfeet fit in the SFAB Type I. However, some babies with clubfeet have very fat and puffy feet. The feet are broad and short and they do not fit easily in the regular SFAB. For these babies there is the SFAB Type II. All the steps in the production of the SFAB remain the same, but the shoe patterns are slightly different to accommodate this type of foot.

1 cm

Size **6** SFAB patterns

TYPE II



1 = upper leather

2 = heel strap

3 = tongue


4 = heel counter

HPM
103/2004


1 cm

Size **6** SFAB patterns


TYPE II



5



6



7

5 = outersole (plywood)

6 = sock lining


7 = sole (rubber)

RP
03/2009

1 cm

Size **7** SFAB patterns


TYPE II




1



2



3



4

1 = upper leather

2 = heel strap

3 = tongue




4 = heel counter

SM
03/2009

1 cm

Size **7** SFAB patterns

TYPE II



5 = outersole (plywood)

6 = sock lining




7 = sole (rubber)

Handwritten signature and date: 03/2009

1 cm

Size **8** SFAB patterns

TYPE II



1 = upper leather

2 = heel strap

3 = tongue




4 = heel counter

SA
03/2009

1 cm

Size **8** SFAB patterns

TYPE II



5 = outersole (plywood)




6 = sock lining

7 = sole (rubber)

SA
03/2009

1 cm

Size **9** SFAB patterns TYPE II



1 = upper leather

2 = heel strap

3 = tongue




4 = heel counter

MS
03/2009

1 cm

Size **9** SFAB patterns

TYPE II



5 = outersole (plywood)

6 = sock lining




7 = sole (rubber)

SA
03/2009

1 cm

Size **10** SFAB patterns

TYPE II



1 = upper leather

2 = heel strap

3 = tongue




4 = heel counter

MA
02/2009

1 cm

Size **10** SFAB patterns

TYPE II



5 = outersole (plywood)

6 = sock lining




7 = sole (rubber)

SM
03/2009

1 cm

Size **11** SFAB patterns

TYPE II



1 = upper leather

2 = heel strap

3 = tongue




4 = heel counter

HS1
03/2009

1 cm

Size **11** SFAB patterns

TYPE II



5 = outersole (plywood)

6 = sock lining




7 = sole (rubber)

LP
03/2009

1 cm

Size **12** SFAB patterns

TYPE II



1 = upper leather

2 = heel strap

3 = tongue




4 = heel counter

LSM
03/2009

1 cm

Size **12** SFAB patterns

TYPE II



5 = outersole (plywood)

6 = sock lining





7 = sole (rubber)

LS1
03/2009

1 cm

Size **13** SFAB patterns

TYPE II



1 = upper leather

2 = heel strap

3 = tongue




4 = heel counter

Handwritten: LSP
03/2009

1 cm

Size **13** SFAB patterns

TYPE II



5 = outersole (plywood)

6 = sock lining





7 = sole (rubber)

SP
03/2009

1 cm

Size **14** SFAB patterns

TYPE II



1 = upper leather

2 = heel strap

3 = tongue


4 = heel counter

LSM
03/2009


1 cm

Size **14** SFAB patterns


TYPE II



5



6



7

5 = outersole (plywood)

6 = sock lining

7 = sole (rubber)

HP
03/2009

The Steenbeek Foot Abduction Brace (SFAB) is a device to maintain the extreme foot abduction (70 degrees) in children with treated clubfeet. The SFAB forms an integral part of the Ponseti Method of Clubfoot treatment, which targets new-born babies with Clubfeet.

Correction of the Clubfoot is achieved by applying serial casts of p.o.p., in combination with a minor surgical procedure (percutaneous tendo achilles tenotomy for 90% of the babies) under local anaesthesia. The SFAB is necessary to maintain this correction.

Global HELP Organization is accumulating publications that are available on our web-site for free and can be printed at minimal cost. Please visit our web-site at www.global-help.org or at www.orthobooks.org.

Publications


All publications from our web-site are free. Publications are available in several forms.

PDF for Web site and CD libraries

All publications are available in PDF format. PDF files can be downloaded from our web-site by clicking on the publication title or image. These files are then copied on your computer and can be printed on personal printers in color or black-and-white. CD libraries include many publications in one convenient compact disc. Some CD libraries are restricted to be used only in developing countries.

Printed publications

Some publications are available for “on demand” printing. Order single or multiple printed copies on our web-site at www.global-help.org. These publications are available for only the cost of printing and mailing.



Web site addresses:

www.global-help.org
www.orthobooks.org

**Send comments, questions
or photographs of a HELP
publication in use to:**
questions@global-help.org

Copyright © 2009 Global HELP Organization
Originally printed in 2004, “Steenbeek Brace For Clubfoot” was updated in 2006.

ISBN 978-1-60189-068-9



9 781601 890689



GLOBAL HELP
HEALTH EDUCATION USING LOW-COST PUBLICATIONS