



Support Essentials

User Guide

User Configuration for EMIS Web


Introduction

This guide shows you how to configure users in **EMIS Web** to enable access to **DXS Point-of-Care** and ensure correct login permissions.


Please note: user configuration is required for each user who needs access to DXS Point-of-Care.

Instructions


1. Log into **EMIS Web** and open the '**EMIS**' menu.




2. Hover over 'System Tools'.




3. And select 'EMAS Manager'.




4. In the left-hand menu, select 'Partner API'.




5. A list of available APIs will be displayed.




6. Scroll down and select 'DXS Point-of-Care'.




7. Click 'Edit Users'.




8. A list of users will appear.




9. You can sort the list by username by clicking the 'User' column.



10. Tick the checkbox next to the user who requires access to DXS Point-of-Care.




11. Enter 'Password1' in both password fields, then click 'OK' to save.




12. Ensure 'DXS Point-of-Care' is still highlighted, then select 'Login Access'.




13. Another user list will appear.




14. You can sort this list by clicking the 'User' column.



15. Tick both 'Auto Login' and 'Allow Login' for the relevant user.



16. Click 'OK' to save.




17. This enables access and registers the user for auto-login.

Installation Note:

18. When DXS Point-of-Care is installed on a workstation:

- The user will be prompted to **activate** DXS Point-of-Care.
- A Windows Firewall prompt may appear – select '**Allow access**' if prompted.



- If DXS Point-of-Care does not load:
 - Ensure it has been activated on that workstation.
 - Restart the device.

Support



0800 028 0004 (option 1)



support-uk@dxs-systems.com



dxs-systems.co.uk

Training



0800 028 0004 (option 2)



training@dxs-systems.co.uk



training.dxs-systems.com

