



STABILITY INDICATING RP-HPLC METHOD FOR THE DETERMINATION OF PIPERACILLIN AND TAZOBACTAM AND THEIR RELATED SUBSTANCES IN BULK AND PHARMACEUTICAL FORMULATION

Khyathi Navle^{I*}, B. Lakshmi Prasanna^{II}, Hema Sree Dometti^I and Shankar Moodu^I

^{I, II}Department of Pharmaceutical Analysis, Chalapathi Institute of Pharmaceutical Sciences, Guntur – 522034, AP, India.

Article Received on
05 June 2017,

Revised on 26 June 2017,
Accepted on 16 July 2017,

DOI: 10.20959/wjpps20178-9829

*Corresponding Author

Khyathi Navle

Department of
Pharmaceutical Analysis,
Chalapathi Institute of
Pharmaceutical Sciences,
Guntur – 522034, AP,
India.

ABSTRACT

A simple, sensitive and precise high performance liquid chromatographic method was developed for the determination of related substances of Piperacillin and Tazobactam in bulk and pharmaceutical formulation. The Impurities were well separated on a Hypersil ODS column (250mm X 4.6mm, 5 µm) by the gradient program using Phosphate buffer (pH 6.0) and Acetonitrile at a flow rate of 1.0 mL /min with detection wavelength at 220 nm. The developed method separates Piperacillin and Tazobactam from its related substances by 45 minutes run time. The calibration curve shows linearity in the concentration range of 0.3-75.5 µg/ml for Piperacillin, 0.3-3.3 µg/ml for Tazobactam, 0.4-3.1 µg/ml for Aminosulfinic acid, 0.4-17.2 µg/ml for Ampicillin, 0.4-7.9 µg/ml for D-phenylglycyl

EDP, 0.4-15.8 µg/ml for Piperacillin sulfoxide with a correlation coefficient 0.9997, 0.9993, 0.9998, 0.9999, 0.9992, 0.9996 respectively. The method was validated for specificity, linearity, accuracy, precision, robustness and solution stability. Recovery was found to be in the range of 75-110% and RSD is <2%. The proposed method was successfully applied for the determination of Piperacillin and Tazobactam and their related substances in bulk and pharmaceutical formulations.

KEYWORDS: Piperacillin(PIPER), Tazobactam(TAZO), HPLC, Related substances, Method Development.

INTRODUCTION

Piperacillin is a broad spectrum antibiotic, chemically it is [2S-[2a,5a,6b(S*)]]-6-[[[(4-ethyl-2,3-dioxo-1-piperazinyl) carbonyl] amino] phenyl- acetyl] amino-3,3 dimethyl-7-oxo-4-thia-1-azabicyclo-[3,2,0]heptanes-2-carboxylic acid. It exerts bactericidal activity by inhibition of both septum and cell wall. Whereas Tazobactam is a β -lactamase inhibitor and its chemical formula is (2S,3S,5R)-3-methyl-7-oxo-3-(1H-1,2,3-triazol-1-ylmethyl)-4-thia-1-azabicyclo-[3.2.0]heptanes-2-carboxylic acid-4,4-dioxide. In vitro, piperacillin and tazobactam combination is active against a variety of gram-positive and gram-negative aerobic and anaerobic bacteria. Ampicillin(<0.2% w/W), D-phenylglycyl EDP (<0.1% w/W), Piperacillin sulfoxide(<0.2% w/W) are the major impurities of piperacillin and Aminosulfonic acid(<1.0% w/W) is an impurity of tazobactam.

Several analytical methods were reported for the simultaneous estimation of Piperacillin and Tazobactam include second-order derivative spectrophotometry, UV-spectrophotometric absorption correction method, HPLC. Literature survey revealed that very few methods were reported for the simultaneous estimation of PIPER and TAZO and their related substances. So, an attempt has been made to develop an accurate, precise and economic RP-HPLC method for the determination of related substances(1-6).




Fig.1 Structure of Piperacillin




Fig.2 Structure of Tazobactam




Fig.3 Structure of Ampicillin




Fig.4 Structure of Aminosulfonic acid




Fig.5 Structure of D-phenylglycyl EDP




Fig.6 Structure of Piperacillin sulfoxide

EXPERIMENTAL

MATERIALS AND METHOD

Chemicals

Working reference standards of Piperacillin, Tazobactam acid and reference standards of ampicillin, aminosulfonic acid, D-phenyl glycel EDP, Piperacillin sulfoxide were gift samples from Aurabindo Pharma Pvt., Ltd., Hyderabad. Orthophosphoric acid was purchased from Lancaster, 40% Tetrabutylammonium hydroxide aqueous solution, Acetonitrile, Water, sodium hydroxide, hydrochloric acid, hydrogen peroxide were purchased from Merck

specialities private limited, Mumbai. Piperacillin and tazobactam for injection- 4.5 gm was commercially purchased.

Equipment

Alliance waters 2695 HPLC separation module provided with UV/Visible detector was used. Data acquisition was carried out by using Empower software. Separation of compounds was performed on Hypersil ODS C₁₈ (250mm×4.6mm, 5 μ).

Chromatographic parameters

Mobile phase consists of solution A- a mixture of phosphate buffer: acetonitrile (98:2%v/V) was prepared, sonicated for 5 min and filtered through membrane filter (0.45 μm) then P^H adjusted to (6±0.05) using NaOH solution. Solution B- Acetonitrile was alone used as solution B. Flow rate was 1.0 ml/min and injection volume was 20 μl. Run time was 45 min with a detection wavelength of 220 nm. Column temperature was maintained at 30⁰C.

Preparation of Diluent

A mixture of phosphate buffer solution and Acetonitrile in the ratio of 98:2%v/V with P^H adjusted to (6±0.05) was used.

Preparation of stock and standard solution

Stock solutions of Piperacillin and Tazobactam were prepared by dissolving 10 mg of PIPER in 10 ml of diluent and 6 mg of TAZO in 10 ml of diluents respectively. From the above stock solutions 5 ml of PIPER and 3ml of TAZO were solutions were spiked in 100 ml volumetric flask and volume was made upto the mark with diluent to prepare a standard solution with concentration range of 18μg/ml for TAZO and 50μg/ml for PIPER.

Preparation of Test solution

285 mg of dry powder (for injection) was weighed and transferred to 50 ml volumetric flask, to this 5 ml of acetonitrile was added and sonicated. Volume was made upto 100 ml with diluents and filtered through 0.45 μm or finer porosity membrane filter.

VALIDATION OF THE HPLC METHOD

System suitability

System suitability is commonly used to verify resolution, column efficiency, and repeatability of the chromatographic system to ensure its adequacy for a particular analysis. Standard solution of Piperacillin and Tazobactam(system suitability solution) was prepared and

injected into HPLC system in five replicates and chromatograms were recorded. The system suitability parameters, %RSD and calculated data was presented in Table-1.

Precision

The precision of an analytical method is the degree of agreement among individual test results when the method is applied repeatedly to multiple sampling of homogeneous sample. Six preparations individually using single batch of PIPER and TAZO for injection spiked with related substances at 100% of their specification limit were prepared and injected into HPLC. Percentage RSD for peak areas of all components were calculated. Less than 10% RSD indicates the method was precise and data represented in Table-2.

Accuracy

Accuracy of an analytical method is the closeness of test results obtained by that method to the true value. To ascertain accuracy solutions were prepared in triplicate using PIPER and TAZO injection spiked with related substances at levels 50%, 100% and 150% of specification as per proposed method and injected each solution into HPLC. The percentage recovery and RSD were calculated and the data was tabulated in Table-3.

Limit of Quantification and Limit of Detection

Limit of Quantification and Limit of Detection values of each of related substances were predicted from a separate linearity data performed at a lower concentration (below 50% specification level). Each predicted concentration was verified by preparing the solutions and injecting each solution six times into HPLC. Percentage RSD for PIPER, TAZO and related substances were calculated and data was given in Table-4.

Linearity

To determine the linearity a series of solutions were prepared using PIPER and TAZO working standards and their related substances at concentration levels from LOQ to 150% of specification and each solution was injected to HPLC. Correlation coefficient was calculated by plotting calibration curve using peak area ratio Vs concentration of the standard solution. Data was given in Table-5.

Specificity

All the solutions of related substances were prepared and injected individually to ensure that there is no interference of related substances at the retention times of PIPER and TAZO. Test

solution(control sample) and test sample spiked with known related substances were prepared and injected into HPLC. Data was recorded and given in Table-6.

Robustness

The robustness of an analytical method is a measure of its capacity to remain unaffected by small but deliberate variations in method parameters and provides an indication of its reliability during normal usage. Robustness was done by changing the column temperature ($\pm 5^{\circ}\text{C}$), flow rate ($\pm 10\%$), Changing the wavelength ($\pm 5\text{ nm}$). Standard solution and test solution spiked with known related substances at specification level were prepared and injected into HPLC at different variable conditions. The data was recorded and given in Table-7.

RESULTS AND DISCUSSION

Various HPLC, UV, Ion pair HPLC were developed and published for estimation of Piperacillin, Tazobactam and their related substances in bulk and pharmaceutical formulations. Present work was aimed to develop a simple, accurate and precise RP-HPLC method for the determination of related substances in Piperacillin and Tazobactam injection in a single run with less run time so that the method can be utilised in routine quality control laboratories.

Method development

Based on the results obtained from various trials detection wavelength was selected at 220 nm since PIPER, TAZO and their related substances showing maximum absorbance at around 200-240 nm. Hypersil ODS 250 mm \times 4.6 mm \times 5 μ column was selected as it shows better aqueous mobile phase compatibility. By using this column no merging of peaks was observed. The choice of buffer was made based on the P^{ka} values of Piperacillin and Tazobactam. Phosphate buffer of P^{H} (6 ± 0.05) was selected as it improves peak shape. A mixture of Phosphate buffer and acetonitrile (solution A) and pure acetonitrile(solution B) was selected as mobile phase because this separate all the peaks. on the basis of solubility of all the impurities as well as PIPER and TAZO, the diluents was selected to be 98:2%v/v:buffer:acetonitrile. Column oven temperature was set at 30°C to activate the column and to achieve the early separation. Hence the optimized HPLC parameters were as follows: flow rate 10ml/min, injection volume 20 μl and a gradient program with mobile phase which consists of solution A and solution B. The developed method separates

piperacillin, Tazobactam and their related substances by 45 min run time. No interference due to diluents and excipients was observed.

Method Validation

The method was validated as per ICH and USP guidelines.

System suitability

The percentage RSD for peak areas of six replicate injections of the standard solution should be less than 5.0%. USP plate count for piperacillin peak should not less than 5000 and USP tailing for PIPER should not more than 2.0 indicate the system was suitable for the proposed method.

Precision

System precision was calculated by using six replicates of standard solution. The percentage RSD should not be more than 5.0%. Method precision was carried out by using single batch of PIPER and TAZO for injection spiked with related substances at specification level. The percentage RSD should be less than 10.0%. The data was given in Table-2.

Accuracy

Recovery should be within 90.0% to 110.0% for related substances of 0.3 to 1.0% specification and 85.0% to 115.0% for related substances of 0.2 % specification indicates that the method is accurate.

Limit of Quantification(LOQ) and Limit of Detection(LOD)

LOQ and LOD values of each related substances were predicted from a separate linearity data. The percentage RSD should be within the acceptance limit of <10.0 for LOQ and 33.0% for LOD. The determined LOQ and LOD values for PIPER, TAZO and their related substances are presented in Table-4.

Linearity

Linearity curves of Piperacillin, Tazobactam and their related substances were calibrated by using a series of solutions at concentration levels of LOQ to % of specification. Correlation coefficient which represents linearity should meet acceptance limit of >0.990. The regression values are shown in Table-5, with linearity curves of Piperacillin, Tazobactam and their related substances are represented in Figure 7-12.

Specificity

Specificity of the method was demonstrated by injecting diluents, test sample spiked with known related substances so that peaks of all the components should meet the acceptance criteria of peak purity < purity threshold. Forced degradation studies were carried out to conform the proposed method has stability indicating capability and selectivity. Piperacillin and Tazobactam for injection was stressed with acid (2M HCl), base (0.05 M NaOH), peroxide (10% H₂O₂), thermal (105⁰C), photolytic (10 K Lux) and humidity (92% RH) conditions respectively. Peak purity for all the mentioned conditions was established and peak purity should meet acceptance criteria of peak purity < purity threshold. From all the peak purity data of PIPER and TAZO peaks at every degradation sample shows that the peaks are homogenous and there are no co-eluting peaks conforms that the method is stability indicating and specific.

Robustness

The robustness of an analytical method is a measure of its capacity to remain unaffected by small but deliberate changes. The proposed method is robust as it meets the acceptance criteria of USP plate count for PIPER <5000 and USP tailing for PIPER peak is <2.0.

Table -1 Results of System Suitability

System suitability parameters	Obtained value	Acceptance criteria
USP theoretical plates for Piperacillin peak	320802	USP plate count NLT 5000
Percentage RSD for PIPER from six replicate injections of standard solution	1.0%	Percentage RSD NMT 5.0%
Percentage RSD for PIPER from six replicate injections of standard solution	0.3%	Percentage RSD NMT 5.0%
USP tailing for Piperacillin peak	1.1	Tailing NMT 2.0

Table -2 Results for Precision and Intermediate Precision

Substance	Method precision		Intermediate precision	
	% Impurity	% RSD	% Impurity	% RSD
Ampicillin	0.236	0.8	0.221	0.5
D-phenylglycyl EDP	0.109	0.9	0.111	0.9
Piperacillin sulfoxide	0.258	0.4	0.231	0.0
Aminosulfinic acid	0.481	1.9	0.419	1.0

Table -3 Results for Accuracy

Substance	% Recovery (Mean of Three Replicates)		
	At 50%	At 100%	At 150%
Ampicillin	102.9	103.2	103.1
D-phenylglycyl EDP	103.1	103.4	102.1
Piperacillin sulfoxide	101.9	103.0	103.3
Aminosulfinic acid	98.5	96.1	96.0

Table -4 Results for LOD and LOQ




Substance	LOD($\mu\text{g/ml}$)	%RSD	LOQ($\mu\text{g/ml}$)	%RSD
Ampicillin	0.208	13.2	0.416	5.8
D-phenylglycyl EDP	0.230	12.9	0.460	9.0
Piperacillin sulfoxide	0.215	13.3	0.430	5.0
Aminosulfinic acid	0.230	13.8	0.460	7.9

Table -5 Results for Linearity

Substance	Linearity range ($\mu\text{g/ml}$)	Correlation coefficient R^2
Piperacillin	0.49 to 77.76	0.9997
Tazobactam	0.39 to 3.37	0.9993
Ampicillin	0.41 to 17.24	0.9999
D-phenylglycyl EDP	0.45 to 7.93	0.9992
Piperacillin sulfoxide	0.42 to 15.89	0.9996
Aminosulfinic acid	0.40 to 3.13	0.9998

Table -6 Results for Forced Degradation Studies

Mode of Degradation	Conditions	Peak purity of Piperacillin		Peak purity of Tazobactam	
		Purity angle	Purity threshold	Purity angle	Purity threshold
Acid degradation	2M HCl/initial	0.050	0.276	0.097	0.097
Base degradation	0.05M NaOH/initial	0.068	0.274	0.080	0.080
Peroxide degradation	10% H_2O_2 /85 $^\circ\text{C}$ /5 min	0.070	0.266	0.101	0.101
Thermal degradation	105 $^\circ\text{C}$ /120 Hours	0.107	0.296	0.524	0.524
Photolytic degradation	10K Lux/240 Hours	0.046	0.258	0.059	0.059
Humidity degradation	92% RH/ 25 $^\circ\text{C}$ /48 Hours	0.029	0.258	0.027	0.027

**Fig-7. Linearity Plot of Piperacillin****Fig-8. Linearity Plot of Tazobactam****Fig-9. Linearity Plot of Ampicillin**




Fig-10. Linearity Plot of D-Phenylglycyl EDP




Fig-11. Linearity Plot of Piperacillin Sulfoxide




Fig-12. Linearity plot of Aminosulfinic acid




Fig No-13 Chromatogram of Reference Solution (spiked with known related substances)

S. No	Name	RT(min)
1	Aminosulfinic acid	2.067
2	Tazobactam	6.651
3	Ampicillin	14.319
4	D-phenylglycyl EDP	17.934
5	Piperacillin sulfoxide	24.973
6	Piperacillin	28.00




Fig No-14 Chromatogram of System suitability

CONCLUSION

The present study describes development and validation of a RP-HPLC method for determination of Piperacillin and Tazobactam and their related substances in bulk and pharmaceutical formulation. The developed method is simple, accurate, sensitive, precise and can be effectively utilized for routine analysis in quality control laboratories, research institutions and clinical pharmacokinetic studies as this method separates PIPER and TAZO from their related substances by 45 min run time. From overall observations it was concluded that the developed method was more simple, accurate, sensitive, precise and found to be within prescribed limits.

ACKNOWLEDGEMENTS

Thankful to Aurabindo Pharma Pvt. Ltd., India for providing gift samples.

REFERENCES

1. Willard HH, Merritt LL, Dean JA and Settle FA. *Instrumental Methods of Analysis*. 7th ed. Delhi: CBS Publishers and Distributors. 2001; 3.
2. Skoog DA, West DM and Holler FJ. *Fundamentals of Analytical Chemistry*. 7th ed. Philadelphia: Saunders College Publishing. 1996; 1-3.
3. Sharma BK. *Instrumental Methods of Chemical Analysis*. 21st ed. Meerut: Goel Publishing House, 2002; 3-5.
4. Skoog DA, Holler FJ, Timothy A and Nieman NW (). *Principle of Instrumental Analysis*. 5th ed. Bangalore: Eastern Press. 2004; 1-2, 678-688, 695-696.
5. Scott RPW. *Technique and Practice of chromatography*. Marcel Dekker: Vol. 70. New York: 2003; 1-12.
6. Jeffery GH, Basset J, Mendham J and Denney RC. *Vogel's textbook of Quantitative Chemical analysis*. 5th ed. England: Longman Publication. 1996; 647-649.
7. Connors KA. *A textbook of Pharmaceutical Analysis*. 8th ed. New York: Wiley-Interscience. 1999; 408-421.
8. Hamilton RJ and Sewell PA. *Introduction to HPLC*. 2nd ed. London: Chapman and Hall. 1982; 189.
9. Snyder LR, Kirkland JJ and Glajch JL. *Practical HPLC Method Development*. 2nd ed. New York: Wiley. 1997; 1-20.
10. Sethi PD. 'High Performance Liquid Chromatography', *Quantitative Analysis of Pharmaceutical Formulations*, 1st ed. New Delhi: CBS Publishers and Distributors. 2001; 3-11, 116-120.
11. Munson JW. *Pharmaceutical Analysis: Modern Methods (Part B)*. New York: Marcel Dekker. 2001; 51-54,120,175.
12. Scott RPW. *Liquid Chromatography for the Analyst*. New York: Marcel Dekker. 1993; 67: 15-23, 265-272.
13. Chatwal GR and Anand SK. *Instrumental Methods of Chemical Analysis*. 5th ed. Delhi: Himalaya Publishing House. 2004; 2.599-2.605.
14. Ewing's. *Analytical Instrumentation Handbook*, 3rd ed. New York. Cazes J, Marcel Dekker. 2005; 995-998.

15. Ahuja, S.; Mills, K.; Handbook of isolation and characterization of impurities in pharmaceuticals, 1st ed.; Elsevier: California, 2003.
16. Sanjay B. Bari, B.R. Kadam, Y.S. Jaiswal, A. Shirkhedkar, Eurasian J of Anal Chem. 2007; 7: 135.
17. Gorog, S.; Mary, J. Journal of chromatography 2000; 5: 154.
18. ICH. impurities of Products Q3B (R2), International Conference on Harmonization, IFPMA, Geneva 2003.
19. Ahuja S Impurities Evaluation of Pharmaceuticals. Marcel Dekker, New York, 1998; 142.
20. Gorog S Identification and Determination of Impurities in Drugs. Elsevier Science Publishing Company, Amsterdam, 2000; 154.
21. ICH Harmonised Tripartite Guideline (Nov 2005) Validation of Analytical Procedures: Text and Methodology Q2 (R1).
22. ICH, (October 1993). Q1A Stability testing of New Drug Substances and Products in Proceedings of the International Conference on Harmonization, Geneva.
23. The United States Pharmacopoeia, The Official Compendia of Standards, (2006), 29th ed., Rockville, MD, USP convention Inc.
24. BJNN Sai*, S Mallikarjuna, DHH Theja, J Ravindra Reddy, Y Padmanabha Reddy, "Uv – Spectrophotometric Absorption Correction Method For The Estimation Of Piperacillin And Tazobactam In Injectables", Inventi Impact: Pharm Analysis & Quality Assurance, 2012.
25. C.H. Xia, Y.Q. Xiong and G.J. Wang, An improved highperformance liquid chromatographic method with a solidphase extraction for the determination of piperacillin and tazobactam: application to pharmacokinetic study of different dosage in Chinese healthy volunteers, Biomed. Chromatogr., 2007; 21: 680-686.
26. P.N.S. Pai, G.K. Rao, M.S. Murthy and H. Prathibha, Simultaneous estimation of piperacillin and tazobactam in injection formulations, Ind. J. Pharm. Sci., 2006; 68: 799-801.
27. Arzuaga, A. Isla, A.R. Gascon, J. Maynar, A. Martin, M.A. Solinis, D. Toral and J.L. Pedraz, Quantitative and stability of piperacillin and tazobactam in plasma and venous hemofiltration by HPLC, Biomed. Chromatogr., 2005; 19: 570- 578.
28. M.L. Qi, R.L. Chen, R.H. Cong and P. Wang, J. Liq. Chrom. Relat. Tech., A validated method for simultaneous determination of piperacillin sodium and sulbactam sodium in sterile powder for injection by HPLC, 2003; 26: 665-676.

29. R. Trittler, M. Ehrlich, T.J. Galla, R.E. Horch and K. Kummerer, New and rapid fully automated method for determination of tazobactam and piperacillin in fatty tissue and serum by column switching liquid chromatography, *J. Chromatogr. B*, 2002; 775: 127-132.
30. V. Augey, P.Y. Grosse, G. Albert, M. Audran and F. Bressolle, High performance chromatographic determination of tazobactam and piperacillin in human plasma and urine, *J. Chromatogr. B*, 1996; 682: 125-136.
31. A.Lakshmana Rao, K.Sai Krishna¹, Ch.Kiran Kumar, T.Raja. Simultaneous Determination Of Piperacillin And Tazobactam In Bulk And Pharmaceutical Dosage Forms By RPHPLC. *International Journal of Pharmacy and Pharmaceutical Sciences*. 2011; 3(2): 134-136.
32. M. Gandhimathi, M. Saravanakumar, T.K. Ravi Validated Ion Pair HPLC Method for Simultaneous Estimation of Ceftriaxone Sodium and Tazobactam Sodium in Dosage Form. *International J. of Pharma and Bio Sciences* 2010; 1(4): 17-22.
33. Wang P., Yang J., A liquid chromatographic method for simultaneous determination of amoxicillin sodium and Tazobactam in combination formulation, *J. Pharm. Biomed. Anal.*, 2004; 36: 565- 569.