

THE FUTURE IS CLEAR

Cassit Product Catalogue

December 2021



ABOUT CASSIT


Cassit develops and manufactures innovative transparent splints and casts for people with hand injuries or disabilities such as fractures, sprains, carpal tunnel syndrome, CP, stroke and more.

Our splints are made of unique thermoplastic material, comes with Neoprene straps and can be customized in 3 minutes to anyone, everywhere, easily, with only hot water and a towel.




1

SOFTEN



2

CUSTOMIZE



3

INSERT STRAPS

OUR BENEFITS



STRONG & FLEXIBLE



REMOLDABLE




LIGHTWEIGHT



SUSTAINABLE



WASHABLE



COMFORTABLE & AERATED



TRANSPARENT



EASY APPLICATION




WRIST ULNAR SPLINT

2.5mm


Customizable wrist splint with stiff lateral support. Provides appropriate support to the forearm, wrist and palm while enabling fingers movement. Straightens deformations to the available range of motion and prevents further deterioration.

Common Uses: Suspected nondisplaced Wrist Fracture, sprains, arthritis, medium-high muscle tone, patients with cerebral palsy or after stroke with an ability to use fingers. Post-surgical wrist fracture. Reduction of pain and inflammation at origin of the Extensor Carpi Radialis Brevis (ECRB).

Straps colors options:



CN	Size	Thickness	Measurements
10023	Medium	2.5 mm	223*116 mm
10025	Large	3 mm	262*135 mm



WRIST ULNAR SPLINT


WRIST THUMB SPLINT

2.5mm


Customizable wrist and thumb splint, reduces thumb motion and immobilize the wrist. Provides stabilization and protection of both the base of the thumb and the wrist while enabling fingers movement.

Common Uses: CMC arthritis, De Quervain's tenosynovitis, treatment of sprain of the thumb or wrist, prevention of wrist ulnar deviation and thumb flexion, reduction of pain, reduction of inflammation.

Straps colors options:



CN	Size	Thickness	Measurements
10049	Large	3 mm	265*172 mm



WRIST THUMB SPLINT



WRIST SPLINT

2.5mm

Customizable functional wrist splint. Brings wrist to optimal functional position while enabling fingers movement and grasping tools


Common Uses: Carpal Tunnel Syndrome, rheumatoid arthritis, reduction of wrist motion, reduction of pain, immobilization of wrist to prevent flexion.

4mm


Customizable functional wrist splint. Brings wrist to optimal functional position while enabling fingers movement and grasping tools

Common Uses: Sprains, Suspected nondisplaced Wrist Fracture, Acute stages of rheumatoid arthritis, Proper wrist positioning to increase maximum function, provide stability of wrist, allows full MP flexion and thumb opposition.

Straps colors options:



CN	Size	Thickness	Measurements
10017	Large	4 mm	237*117 mm
10015	Large	2.5 mm	237*117 mm
10011	Small	2.5 mm	182*91 mm



WRIST SPLINT




RESTING SPLINT


Customizable full length hand splint, from the forearm and up to the end of the fingers. Provides full support, protection and rest to the wrist and fingers in a neutral position. Ideal for the treatment of contractures, hand, and finger deformities for patients with spasticity.

Common Uses: Paralyzed hand and wrist, wrist drop, spastic hand after relief treatment, post-surgery- Extensor tendons, preventing loss of motion after an acute trauma such as burns, pain relief, preservation of ROM.

Straps colors options:



CN	Size	Thickness	Measurements
10003	Medium	3 mm	315*126 mm
10005	Large	4 mm	371*148 mm




RESTING SPLINT

FINGER SPLINT


Customizable finger splint, providing support and stability to the middle and upper finger joints. Immobilize the fingers in a straight or bent position and stabilize the phalanx joints.

Common Uses: Suspected Finger Fractures, sprains, cuts, dislocations PIP and DIP, drop finger, baseball finger, arthritis.

Straps colors options:



CN	Size	Thickness	Measurements
10032	Medium	1.8 mm	111*37 mm
10042	Large	1.8 mm	111*37 mm
10041	XL	1.8 mm	116*38 mm



FINGER SPLINT
Wide / Narrow




UNIVERSAL CUFF

Sponge tube with customizable plastic handle to assist grasping. An open clip adjusted to the palm shape to stabilize a tool in the hand and prevent it from falling. The Universal Cuff is the perfect aid to help gain complete control over the use of various tools such as eating utensils, toothbrush and more.

Common Uses: People with weakened/ loose grip, people with neurological-related disability, low fine motor skills, tremor, Parkinson's Disease.

CN	Size	Thickness	Measurements
10066	Small	3 mm	125*19 mm
10068	Large	3 mm	172*35 mm



UNIVERSAL CUFF

WRITING LOOP

Writing loop for wrapping around two fingers and the thumb. Suitable for learning proper grip/pinch of the writing tool.

Common Uses: People with motor disabilities of fingers, children with fine motor skills difficulties, children and adults with Cerebral Palsy, CVA or weak grip/pinch.

CN	Size	Thickness	Measurements
10501	Small	2 mm	40 mm
10503	Medium	2 mm	50 mm
10505	Large	2 mm	62 mm




Writing Loop




ACCESSORIES

CN	Prodcut Name	Size	Thickness	Measurements
10601	Wide Peripheral Strap	-	3 mm	180*50 mm
56001	Extra Straps Set Mix	-	1.5-3 mm	Various
10703	Heat Resistant Gloves	Medium	1.5-2 mm	22*130 mm
10705		Large	1.5-2 mm	25*150 mm
10707		Extra Large	1.5-2 mm	270*160 mm
10801	Cassit Towel	-	3 mm	500*300 mm
16013	Flat-End Forceps	-		234*63.5*26 mm
60002	Customization Tub	-		



Cassit Towel

Perfect size for protecting patients during customization, made of 100% cotton.



Customization Tub


Customization tub - reaches 100 degrees Celsius in 10 min. Lightweight and compact tub that fits all sizes of Cassit splints.

Supplied with a thick base to keep the splint away from the heating element.

Temperature range: 90°C - 100°C


Inner tub dimensions: 370 L x 300 W x 78 H mm

Outside dimensions: 500 L x 315 W x 205 H mm



Heat Resistant Gloves

Heat resistant gloves made of high quality Neoprene.



Flat-End Forceps

Flat-end forceps for gripping and lifting the splint while hot.



VISIT US ONLINE

www.cassit-ot.com

CONTACT US

CASSIT ORTHOPEDICS Ltd

26 Elifelet st. Tel Aviv, Israel

TEL: +972-542-999444

Info@cassit-ot.com

FOLLOW US

 [CassitIsrael](#)

 [cassit_ot](#)

