


VIVI CAP1 - Insulin Temperature Shield


By TempraMed



Keeps Insulin cool
everyday, everywhere
Years of use without
maintenance



- Protects insulin from overheating¹ or freezing
- FDA registered, scientifically validated device
- No preparation needed - just "click-n-go"
- Includes insulin temperature indicator
- Reusable insulin pen cap
- Small and lightweight - pocket fit
- Maintenance -free everyday use for years²
- Available for prefilled and refillable insulin pens



In testing³ VIVI CAP1 demonstrated superior performance over water evaporation based device x4 its weight.



Learn more at www.my-vivi.com

VIVI CAP1™ is manufactured and distributed by TempraMed Inc. FDA registration number 3012381649
19481 SE 14th St. Sammamish, WA 98075, USA, info@TempraMed.com, www.tempramed.com

¹ VIVI CAP1 keeps insulin cool within safe temperatures below 29°C (84.2°F) for minimum of 12 hours even in a constant environmental temperature of 37.8°C (100°F), when the device is placed inside a personal handbag.

² VIVI CAP1 is 'reactivated'/'filled up' to provide safe storage conditions, without any special user intervention, whenever the ambient temperature gets below 27°C (80.6°F), for subsequent periods of 12 hours and for up to 365 days of daily use, under limited warranty. Reactivation time varies and depends on the use conditions and ambient temperature. If needed VIVI CAP1 can be placed in the refrigerator without the insulin pen, to accelerate reactivation.

³ Pflutzner A. et al, Expert Opin Drug Deliv. 2017 Jun 14 (6) 709-714

