



ref. 54118M

COVID-19 ANTIGEN TEST KIT



Requirements

- Directive 98/79/EC on in vitro diagnostic medical devices

Packaging


Unit color box

Languages included: english, deutch, italian, spanish,

500 pcs per box (59,5x49x5x35cm)

Deepblue SL030101NST-1

Content



Use

In Vitro qualitative detection of the SARS CoV-2 antigen in human nasal swab specimen. Intended for personal use by untrained layman as rapid test method for novel coronavirus infection

Technical details

Sensitivity	96.4%
Specificity	99.8%
Storage	4°C – 30°C

Reference RT-PCR Assay				
		POS	NEG	Total
DEEP BLUE SARS-CoV-2 Ag Test	POS	106	1	107
	NEG	4	409	413
	TOTAL	110	410	520